

DISPODIPPER LABOPLAST® AND STERIPLAST®



Sampling tube manufactured in a clean room specially designed for disposable use. Due to the blue color, it is also perfect for use in the food, animal feed and pharmaceutical industries. The scaled sample container can be used for an exact measurement. After sampling, the container can be closed and is leakproof. The closed handle structure ensures secure handling.

The sample container is available with 50 ml, 100 ml or 250 ml (1.69 oz., 3.38 oz. or 8.45 oz.) content. Depending on the design, container openings can be sampled from a diameter of 52 mm (2" barrel opening), 65 mm (2.56 in.) or 85 mm (3.34 in.). This makes DispoDipper ideal for extracting samples from barrels, canisters and other containers.

DispoDipper is completely manufactured in one piece and can be used for sampling and subsequent transport of the sample in one device. In this way, cross-contamination or sample impurities are almost impossible.

For sampling, first unpack the sampler and detach the cover at the predetermined breaking point from the handle. Afterwards extract the sample using the sample container and securely close with the cover. After sampling, the handle at the predetermined breaking point can be detached from the sample container. This allows for easy and space-saving transport of the sample.



LaboPlast

SteriPlast

- ▶ PP, blue transparent
- ▶ Individually packaged and optionally sterilized by gamma rays
- ▶ Production and packaging according to clean room class 7
- ▶ Complies with EU food requirements and FDA regulations
- ▶ Scaled sampling tube volume in ml and oz.
- ▶ Immersion depth 270 mm (10.63 in.)

Capacity ml (oz.)	Model	Dia. mm	Height with cap mm	Total length mm	Pieces per pack	Price \$	Item no.
50 (1.69)	standard	40 (1.57)	98 (3.86)	417 (16.41)	20	\$28.46	76183-154
50 (1.69)	sterilized	40 (1.57)	98 (3.86)	417 (16.41)	20	\$35.96	76183-094
100 (3.38)	standard	53 (2.08)	98 (3.86)	415 (16.34)	20	\$40.81	76299-420
100 (3.38)	sterilized	53 (2.08)	98 (3.86)	415 (16.34)	20	\$51.04	76299-432
250 (8.45)	standard	72 (2.83)	116 (4.57)	429 (16.89)	20	\$46.33	76299-434
250 (8.45)	sterilized	72 (2.83)	116 (4.57)	429 (16.89)	20	\$57.87	76299-436