



Anachemia
A VWR Company

Product Specification

Material	74796-300
Material description	Potassium persulfate, crystals
Grade	ACS
CAS Number	7727-21-1
Molecular formula	K ₂ S ₂ O ₈
Molecular mass	270.32

Characteristics	Specifications
Assay	Min. 99.0 %
Chlorine compounds (as Cl)	Max. 0.001 %
Heavy metals (as Pb)	Max. 0.001 %
Insoluble matter	Max. 0.005 %
Fe (Iron)	Max. 5 ppm
Mn (Manganese)	Max. 2 ppm

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Michel Boisvert, Head of the Laboratory - Lachine
VWR Anachemia Canada, 255 rue Norman, Lachine, H8R 1A3,
QC, Canada

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.