


# PurFlock® Ultra Flocked Swab

Ultrafine tip, polystyrene handle with breakpoint at 100mm



## Description

An ideal choice in flock swabs with excellent absorption and elution of specimens when HydraFlock's high levels of absorbency aren't necessary. PurFlock Ultra's high purity characteristics makes it an ideal swab for various rapid diagnostic tests.

Patented: US 8420385

Patents pending: Canada 2,744,678, Europe 10251389.2

## Dimensions

6.02" Overall Length (152.908 mm)

### Tip

■ 0.140" Dia. (3.556 mm)

■ 0.66" Length (16.76 mm)

### Handle

■ 0.1" Dia. (2.54 mm) tapered to 0.036" (0.914 mm)

■ 5.99" Length (152.146 mm)

## Standards

■ CE marked for medical use

■ ISO 9001:2008

■ ISO 13485:2003

■ FDA GMP's

■ Certificate of sterility and/or conformance available with each shipment upon request

## Qualities

■ Medical grade quality

■ Highly absorbent tip

■ Soft smooth tip

## Applications

■ For nasopharyngeal and urethral cell collection

■ Rapid diagnostic testing

■ Specimen collection

■ Forensic evidence collection

| Ref #                               | Packaging                                                                            |
|-------------------------------------|--------------------------------------------------------------------------------------|
| <b>3317-U</b><br><i>Non-sterile</i> | 10/100<br>(10 bags per case, 100 per bag)                                            |
| <b>25-3317-U</b><br><i>Sterile</i>  | 10/50/1<br>(10 boxes per case, 50 individually wrapped packages, 1 swab per package) |



**PurFlock®  
Ultra**

Puritan Medical Products Company LLC  
Guilford, Maine 04443-0149 USA

Tel 800-321-2313 (US and Canada Only)  
207-876-3311

Fax 800-323-4153  
207-876-3130

Web [www.puritanmedproducts.com](http://www.puritanmedproducts.com)

Email [sales@puritanmedproducts.com](mailto:sales@puritanmedproducts.com)

**1-800-321-2313**  
**[puritanmedproducts.com](http://puritanmedproducts.com)**