


Component Schematic

96-well Microplate (PP) Deep Well (Round) raised rim



L	127.90 ± 0.10	Hw	39.00 ± 0.05
W	85.80 ± 0.10	P	8.95 ± 0.05
H	42.10 ± 0.05	D	6.60 ± 0.05
Hs	3.25 ± 0.10	C	12
X1	14.70 ± 0.05	R	8
Xn	113.20 ± 0.05		
Y1	11.50 ± 0.05		
Yn	74.30 ± 0.05		

Volume 1,300µl

Well form
Material

U Bottom
PP