





NA Cat	EU Cat#	Short Description English	Tube Capacity	Number of places	Slot Size	Width (mm)	Length (mm)	Height (mm)	Working Temperature	Vents	Product Color
89235-996	479-2419	VWR FREEZER BOX 50TUBE MIX PK5	1-2ml	50	12.5	71.6	139.8	56.7	-80°C to 121°C	NA	Mix
89235-998	479-2420	VWR FREEZER BOX 50TUBE BLU PK5	1-2ml	50	12.5	71.6	139.8	56.7	-80°C to 121°C	NA	Blue
89236-000	479-2421	VWR FREEZER BOX 50TUBE GRN PK5	1-2ml	50	12.5	71.6	139.8	56.7	-80°C to 121°C	NA	Green
89236-002	479-2422	VWR FREEZER BOX 50TUBE NTL PK5	1-2ml	50	12.5	71.6	139.8	56.7	-80°C to 121°C	NA	Natural
89236-004	479-2423	VWR FREEZER BOX 50TUBE ORG PK5	1-2ml	50	12.5	71.6	139.8	56.7	-80°C to 121°C	NA	Orange
89236-006	479-2424	VWR FREEZER BOX 50TUBE PRL PK5	1-2ml	50	12.5	71.6	139.8	56.7	-80°C to 121°C	NA	Purple
89236-008	479-2425	VWR FREEZER BOX100TUBE MIX PK5	1-2ml	100	12.5	140.6	140.6	57.2	-80°C to 121°C	NA	Mix
89236-010	479-2426	VWR FREEZER BOX100TUBE BLU PK5	1-2ml	100	12.5	140.6	140.6	57.2	-80°C to 121°C	NA	Blue
89236-012	479-2427	VWR FREEZER BOX100TUBE GRN PK5	1-2ml	100	12.5	140.6	140.6	57.2	-80°C to 121°C	NA	Green
89236-014	479-2428	VWR FREEZER BOX100TUBE NTL PK5	1-2ml	100	12.5	140.6	140.6	57.2	-80°C to 121°C	NA	Natural
89236-016	479-2429	VWR FREEZER BOX100TUBE ORG PK5	1-2ml	100	12.5	140.6	140.6	57.2	-80°C to 121°C	NA	Orange
89236-018	479-2430	VWR FREEZER BOX100TUBE PRL PK5	1-2ml	100	12.5	140.6	140.6	57.2	-80°C to 121°C	NA	Purple
89236-020	479-2431	VWR FREEZER BOX 81TUBE MIX PK5	1-2ml	81	12.5	126.5	126.5	56.7	-80°C to 121°C	NA	Mix
89236-022	479-2432	VWR FREEZER BOX 81TUBE BLU PK5	1-2ml	81	12.5	126.5	126.5	56.7	-80°C to 121°C	NA	Blue
89236-024	479-2433	VWR FREEZER BOX 81TUBE GRN PK5	1-2ml	81	12.5	126.5	126.5	56.7	-80°C to 121°C	NA	Green
89236-026	479-2434	VWR FREEZER BOX 81TUBE NTL PK5	1-2ml	81	12.5	126.5	126.5	56.7	-80°C to 121°C	NA	Natural
89236-028	479-2435	VWR FREEZER BOX 81TUBE ORG PK5	1-2ml	81	12.5	126.5	126.5	56.7	-80°C to 121°C	NA	Orange
89236-030	479-2436	VWR FREEZER BOX 81TUBE PRL PK5	1-2ml	81	12.5	126.5	126.5	56.7	-80°C to 121°C	NA	Purple

50 Places Box



81 / 100 Places Box



Application Temperature: -80°C to +121°C

Material Recycle Code



Description: Freezer Box

Drawing:
Not to Scale



Date: 10/16/2023

Revision: 1