

3200 Secure II Standard Features

Biology

From the highest quality living and preserved specimens to the widest selection of media, microscope slides, and more, we have more options for your biology classroom than anywhere else.



- New Hospital Grade 84" Long x 36" wide X 6" thick foam mattress
- Three-sided backrest (fowler)/knee gatch controls
- Caregiver siderail controls for hi/lo, head up/head down & knee up/knee down positions
- Patient siderail controls for backrest and knee positions
- Caregiver footboard control panel for hi/lo, backrest (fowler), knee gatch, cardiac chair-in and chair-out, Trendelenburg and reverse trend positions as well as patient siderail control lockouts. (Optional: Integrated scale controls and display and Chaperone® Bed exit controls)
- Removable footboard and headboard. (Headboard can be upgraded for emergency CPR.)
- One-hand, fold-down head and foot end siderails
- Open architecture – flat litter deck
- StayPut Technology – Retractable bed frame
- Automatic contour positioning
- 6" (cm) casters
- Four-wheel, steel brake with swivel lock system
- Centrally-located activation pedals
- Locking steer control caster for improved steering and mobility with activation pedal at head end of bed.
- Downward motion interrupt
- Night light
- Instant CPR release
- Four foley bag hooks
- Foot prop rod
- Trendelenburg and Reverse Trendelenburg controls
- Foot end digital trend and backrest (fowler) angle display)
- Four IV pole and traction equipment socks
- Backrest (fowler) angle gauge
- Hospital grade power plug
- DC motor

Specifications:

Overall Length: 93" (236 cm)

Overall Width:

- Siderails Up: 42.5" (108 cm)
- Side
- rails Down: 40" (102 cm)

Weight Capacity: 500 lbs (228 kg)

Height Range (to litter top):

- High: 30" (76 cm)
- Low (standard with 6' casters: 16" (41 cm)

Litter Positioning:

- Backrest: 0° – 60°
- Knee Gatch: 0° – 40°
- Trend/reverse Trend: $\pm 12^\circ$
- Retraction: 14" (36 cm)

Patient Surface: 35" x 84" (89 cm x 213 cm)

Caster Diameter:

- Standard: 6" (15 cm)
- Optional: 8" (20 cm)

