



Technical Data Sheets

Kimberly-Clark* Classic* Surgical Mask

Pleat style face mask with ties.

48201

62369

EN14683 Type I



Catalog Numbers:

Classic* Surgical Mask, Blue

- Green

Code # **48201**
(formerly 47123-010,
80-250, 80-410,
37550, 375501,
48200)
Code # **62369**
(formerly 48205)

Performance Data:

Sub-Micron Particle Filtration Efficiency (0.1 micron):	$\geq 98\%$ ¹
Differential Pressure, or "Delta P" (Breathability):	$< 20 \text{ Pa/cm}^2$ ²
Bacterial Filtration Efficiency (BFE in vitro, 2.9 micron):	$\geq 96\%$ ²

Composition:

Polypropylene, Polyester, Cellulose, Aluminum (nose wire)

Notes:

Natural Latex Free, Non-Sterile, Single Use Only

Kimberly-Clark Corporation
Health Care
1400 Holcomb Bridge
Roswell, Georgia USA 30076

Belgicastraat 13
B-1930 Zaventem
Brussels, Belgium
Tel: +32 (0) 2 711 26 54
Fax: +32 (0) 2 711 26 90

Commitment to Excellence

If, for any reason, our products do not meet your expectations, please let us know your comments or suggestions for improvement. Your input will result in a concerted effort on our part to meet your requirements. Our goal is to provide quality products that completely meet your needs time after time.



For Further information contact:
Partners in Quality* at 1-800-KCHELPS
(1-800-524-5377)

¹ Nelson Laboratories, 6280 South Redwood Rd., Salt Lake City, UT 84123. Latex Particle Challenge Test. Report No. 286242

² Nelson Laboratories, Differential Pressure Test. Report No. 286239

* Registered Trademark or Trademark of Kimberly-Clark Worldwide, Inc. or its affiliates. ©2006 KCWW