

PRODUCT DATA SHEET

PART NUMBER:414004-266

PART DESCRIPTION:Aspirating Pipet-5ml

MATERIAL

Polystyrene

PRODUCT REQUIREMENTS

Non-pyrogenic, Gamma Irradiated

Dimensions(mm)

Tip:

Length=37.6

Interior Dia.=1.53

Outside Dia.=3.3

PIPET SIZE:5ml

INCREMENTS:N/A

GRADUATIONS:N/A

Mouthpiece:

Length=N/A

Interior Dia.=N/A

Outside Dia.=N/A

Plug Length=N/A

OVERFILL CAPACITY:N/A

COLOR CODE:N/A

Extrusion

Interior Dia.=6

Outside Dia.=8


Length=N/A

GENERAL DIMENSIONS

A=347

B=3.3

C=8



PACKAGING REQUIREMENTS

Individually wrapped packaging,50/bag,200/case

INSPECTION REQUIREMENTS

Representative production samples are collected and inspected in accordance with current applicable product specifications.