

Product Name:	FRV-30
Description:	<p> VWR Catalog No: 10819-766 VWR Standard Flammable Storage Refrigerator 30 Cu. Ft. Capacity One solid swing door, right hinged, not reversible Mechanical thermostat Adjustable operating temperature range: 1°C – 10°C Five adjustable shelves Keyed door lock Magnetic door gasket for positive seal Heavy duty door handle Leveling legs Probe access port (3/8") White textured cabinet and door High density urethane foam cabinet and door insulation CFC-free refrigerant (R404a), foam insulation and packaging Two year parts and labor warranty Exterior dimensions: 35 ¾" W x 35 ½" D x 72 ½" H Cycle defrost 115V, 60Hz, 6 Amps, 1/3 HP UL/C-UL listed Units meet NFPA guidelines 45 and 70 Shipping weight: 381 lbs. </p>