



Interior

Front View



Side View


115VAC/60Hz/1PH
NEMA-5-15P

Cat# 76210-394/HCUCBI0420LH

Cu.Ft.	Defrost	Door	Int Door	Shelves	Interior Dimensions			Refrigerant	H.P.	Amps	Weight
					W"	D"	H"				
4.2	Cycle	1 Solid LH	0	2 Fixed	19	17 1/4	24 1/2	R600a	1/8	1.5	100 lbs.

Note: Interior dimensions are estimates and subject to change without notice