

PALL AcroPrep™ Advance 350
DNA Binding
PN 8032 Lot No. FZ3640 (18) 82617

PALL AcroPrep™ Advance
PN 8071 Lot No. FZ3640

PALL AcroPrep™ Advance
PN 8132 DNA Binding
Lot No. FZ5120 (18) 82544

FOR A LIMITED TIME

Get 25% off AcroPrep™ Advance 96-well Filter Plates

- AcroPrep Advance 96-well Filter Plates can be used across a wide range of laboratory applications including; multiplexing, lysate clearance, protein purification, DNA purification, ultrafiltration and solvent filtration.
- Available in well volumes of 350 µL, 1 mL and 2 mL deep well formats.
- Range includes Pall's AcroPrep Advance 96-well Long Tip Filter Plate for Nucleic Acid Binding. This plate incorporates a silica-based quartz glass fiber media to allow for the efficient binding of DNA and RNA, while providing smooth flow and rapid processing of samples.

Ordering Information for AcroPrep Advance 96-Well Filter Plates

Protein Purification

VWR Cat. No.	Pall Part No.	Description	Pkg
MP89135-684	8071	350 µL, Mustang® Q membrane	10/pkg
MP97065-368	8022	350 µL, Mustang Q membrane	10/pkg
MP89135-686	8072	350 µL, Mustang S membrane	10/pkg
MP89135-704	8171	1 mL, Mustang Q membrane	5/pkg
MP89135-706	8172	1 mL, Mustang S membrane	5/pkg

Nucleic Acid Binding (NAB)

MP10158-728	8133	1 mL	5/pkg
-------------	------	------	-------

DNA Purification

MP89137-182	8032	350 µL, DNA binding	10/pkg
MP89135-690	8132	1 mL, DNA binding	5/pkg

Lysate Clearance – Supor® Membrane

VWR Cat. No.	Pall Part No.	Description	Pkg
MP97052-116	8075	350 µL, 3 µm glass fiber/0.2 µm	10/pkg
MP97052-114	8040	350 µL, 3 µm glass fiber/1.2 µm	10/pkg
MP97052-130	8175	1 mL, 3 µm glass fiber/0.2 µm	5/pkg
MP97052-132	8275	2 mL, 3 µm glass fiber/0.2 µm	5/pkg

Ultrafiltration

MP97052-106	8033	350 µL, Omega™ 3K MWCO	10/pkg
MP97052-108	8034	350 µL, Omega 10K MWCO	10/pkg
MP97052-110	8035	350 µL, Omega 30K MWCO	10/pkg
MP97052-112	8036	350 µL, Omega 100K MWCO	10/pkg
MP89135-696	8163	1 mL, Omega 3K MWCO	5/pkg
MP89135-698	8164	1 mL, Omega 10K MWCO	5/pkg
MP89135-700	8165	1 mL, Omega 30K MWCO	5/pkg
MP89135-702	8166	1 mL, Omega 100K MWCO	5/pkg

Solvent Filtration – wwPTFE Membrane

VWR Cat. No.	Pall Part No.	Description	Pkg
MP76308-642	8582	350 µL, 0.2 µm, short tip	10/pkg
MP76308-644	8584	350 µL, 0.45 µm, short tip	10/pkg
MP76308-646	8682	1 mL, 0.2 µm, short tip	5/pkg
MP76308-648	8684	1 mL, 0.45 µm, short tip	5/pkg
MP76308-650	8782	2 mL, 0.2 µm, short tip	5/pkg
MP76308-652	8784	2 mL, 0.45 µm, short tip	5/pkg
MP76423-864	8586	350 µL, 0.2 µm, long tip	10/pkg
MP76423-866	8686	1 mL, 0.2 µm, long tip	5/pkg
MP76423-868	8588	350 µL, 0.45 µm, long tip	10/pkg
MP76423-870	8688	1 mL, 0.45 µm, long tip	5/pkg

Solvent Filtration – PTFE Membrane

MP89137-184	8047	350 µL, 0.2 µm, short tip	10/pkg
MP89137-186	8048	350 µL, 0.45 µm, short tip	10/pkg
MP89135-692	8147	1 mL, 0.2 µm, short tip	5/pkg
MP89135-694	8148	1 mL, 0.45 µm, short tip	5/pkg
MP89135-710	8247	2 mL, 0.2 µm, short tip	5/pkg
MP89135-712	8248	2 mL, 0.45 µm, short tip	5/pkg

Aqueous Filtration

VWR Cat. No.	Pall Part No.	Description	Pkg
MP97052-096	8019	350 µL, 0.2 µm Supor membrane	10/pkg
MP97052-098	8029	350 µL, 0.45 µm Supor membrane	10/pkg
MP97052-100	8039	350 µL, 1.2 µm Supor membrane	10/pkg
MP97052-102	8027	350 µL, 30-40 µm PP/PE non-woven media	10/pkg
MP97052-104	8031	350 µL, 1 µm glass fiber	10/pkg
MP97052-124	8119	1 mL, 0.2 µm Supor membrane	5/pkg
MP97052-126	8129	1 mL, 0.45 µm Supor membrane	5/pkg
MP97052-128	8130	1 mL, 1.2 µm Supor membrane	5/pkg
MP89135-688	8131	1 mL, 1 µm glass fiber	5/pkg
MP30621-080	8127	1 mL, 30-40 µm PP/PE non-woven media	5/pkg
MP89135-708	8231	2 mL, 1 µm glass fiber	5/pkg
MP30621-082	8227	2 mL, 30-40 µm PP/PE non-woven media	5/pkg

Multiplexing

MP97052-118	8049	For multiplex assays	10/pkg
-------------	------	----------------------	--------

Receive up to 25% off of VWR published list price AcroPrep Advance 96-well filter plates between January 1, 2022 to June 30, 2022. Products Ordered must be from list associated for this promotion. Promotion ends June 30, 2022. Visit vwr.com/pall for information on qualifying products. For further details contact your VWR or Pall rep.

Promotional part numbers are available for a limited time, must be used at the time of purchase, and may be ordered via the “ORDER ENTRY” link on vwr.com. (Promotional part numbers cannot be found via Search on vwr.com.) For assistance with ordering, please contact VWR Customer Service at 1 800 932 5000.



Ordering:

vwr.com

© Copyright 2022, Pall Corporation. Pall, , AcroPrep, Omega, Seitz and Supor are trademarks of Pall Corporation. ® Indicates a trademark registered in the USA.

GN21.0705

1/22


VWR.COM

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by country, federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. Offers valid in countries listed above, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2022 Avantor, Inc. All rights reserved.

0122 Lit. No. 160249WREV Promo Code 5784