



VWR® 25cm Cell Scraper

Description: 25cm, Cell Scrapers, Sterile, Individually wrapped

Purpose: The process of scraping off and collecting cells more easier and effective

Materials


Handle: ABS (Acrylonitrile Butadiene Styrene)

Color: Green

Blade: TPE (Thermo Plastic Elastomer)

Color: Natural

Dimensions (Unit: mm)



Features

- Temperature range: -20°C to +60°C
- Autoclavability: No
- DNase/RNase free and Non-pyrogenic
- Sterilized by irradiation SAL10⁻⁶ (ISO11137)
- Shelf life: 3 years after month of production
- Manufactured in a class 100,000 room environment
- Manufactured under ISO13485 and ISO9001 quality management system

VWR® 25cm Cell Scraper

| VWR EU ENA. No. | VWR NA Cat. No. | Description | Sterile | Qty/Tray | Qty/Pk | Qty/Case |
|--------------------|--------------------|--|---------|----------|--------|----------|
| 734-2602 | 10062-904 | Cell scrapers 25cm length, scraper blade 20mm, sterile/1 | Yes | 1 | 100 | 100 |
| 734-2603 | 10062-906 | Cell scrapers 25cm length, lifter blade 20mm, sterile/1 | Yes | 1 | 100 | 100 |



VWR® 39cm Cell Scraper

Description: 39cm, Cell Scrapers, Sterile, Individually wrapped

Purpose: The process of scraping off and collecting cells more easier and effective

Materials


Handle: ABS (Acrylonitrile Butadiene Styrene)

Color: Green

Blade: TPE (Thermo Plastic Elastomer)

Color: Natural

Dimensions (Unit: mm)



Features

- Temperature range: -20°C to +60°C
- Autoclavability: No
- DNase/RNase free and Non-pyrogenic
- Sterilized by irradiation SAL10⁻⁶ (ISO11137)
- Shelf life: 3 years after month of production
- Manufactured in a class 100,000 room environment
- Manufactured under ISO13485 and ISO9001 quality management system

VWR® 39cm Cell Scraper

| VWR EU ENA. No. | VWR NA Cat. No. | Description | Sterile | Qty/Tray | Qty/Pk | Qty/Case |
|--------------------|--------------------|--|---------|----------|--------|----------|
| 734-2604 | 10062-908 | Cell scrapers 39cm length, scraper blade 30mm, sterile/1 | Yes | 1 | 100 | 100 |
| 734-2605 | 10062-910 | Cell scrapers 39cm length, lifter blade 30mm, sterile/1 | Yes | 1 | 100 | 100 |