

PDS No. 73117x	PRODUCT DATA SHEET				Page 1 of 1
Revision 02	Disposable Inoculation Loop, 10 µl, PS				 greiner bio-one
	Greiner Item-No. 73117x				
Valid for Item-No.:	731170 (sterile)	731171 (sterile)	731175 (sterile)		

1.	Description / Specification	
1.1	Description	PS flexible disposable inoculation loop, 10 µl
1.2	Dimensions	Length: 200 mm
1.3	Volume	Volume (wetting): 10 µl
1.4	Material / Resin	PS (Polystyrene)
1.5	Colour	Blue
1.6	Sterilization	SAL 10 ⁻³
1.7	Quality Control	- <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	Hydrophobic
2.2	Temperature range	-20°C to +60°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: www.gbo.com/bioscience →Products →Literature →Technical Information →Chemical Resistance of Resins
2.6	Shelf life	5 years after month of production
2.7	Other Information	-

3.	Packaging	731170	731171	731175
3.1	Pieces / Bag	50	1	10
3.2	Pieces / Box	2000	600	3000
3.3	Lot-No.	PPYWWX (working place, year, week, index)	PPYWWX (working place, year, week, index)	PPYWWX (working place, year, week, index)
3.4	Other Information	-	-	

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 01	Date 14 November 2006	Date 15 November 2006	Date 20 November 2006	
Date 15.03.05	Name S. Kaelberer	Name Dr. L. Marchetti	Name A. Schulz	