

Certificate Of Quality

Product Information

Product Number: 10600023
Product Name: Amersham™ Hybond™ P 0.45 PVDF
Lot Number: A28278266
Converting Lot: A28843305
Format: 300 mm x 4 m

Tests

Nominal Pore Size: 0.45 µm

Test	Range
Thickness [µm]	165 – 190
Flow rate with 2-Propanol, Ø 50mm [s/100ml]	50 – 140
Bubble point with 2-Propanol [bar]	0.9 – 1.1
Tensile Strength [N/15mm]	≥ 3.5

Conformance & Quality systems statement

This is to certify that this product conforms to GE Healthcare Life Sciences specifications.

All products are manufactured via a Quality System certified to **ISO 9001** and **EN ISO 13485** and tested in accordance with documented quality procedures and approved as a result of meeting the required specification.

Quality data are stored.

Electronic signature

This document has been electronically produced and is valid without a signature.

Dassel, 05-Apr-2019 / AHe

Version AB



imagination at work

www.gelifesciences.com

GE Healthcare UK Limited
Amersham Place Little Chalfont
Buckinghamshire HP7 9NA UK