


AC-1212xx	CUSTOMER DRAWING	Page 1 of 1
Revision 03	Cryo.s™ 2 ml, round bottom	 greiner bio-one
	Item-No. 121 2xx	

NOT TO SCALE



Cat.-No.	Description	A	B	C	D	E	F
121 261	Cryo.s™ 2 ml, natural	Ø 10.6	Ø 12.5	Ø 12.6	38.0	47.85	17.0
121 263	Cryo.s™ 2 ml, natural	Ø 10.6	Ø 12.5	Ø 12.6	38.0	47.85	17.0
121 277	Cryo.s™ 2 ml, green	Ø 10.6	Ø 12.5	Ø 12.6	38.0	47.85	17.0
121 278	Cryo.s™ 2 ml, yellow	Ø 10.6	Ø 12.5	Ø 12.6	38.0	47.85	17.0
121 279	Cryo.s™ 2 ml, blue	Ø 10.6	Ø 12.5	Ø 12.6	38.0	47.85	17.0
121 280	Cryo.s™ 2 ml, red	Ø 10.6	Ø 12.5	Ø 12.6	38.0	47.85	17.0

All dimensions in mm

Customer drawing subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 02	Date November 20, 2008	Date November 20, 2008	Date November 21, 2008	
Date November 17, 2006	Name F. Falter	Name F. Sonntag	Name F. Sonntag	