

Boreal2 AP Stereomicroscope

| | |
|----------------------------|---|
| | 470014-538 |
| Eyepiece | Paired 10x, 20mm |
| Head | Fixed 45° inclined head |
| Magnification | 2X 4X |
| Stage | Frosted glass and reversible black & white |
| Working distance | 89mm |
| Focusing | Coarse focus with tension adjustment and slip clutch |
| Top illumination | Corded LED |
| Bottom illumination | Corded LED |
| Replacement bulb | (T) Item #470330-980 (B) #470330-978 |
| Product includes | Removable soft eyecups and vinyl dust cover |
| Dimensions | 12" (H) x 7" (W) x 9" (D) |
| Warranty | Boreal Lifetime (includes shipping, excludes bulb & fuse) |

Boreal



**Lifetime
Warranty**

Boreal Microscope Warranty; Best in the Industry

Life-Time Warranty; No-Cost, No Hassle Warranty

Our extraordinary, all inclusive lifetime warranty covers all Boreal Microscope repairs beginning from the day you first receive your Boreal Microscope. It covers any defect in materials or workmanship for the life of the microscope with no restrictions on ownership. If anything goes wrong with your Boreal Microscope (with exception of light bulbs, fuses, routine cleaning and normal wear), we repair or even replace it, if necessary (at our discretion).

All repairs are done by our in-house staff of trained technicians. All repairs are performed at no cost to you. There is no charge for parts or labor. If your Boreal Microscope needs to be shipped back to us, we provide everything you will need to do so. We even cover the cost of shipping and handling - both from you to us, and back. There is never a charge for any of these services.