

| | |
|---------------|--|
| Product Name: | PR051WWW/0 |
| Description: | <p>Corepoint Scientific™ White Diamond Series Undercounter Refrigerator Built-In 4.6 Cu. Ft. Capacity</p> <p>One solid swing door, right hinged*</p> <p>Microprocessor temperature controller</p> <p>Digital temperature display</p> <p>Audible and visual high and low temperature alarms</p> <p>Remote alarm contacts</p> <p>2 probes; 1 in air and 1 in glycol bottle</p> <p>Adjustable operating temperature range: 1°C – 10°C</p> <p>Three shelves (two adjustable/one fixed) with guard rail on back to prevent product from being pushed against the back wall</p> <p>Back wall evaporator cover</p> <p>Forced draft circulation</p> <p>Recessed door handle</p> <p>Keyed door lock</p> <p>Magnetic door gasket for positive seal</p> <p>Leveling legs</p> <p>Probe access port (3/8")</p> <p>Powder coat exterior</p> <p>High density urethane foam cabinet and door insulation</p> <p>HFC-free refrigerant (R600a), foam insulation and packaging</p> <p>Three year parts and labor warranty, plus an additional two year compressor parts only warranty</p> <p>Pyxis®, Omnicell® and AcuDose RX® compatible</p> <p>Exterior dimensions: 23 ¾ "W x 24 ½" D x 33 3/8 "H</p> <p>Cycle defrost</p> <p>115V, 60 Hz, 1/8 HP, 1.74 Amps</p> <p>UL/C-UL listed</p> <p>Shipping weight: 137 lbs.</p> <p>*This product is available left hinged, product number PR051WWWLH/0</p> |