



Anachemia
A VWR Company

Product Specification

Material	62786-460
Material description	Nitric acid
Grade	ACS
CAS Number	7697-37-2
Molecular formula	HNO ₃
Molecular mass	63,01

Characteristics	Specifications
Assay	68,0 --> 70,0 %
Appearance	Clear colourless liquid
Free from suspended matter or sediment	Passes test
Color	Max. 10 APHA
Heavy metals (as Pb)	Max. 0,2 ppm
Residue after ignition	Max. 5 ppm
Cl (Chloride)	Max. 0,5 ppm
SO ₄ (Sulphate)	Max. 1 ppm
As (Arsenic)	Max. 0,01 ppm
Fe (Iron)	Max. 0,2 ppm

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Alicia Bourdeau, Head of the Laboratory - Rouses Point
VWR Chemicals, LLC, 3 Lincoln Blvd, Rouses Point, NY
12979-1045, USA