



VWR® Carboys

Polypropylene Carboys

Use for storing and dispensing solutions and media. Ideal for sterile water storage, polypropylene (PP) carboys can be autoclaved before refilling to stop recurring bacterial growth. Excellent chemical resistance. Graduated to identify the volume dispensed. Carboys manufactured with 100% virgin Medical grade resin meets US FDA & USP Class VI Requirements.

HDPE Carboys

The high-density polyethylene (HDPE), disposable carboys are ideal for storing and transfer of sterile fluids and pharmaceutical/biotech reagents. Carboys have both wide mouth and narrow mouth option to suite the application. Carboys have excellent chemical resistance and are well suited for use from –100°C to 100°C. Sterile Carboy available for high purity application. Graduated to identify the volume dispensed. Carboys manufactured with 100% virgin Medical grade resin meets US FDA & USP Class VI Requirements.



- **Carboy:** Refer to table for material type, 100 % Virgin, FDA 21 CFR Compliance, USP Class VI passed.
- **Closure:** Polypropylene, 100% Virgin, FDA 21 CFR Compliance, USP Class VI passed with food-grade color..
- **BSE/TSE Statement:**—These products are manufactured without the use of raw materials of animal origin.
- **Latex Free Statement:** No materials containing latex or natural rubber are used in the manufacturing, handling and packaging processes for these products.



VWR® Carboys

VWR® Rectangular Carboy with Spigot, Wide Mouth						
Capacity	Material	Cap Size	Neck ID	Sterility	Cat. No.	UOM
10 L (2.6 gal.)	HDPE	100 mm	85 mm	Nonsterile	10755-094	Each
10 L (2.6 gal.)	PP	100 mm	85 mm	Nonsterile	10755-096	Each
20 L (5.3 gal.)	HDPE	100 mm	85 mm	Nonsterile	10755-106	Each
20 L (5.3 gal.)	PP	100 mm	85 mm	Nonsterile	10755-108	Each
VWR® Round Carboy, Wide Mouth, Without Spigot						
10 L (2.6 gal.)	LDPE	100 mm	85 mm	Nonsterile	10755-100	Each
10 L (2.6 gal.)	PP	100 mm	85 mm	Nonsterile	10755-102	Each
20 L (5.3 gal.)	PP	100 mm	85 mm	Nonsterile	10755-122	Each
VWR® Round Carboy, Narrow Mouth, Without Spigot						
10 L (2.6 gal.)	HDPE	83 mm	65 mm	Nonsterile	10755-098	Each
10 L (2.6 gal.)	LDPE	83 mm	65 mm	Nonsterile	10755-084	Each
10 L (2.6 gal.)	PP	83 mm	65 mm	Nonsterile	10755-086	Each
20 L (5.3 gal.)	HDPE	83 mm	65 mm	Sterile	10755-121	Each
20 L (5.3 gal.)	HDPE	83 mm	65 mm	Nonsterile	10755-110	Each
20 L (5.3 gal.)	LDPE	83 mm	65 mm	Nonsterile	10755-112	Each
20 L (5.3 gal.)	PP	83 mm	65 mm	Nonsterile	10755-114	Each
50 L (13.2 gal.)	LDPE	83 mm	65 mm	Nonsterile	10755-124	Each
50 L (13.2 gal.)	PP	83 mm	65 mm	Nonsterile	10755-126	Each
VWR® Heavy-Duty Round Carboy, Narrow Mouth, Without Spigot						
10 L (2.6 gal.)	PP	83 mm	65 mm	Nonsterile	10755-092	Each
20 L (5.3 gal.)	PP	83 mm	65 mm	Nonsterile	10755-104	Each
VWR® Round Carboy, Narrow Mouth, With Spigot						
10 L (2.6 gal.)	LDPE	83 mm	65 mm	Nonsterile	10755-088	Each
10 L (2.6 gal.)	PP	83 mm	65 mm	Nonsterile	10755-090	Each
20 L (5.3 gal.)	LDPE	83 mm	65 mm	Nonsterile	10755-116	Each
20 L (5.3 gal.)	PP	83 mm	65 mm	Nonsterile	10755-118	Each
50 L (13.2 gal.)	LDPE	83 mm	65 mm	Nonsterile	10755-128	Each
50 L (13.2 gal.)	PP	83 mm	65 mm	Nonsterile	10755-130	Each
Accessories						
Stopcock for VWR Carboys					10766-762	Each

Caution

1. Before Autoclaving separate the caps and Spigot/Stopcock threads.
2. Do not stack products, this may cause vacuum and may lead to physical damage/collapse of product.