


AC-1232xx	CUSTOMER DRAWING	Page 1 of 1
Revision 02	Cryo.s™ 1 ml, conical bottom, starfoot base	 greiner bio-one
	Item-No. 123 2xx	



NOT TO SCALE



Cat.-No.	Description	A	B	C	D	E	F
123 261	Cryo.s™ 1 ml, natural	Ø 10.6	Ø 12.5	Ø 12.6	32.0	41.8	17.0
123 263	Cryo.s™ 1 ml, natural	Ø 10.6	Ø 12.5	Ø 12.6	32.0	41.8	17.0
123 277	Cryo.s™ 1 ml, green	Ø 10.6	Ø 12.5	Ø 12.6	32.0	41.8	17.0
123 278	Cryo.s™ 1 ml, yellow	Ø 10.6	Ø 12.5	Ø 12.6	32.0	41.8	17.0
123 279	Cryo.s™ 1 ml, blue	Ø 10.6	Ø 12.5	Ø 12.6	32.0	41.8	17.0
123 280	Cryo.s™ 1 ml, red	Ø 10.6	Ø 12.5	Ø 12.6	32.0	41.8	17.0

All dimensions in mm

Customer drawing subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 01	Date September 12, 2008	Date September 12, 2008	Date September 15, 2008	
Date December 21, 2004	Name F. Falter	Name F. Sonntag	Name F. Sonntag	