

## VWR® MAXIMUM PROTECTION CLEANROOM FACE MASKS

As per the manufacturer, the below product meets the following criteria:

<b>North American Catalog No:</b>	414004-667
<b>Manufacture Lot Number:</b>	200510
<b>Description:</b>	VWR Mask W/Earlps Blue CS500
<b>Date of Manufacture:</b>	May 10, 2020
<b>Country of Origin:</b>	China

### Quality System Compliance

Products are manufactured under the ISO 9001:2015 & ISO13485:2012 standard. Products are Inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP. Inspection records are reviewed and signed off by qualified personnel for product release.

### QC Testing

Representative production samples are collected and inspected in accordance with current applicable product specifications/industry standards-ASTM F2101.

**Bacterial filtration efficiency (BFE):** ASTM F2101 SPEC: >99% Result: >99%

**Pressure Drop (Delta P):** ASTM F2101 SPEC: <3.5 Result: 3.5

### Product Specifications

**Material:** 21TB /25MB /28PP

**BSE/TSE Statement:** The above product contains no ingredients of animal origin and no material derived from, or exposed to animals affected by or under quarantine for Transmitting, Animal Spongiform Encephalopathy (TSE)/Bovine Spongiform Encephalopathy (BSE)

**Non-Pyrogenic:** N/A

**DNase & RNase Free:** N/A

**Cytotoxicity:** N/A

**Latex Free Statement:** The above product is latex free.

**Signed:**



Jamie Ethier  
VP Global Quality  
VWR, Part of Avantor

**Date:** October 12, 2020