

ROTARY EVAPORATOR CONDENSERS

Ace Glass has a comprehensive line of glassware to fit today's popular rotary evaporators. All our rotary evaporator replacement glass is manufactured in the USA, to exacting specifications and standards. Available with and without plastic coating.

Coil Condensers

- Part of glassware sets: G1, G3, G6; A, V
- For all standard and reflux distillations
- Angled or vertical installation
- 35/20 inner ball joint bottom

Dry Ice Condensers



- Part of glassware sets: G5; C, CR
- Condensate trap for low boiling solvents
- Vertical installation
- Cold trap inner body with lid and sealing ring
- 35/20 inner ball joint bottom
- 45/40 bottom joint for use with distribution head

Lower Distribution Head

- Shut-off valve
- Also available with plastic coating

Contact us for...

Custom glass and repair services.



ROTARY EVAPORATOR VAPOR TUBES

These are just a few of the more widely used glassware components to replace your existing glass. All our rotary evaporator replacement glass is manufactured in the USA, to exacting specifications and standards.

Plain Series:

- Plain or with splash guard versions available
- 24/25, 24/40, 29/42 or 45/50 joints

GL Thread Series:

- 45mm or 160mm lengths available
- 29/32, 24/40 or 29/42 joints
- For use with Combi-Clip (*not included*)

Steam Pipe:

- Vapor tube for large scale

Adapter Tube:

- Receiving flask to vapor tube adapter
- 29/32 joint

Also Available...

Contact us for all your rotavap glassware replacement part needs, for both large and bench scale applications.



(From Top to Bottom)
Plain Series, GL Series, Steam Pipe & Adapter Tube



1.800.932.5000 | VWR.COM

Prices and product details are current when published and subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2019 Avantor, Inc. All rights reserved.