

JOINT CLIPS/SLEEVES

Polytetrafluoroethylene (PTFE)

Keck[®] type joint clips snap on and off with ease. These clips have multiple no-pressure uses and can be combined with our PTFE sleeves. PTFE sleeves are for use with glass adapters to *prevent the joint from freezing*, as well as ease of removal for ground glass joints.



Keck[®] Type PTFE Joint Clips

- Won't scratch glass.
- Can be used up to 300°C.
- Can also be used for clamping stopcock plugs and barrels.
- Not for use under pressure.

Joint Size	VWR Cat. No.
14/20	89059-596
24/40	89059-598
29/42	89059-600
45/50	89059-602



PTFE Sleeves

- Precludes contamination from lubricant or leaks.
- Usable to 280°C, and will withstand low to medium vacuum.
- Sold in packs of 3.

0.05mm Wall Thickness

Joint Size	VWR Cat. No.
10/30	10099-430
12/30	89059-894
14/10	89059-872
14/20	89059-582
14/35	89059-896
19/38	89059-900
24/40	89059-876
29/42	89059-906
34/45	89059-880
40/50	10099-432
45/50	89059-940
50/50	89059-942
55/50	89059-944
60/50	10099-434
71/60	10099-436

0.13mm Wall Thickness

7/10	89059-866
10/10	89059-870
10/18	89059-868
14/10	89059-902
14/20	89059-874
24/40	89059-904
29/42	89059-878
34/45	89059-934
45/50	89059-882



PTFE Sleeves w/Gripping Ring


- Sturdy and Reusable.
- Knurled, reinforced gripping ring for easy and safe removal.
- Wall thickness is .4mm (.016").
- Sold in packs of 3.

Joint Size	VWR Cat. No.
14/20	89059-856
24/40	89059-858
29/42	89059-860
34/45	89059-862
45/50	89059-864

STOPPERS

Polytetrafluoroethylene (PTFE)

Wetted surfaces use chemically resistant PTFE. Superior sealing with no need for joint grease. All Ace Glass stoppers come in a wide range of joint sizes, from standard to metric, to fit all of your needs.



PTFE Standard Taper Stopper w/Easy-to-grip Handle

- Easy-to-grip handle helps with removal.

Joint Size	VWR Cat. No.
14/20	89076-968
19/22	89076-970
24/25	89076-972
24/40	89076-974
29/42	89076-974
45/50	89076-978
#13	89076-958
#16	89076-960
#22	89076-962
#27	89076-964
#38	89076-966



Solid PTFE Standard Taper w/Polypropylene Extraction Nut

- Designed to prevent freezing in the joint.
- Extraction nut for easy removal from outer joints.

Joint Size	VWR Cat. No.
10/18	10099-170
10/30	10099-172
12/30	10099-174
14/20	10099-176
14/35	10099-178
19/22	10099-180
19/38	10099-182
24/25	10099-184
24/40	10099-186
29/42	10099-188
34/45	10099-190
45/50	10099-192
55/50	10099-194
#8	10099-196
#9	10099-198
#13	10099-200
#16	10099-202
#19	10099-204
#22	10099-206
#27	10099-208
#32	10099-210
#38	10099-212



PTFE Penny Head Stopper

- Penny head stopper-style handle.

Joint Size	VWR Cat. No.
7/10	89076-942
7/25	89076-944
10/10	89076-946
14/10	89076-948
14/20	89076-950
19/22	89076-952
24/25	89076-954
29/26	89076-956
#8	89076-924
#9	89076-926
#13	89076-928
#16	89076-930
#19	89076-932
#22	89076-934
#27	89076-936
#32	89076-938
#38	89076-940