


Want an easy way to capture your target antigen?



Try the new

Immunoprecipitation Starter Pack

with two media for binding a broad range of antibodies, and easy to follow instructions

- Protein A Sepharose™ 4 Fast Flow and Protein G Sepharose 4 Fast Flow
- 2 x 2 ml – enough for 2 x 80 reactions
- Minimal non-specific adsorption and low ligand leakage
- Step-by-step instructions to guide you through the procedure



Immunoprecipitation is a highly specific and effective technique for capture of antibodies from crude cell lysates. Our new Immunoprecipitation Starter Pack is ideal for starting out in this technique. It contains Protein A Sepharose 4 Fast Flow and Protein G Sepharose 4 Fast Flow, two media with different selectivities which together bind to a wide range of species. Their high specificity and high binding capacity give effective and rapid removal of immune complexes formed between an antigen and its specific antibody in the immunoprecipitation reaction.

Additionally, the media have minimal non-specific adsorption and low ligand leakage.

In combination with techniques like SDS-PAGE and immunoblotting, immunoprecipitation can be used to detect and quantify antigens, determine relative molecular weights, monitor protein turnover and post-translational modifications, and to check for enzyme activity.

Try our new Starter Pack - it comes with comprehensive instructions that not only take you step-by-step through a generic method, but also provide buffer recipes, a troubleshooting guide and references.

Order today

Ordering information

	Pack size	Code No
Immunoprecipitation Starter Pack	2 x 2 ml	17-6002-35
Protein A Sepharose Fast Flow		
Protein G Sepharose Fast Flow		
Data File Immunoprecipitation Starter Pack		18-1141-04

Related products

	Pack size	Code No
Protein A Sepharose Fast Flow	5 ml	17-0974-01
	25 ml	17-0974-04
Protein G Sepharose Fast Flow	5 ml	17-0618-01
	25 ml	17-0618-02
Protein A Sepharose CL-4B	1.5 g	17-0780-01
ExcelGel™ SDS Homogeneous 7.5%	6 gels/pack	80-1260-01
ExcelGel SDS Gradient 8-18%	6 gels/pack	80-1255-53
ExcelGel SDS Buffer Strips	6+6 strips/pack	17-1342-01
Hybond™ ECL™ (6 x 8 cm sheets)	50 sheets	RPN68D
ECL Western Blotting System	1 kit	RPN2108

to order:

Asia Pacific: Tel: +852 2811 8693 Fax: +852 2811 5251 Australasia: Tel: +61 2 9894 5152 Fax: +61 2 9899 7511 Austria: Tel: 01 576 0616 20 Fax: 01 576 0616 27 Belgium: Tel: 0800 73 888 Fax: 03 272 1637
Canada: Tel: 1 800 463 5800 Fax: 1 800 567 1008 Central, East, South East Europe: Tel: +43 1 982 3826 Fax: +43 1 985 8327 Denmark: Tel: 45 16 2400 Fax: 45 16 2424 Finland: Tel: 09 512 3940 Fax: 09 512 1710
France: Tel: 0169 35 67 00 Fax: 0169 35 9677 Germany: Tel: 0761 4903 401 Fax: 0761 4903 405 Italy: Tel: 02 27322 1 Fax: 02 27302 212 Japan: Tel: 81 3 5331 9336 Fax: 81 3 5331 9370 Latin America
Tel: +55 11 3667 5700 Fax: +55 11 3667 87 99 Middle East and Africa: Tel: +30 (1) 96 00 687 Fax: +30 (1) 96 00 693 Netherlands: Tel: 0165 580 410 Fax: 0165 580 401 Norway: Tel: 2318 5800 Fax: 2318 6800
Portugal: Tel: 01 417 7035 Fax: 01 417 7035 Russian Federation: Tel: +7 (095) 232 0250, 956 1137 Fax: +7 (095) 230 6377 South East Asia: Tel: 60 3 724 2080 Fax: 60 3 724 2090 Spain: Tel: 93 594 49 50
Fax: 93 594 49 55 Sweden: Tel: 018 612 1900 Fax: 018 612 1910 Switzerland: Tel: 01 802 81 50 Fax: 01 802 81 51 UK: Tel: 0800 616 928 Fax: 0800 616 927 USA: Tel: +1 800 526 3593 Fax: +1 800 329 3593

Sepharose, ExcelGel, ECL and Hybond are trademarks of Amersham Biosciences or its subsidiaries, Amersham and Amersham Biosciences is a trademark of Amersham plc.

Amersham Biosciences Björkgatan 30, SE-751 84 Uppsala, Sweden. Amersham Biosciences UK Limited Amersham Place, Little Chalfont Buckinghamshire HP7 9NA, England. Amersham Biosciences 800 Centennial Avenue, PO Box 1327, Piscataway, NJ 08855, USA. Amersham Biosciences Europe GmbH Munzinger Strasse 9, D-79111 Freiburg Germany. Amersham Biosciences K.K. Sanken Building, 3-25-1, Hyakunincho, Shinjuku-ku, Tokyo 169-0073, Japan. All goods and services are sold subject to the terms and conditions of sale of the company within the Amersham Biosciences group which supplies them. A copy of these terms and conditions is available on request. © Amersham Biosciences UK Limited 2000 - All rights reserved.