



Quality Certificate

VWR* STIR BARS

This document certifies that the below product meets the following criteria:

North American Catalog No:	58948-952
Description:	VWR BAR PTFE 4X5/8
Material:	Dupont White Teflon #807N-X, Magnet, ALNICO 5
Dimensions:	4 x 5/8" (101.6 X 16mm)
Tolerance	+/- .020

Other:

- Autoclavable at 121 C.
- Materials meet FDA regulations 21CFR177.1550, 21CFR177.1520.
- Teflon used to totally encapsulate the magnetic bar meets US Pharmacopia Class VI.
- Based on the process and materials involved in the manufacturing of this product, this product is considered to be BSE/TSE protein free.
- Based on the information provided by the Manufacturer this product does not contain or come in contact with lead, latex, or other contaminants.
- To the best of the Manufacturer's knowledge this product does not contain or come in contact with any animal origin.
- All items are carefully inspected to conform with ISO 9001:2008, Certificate US94/0161.
- Manufactured in USA.

Signed:

A handwritten signature in blue ink, appearing to read "Ken Crossley", with a stylized flourish at the end.

Ken Crossley
Manager
Quality Assurance

Date : July 16, 2013