

Product Name:	MFS-20-C
Description:	<p>VWR Catalog No: 76285-194</p> <p>VWR Manual Defrost Laboratory Chest Freezer</p> <p>20 Cu. Ft. capacity</p> <p>One solid door</p> <p>Microprocessor temperature controller:</p> <p>Adjustable operating temperature range: -15°C – -25°C</p> <p>Temperature Display & Alarm Module:</p> <p>Sample simulated digital temperature display: .1C resolution, convertible to display temperatures in Celsius or Fahrenheit</p> <p>Audible and visual high/low temperature alarms</p> <p>Sensor error alarm</p> <p>Power failure alarm</p> <p>Resettable 1-point Min/Max History</p> <p>Remote alarm contacts</p> <p>LED Interior light</p> <p>Two removable baskets</p> <p>Swivel casters (2 ¾")</p> <p>Keyed door lock</p> <p>Magnetic door gasket for positive seal</p> <p>White cabinet and lid</p> <p>One probe access port (3/8")</p> <p>High density urethane foam cabinet and door insulation</p> <p>CFC-free refrigerant (R134a), foam insulation and packaging</p> <p>Two year parts and labor warranty plus an additional three year compressor parts warranty</p> <p>Exterior dimensions: 74 3/8"W x 29½"D x 34"H with casters; 31 ¼"H without casters</p> <p>(Height with hinged door open 90° is 61 ¾")</p> <p>Manual defrost</p> <p>115V, 60 Hz, 5 Amps, 1/4 HP</p> <p>Shipping weight: 241 lbs.</p>