

# Avantor® disposable glove guide

01. CLEANROOM GLOVES

02. EXAM OR SAFETY  
GLOVES







# Avantor<sup>®</sup> disposable glove guide

Visit  
[AVANTOR.COM/  
 VWR GLOVES](https://www.avantor.com/vwr-gloves)  
 to search  
 online!

The VWR<sup>®</sup> Collection line of disposable gloves offers a variety of choices depending on your application. Whether your needs are production, laboratory, or industrial focused, we have you covered. Use the guide below to understand the features of our VWR<sup>®</sup> Disposable Gloves and choose the right level of protection from various material, thickness and length options.














## CLEANROOM GLOVES

	Product Line	Packed ISO Class	Static Dissapative	Finger Thickness	Length	Color	Cuff	Texture	Hand Specific/ Ambidextrous	Packaging	Cat. No. Series
<b>Latex</b>											
	VWR <sup>®</sup> CERTICLEAN <sup>®</sup> Class 10 Latex Gloves	ISO Class 4 (FED-STD-209E Class 10/M2.5)	No	8.7	12 in	Natural	Beaded	Full	Ambidextrous	1000/cs	<b>40101-146</b>
	VWR <sup>®</sup> CERTICLEAN <sup>®</sup> Class 100 Latex Gloves	ISO Class 5 (FED-STD-209E Class 100/M3.5)	No	8.7	12 in	Natural	Beaded	Full	Ambidextrous	1000/cs	<b>40101-158</b>
	VWR <sup>®</sup> Cleanroom Latex, Non-Sterile, Class 10	ISO Class 4 (FED-STD-209E Class 10/M2.5)	No	7.3	12 in	White	Beaded	Full	Ambidextrous	1000/cs	<b>76579-674</b>
	VWR <sup>®</sup> Cleanroom Latex, Sterile, Class 10	ISO Class 4 (FED-STD-209E Class 10/M2.5)	No	10.0	12 in	White	Beaded	Fingers	Hand Specific	200pr/cs	<b>76533-186</b>

	Product Line	Packed ISO Class	Static Dissipative	Finger Thickness	Length	Color	Cuff	Texture	Hand Specific/ Ambidextrous	Packaging	Cat. No. Series
<b>Nitrile</b>											
	VWR® CERTICLEAN® Class 10 Nitrile Gloves	ISO Class 4 (FED-STD-209E Class 10/M2.5)	Yes	6.3	12 in	White	Beaded	Fingertips	Ambidextrous	1000/cs	40101-352
	VWR® Maximum Protection Class 10 Nitrile Gloves	ISO Class 4 (FED-STD-209E Class 10/M2.5)	Yes	5.5	12 in	White	Beaded	Fingers	Ambidextrous	1000/cs	89130-846
	VWR® Advanced Protection Class 100 Nitrile Gloves	ISO Class 5 (FED-STD-209E Class 100/M3.5)	Yes	5.5	12 in	White	Beaded	Fingers	Ambidextrous	1000/cs	89130-836
	VWR® Cleanroom Nitrile Non-Sterile Class 10 Non-Tacky	ISO Class 4 (FED-STD-209E Class 10/M2.5)	Yes	6.0	12 in	White	Beaded	Fingers	Ambidextrous	1000/cs	76533-198
	VWR® Cleanroom Nitrile Non-Sterile Class 10 Tacky	ISO Class 4 (FED-STD-209E Class 10/M2.5)	Yes	6.0	12 in	White	Beaded	Fingers	Ambidextrous	1000/cs	76533-210
	VWR® Cleanroom Nitrile, Sterile, Class 10	ISO Class 4 (FED-STD-209E Class 10/M2.5)	Yes	6.5	12 in	White	Beaded	Fingers	Hand Specific	200pr/cs	76533-224
	VWR® Sterile Cleanroom Nitrile Glove, Class 100	ISO Class 5 (FED-STD-209E Class 100/M3.5)	Yes	5.5	12 in	White	Beaded	Fingers	Ambidextrous	1000/cs	76423-448
	VWR® Sterile Class 100 HS Nitrile Gloves	ISO Class 5 (FED-STD-209E Class 100/M3.5)	Yes	5.7	12 in	White	Beaded	Fingers	Hand Specific	200pr/cs	76458-120
	VWR® Cleanroom Nitrile, Non-Sterile, ISO 4, 400mm (16")	ISO Class 4 (FED-STD-209E Class 10/M2.5)	Yes	5.3	16 in	White	Beaded	None	Ambidextrous	1000/cs	76557-792
	VWR® Cleanroom Nitrile, Sterile, ISO 4, 400mm (16")	ISO Class 4 (FED-STD-209E Class 10/M2.5)	Yes	5.3	16 in	White	Beaded	None	Ambidextrous	200pr/cs	76557-804
	VWR® Bagged Examination Nitrile Gloves	ISO Class 7 (FED-STD-209E Class 10,000/M5.5)	Yes	6	12 in	Blue	Beaded		Ambidextrous	1000/cs	77257-062
	VWR® Bagged Examination Nitrile Gloves	ISO Class 7 (FED-STD-209E Class 10,000/M5.5)	Yes	6	12 in	White	Beaded		Ambidextrous	1000/cs	77257-070

## EXAM GLOVES

	Product Line	Feature	Material	Chemo Rated	Finger Thickness	Length	Texture	Cuff	Color	Packaging	Cat. No. Series
	VWR® Powder-Free Nitrile Examination Gloves	Everyday glove, ideal for many applications	Nitrile	Yes	4.7	9 in	Full	Beaded	Blue	100pk/1000/cs	82026-424
	VWR® Premium Thickness Nitrile Examination Gloves	Increased thickness for enhanced protection	Nitrile	Yes	6	9.5 in	Full	Beaded	Blue	100pk/1000/cs	76582-340
	VWR® Black Nitrile Examination Gloves	Exceptional barrier protection, sleek look, fentanyl approved	Nitrile	Yes	5.5	9.5 in	Finger	Beaded	Black	100pk/1000/cs	76307-464
	VWR® Blue Nitrile Examination Glove	Everyday glove, ideal for many applications	Nitrile	Yes	4.7	9 in	Finger	Beaded	Blue	100pk/1000/cs	76518-338
	VWR® Nitrile 200 Powder-Free Examination Gloves	Space-saving and waste-reducing packaging	Nitrile	Yes	3.1	9 in	Finger	Beaded	Blue	200pk/2000cs	89428-748
	VWR® Blue Light Nitrile Examination Glove	Lightweight glove, ideal for many applications	Nitrile	Yes	3	9 in	Finger	Beaded	Blue	100pk/1000/cs	76717-014
	VWR® Ultra Light Black Nitrile Glove	Space-saving and waste-reducing packaging	Nitrile	Yes	2.5	9.5 in	Finger	Beaded	Black	250pk/2500cs	76457-132
	VWR® Aloe Coated Nitrile Examination Gloves	Aloe coated to help with sensitive skin and prevent irritation	Nitrile	Yes	3.5	9.5 in	Finger	Rolled	Blue	200pk/2000cs	76518-176
	VWR® Nitrile Extended Cuff Examination Gloves	Extended cuff for high-risk protection	Nitrile	Yes	6	12 in	Finger	Beaded	Blue	100pk/1000/cs	76582-350
	VWR® Nitrile Maximum Protection Nitrile Examination Gloves	Extra-length cuff for optimal barrier protection	Nitrile	Yes	9	12 in	Full	Beaded	Blue	50pk/500/cs	789142-630
	VWR® Powder-Free Latex Examination Gloves	Natural rubber latex for everyday use, high-barrier protection	Latex	Yes	5.5	9.5 in	Finger	Beaded	Natural	100pk/1000/cs	76319-668