

Control No: 23QC-C-6544  
Date: January 16, 2024

### CERTIFICATE OF CONFORMITY

PRODUCT CODE	LOT NO.	MANUFACTURING DATE	STERILIZATION DATE	EXPIRY DATE
SG3-1838	240110C	2024-01	2024-01	2028-12-31

- I. Visual Inspection : Passed
- II. Dimensional Inspection : Passed
- III. Functional Inspection : Passed
- IV. Sterilization : Passed
- V. Bacterial Endotoxin Test : Passed

The visual, dimensional and functional inspections are checked in accordance with Terumo specifications based on ISO standards.

Sterilization method : Electron Beam Sterilization  
Sterilization dose is determined and audited periodically according to ISO11137.  
Sterilization Assurance Level (SAL) is  $10^{-6}$ .

The bacterial endotoxin limit for medical device is tested in accordance to USP standard.

TERUMO (PHILIPPINES) CORPORATION has established and applied a quality system for manufacturing of medical devices.

Report No.: 176131021-21  
Registration No.: SX 2002969-1  
Furnish requirement: EN ISO 13485:2016  
Valid until: 2025-02-11



Maricris Arcilla  
Quality Control Manager

#### TERUMO (PHILIPPINES) CORPORATION

124 East Main Avenue, Laguna Technopark, Biñan, Laguna  
Tel. No. (049) 541-2111 Fax No. (049) 541-2121