

PRODUCT DATA SHEET

PART NUMBER:414004-267

PART DESCRIPTION:Aspirating Pipet-10ml

MATERIAL

Polystyrene

PRODUCT REQUIREMENTS

Non-pyrogenic, Gamma Irradiated

Dimensions(mm)

Tip:

Length=37.5

Interior Dia.=1.45

Outside Dia.=3.4

PIPET SIZE:10ml

INCREMENTS:N/A

GRADUATIONS:N/A

Mouthpiece:

Length=19.2

Interior Dia.=5.43

Outside Dia.=8

OVERFILL CAPACITY:N/A

COLOR CODE:N/A

Extrusion

Interior Dia.=8

Outside Dia.=9.58

Length=286.3


GENERAL DIMENSIONS

A=9.58

B=8.0

C=3.4

D=343



PACKAGING REQUIREMENTS

Individually wrapped packaging, 50/bag, 200/case

INSPECTION REQUIREMENTS

Representative production samples are collected and inspected in accordance with current applicable product specifications.