

Product: Spin-OUT™ GT-600 **Catalog#:** 786-703, 786-170, 786-171, 786-704 and 786-705

<u>Catalog#</u>	<u>Components</u>	<u>Size</u>
786-703	Spin-OUT™ GT-600, 0.1ml	25/bag
786-170	Spin-OUT™ GT-600, 1ml	10/bag
786-171	Spin-OUT™ GT-600, 3ml	10/bag
786-704	Spin-OUT™ GT-600, 5ml	5/bag
786-705	Spin-OUT™ GT-600, 10ml	5/bag

Material Safety Data Sheet

Section 1: Chemical Identification

Name: Spin-OUT™ GT-600

Catalog #: a component of product Cat. # 786-703, 786-170, 786-171, 786-704 and 786-705 indicated above.

Section 2: Composition/Information on Ingredients Contains: Proprietary resin equivalent to Sephadex G25.

Section 3: Hazards Identification Label precautionary statements. Irritating to eyes, respiratory system, and, skin. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Wear suitable protective clothing.

Section 4: First-Aid Measures

After eye contact: Immediately flush eyes with copious amounts of water for at least 15 minutes.

After skin contact: Immediately wash skin with soap and copious amounts of water.

After Inhalation: None needed for non-volatile items. For items emitting vapor and dust particles: Remove to fresh air. If not breathing give artificial respiration. If breathing is difficult, give oxygen.

After swallowing: Wash out mouth with water provided person is conscious. Call a physician. Wash contaminated clothing before reuse.

Section 5: Fire Fighting Measures

Extinguishing Agents: Water spray, Carbon dioxide, dry chemical powder or appropriate foam.

Special fire fighting procedures: None needed

Unusual fire and explosions hazards: Not known

Section 6: Accidental Release Measures

Human Safety Precautions: Chemical safety goggles, rubber boots and rubber gloves.

Cleaning/Collecting Measures: Sweep up, place in a bag and hold for waste disposal.

Section 7: Handling and Storage Refer to section 8.

Section 8: Exposure Controls/Personal Protection

Personal protective equipment: Chemical safety goggles, compatible chemical-respirator, NIOSH/MSHA-approved respirator, and safety shower and eye bath.

General Protective and Hygienic Measures: Avoid contact with eyes, skin and clothing. Wash thoroughly after handling. Irritant. Keep tightly closed. Store it in a cool dry place.

Section 9: Physical and Chemical Properties

Form: Suspension

Color: Opaque

Odor: None

Section 10: Stability and Reactivity

Thermal Decomposition/Conditions to be Avoided: No decomposition if used according to specifications.

Dangerous Reactions: No dangerous reactions known.

Dangerous Products of Decomposition: No known dangerous decomposition products.

Section 11: Toxicological Information

Primary Irritant Effects: Harmful if inhaled or ingested.

On the skin: Material is irritating to skin and mucous membranes.

On the eye: Causes eye irritation.

Sensitization: No sensitizing effects known.

Section 12: Ecological Information Data not yet available

Section 13: Disposal Considerations

Recommendation: Must not be disposed of together with household garbage. Do not allow product to reach sewage system. Hand over to hazardous waste disposers. Observe all federal, state and local environmental regulations.

Section 14: Transport Information Not dangerous goods. Safe for all transport.

Section 15: Regulatory information Irritating to eyes, respiratory system and skin. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Wear suitable protective clothing.

Section 16: Other Information The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. G-Biosciences shall not be held liable for any damage resulting from handling or from contact with the above product. See reverse side of invoice or packing slip for additional terms and conditions of sale.

G-Biosciences/Genotech St. Louis, MO, USA

Web: www.GBiosciences.com, **email:** technical@GBiosciences.com