

# Material Safety Data Sheet

Pierce® BCA Solid

## 1. Product and company identification

<b>Product name</b>	: Pierce® BCA Solid		
<b>Synonym</b>	: [2,2'-Biquinoline]-4,4'-dicarboxylic acid, disodium salt		
<b>Chemical formula</b>	: C20H12N2O4.2Na		
<b>Supplier</b>	: Thermo Fisher Scientific Pierce Biotechnology P.O. Box 117 Rockford, IL 61105 United States 815.968.0747 or 800.874.3723	<b>Manufacturer</b>	: Thermo Fisher Scientific Pierce Biotechnology P.O. Box 117 Rockford, IL 61105 United States 815.968.0747 or 800.874.3723
<b>Code</b>	: 0023230 0023230B 1879960 1900041 1901041		
<b>MSDS #</b>	: 0228		
<b>Validation date</b>	: 12/20/2010.		
<b>Print date</b>	: 12/20/2010.		
<b>Responsible name</b>	: MSDS (Regulatory Specialist)		
<b>In case of emergency</b>	: <b>CHEMTREC:</b> 800.424.9300 <b>OUTSIDE US:</b> 202.483.7616	<b>Material uses</b>	: Refer to the instruction booklet for proper and intended use. Otherwise, contact supplier for specific applications.
<b>Product type</b>	: Powder.		

## 2. Hazards identification

### Emergency overview

<b>Physical state</b>	: Solid. [Crystalline powder.]
<b>Color</b>	: White to yellowish. [Light]
<b>Odor</b>	: Sweetish.
<b>Signal word</b>	: WARNING!
<b>Hazard statements</b>	: CAUSES RESPIRATORY TRACT, EYE AND SKIN IRRITATION.
<b>Precautionary measures</b>	: Avoid breathing dust. Use only with adequate ventilation. Avoid contact with eyes, skin and clothing. Keep container tightly closed. Wash thoroughly after handling.
<b>OSHA/HCS status</b>	: This material is considered hazardous by the OSHA Hazard Communication Standard (29 CFR 1910.1200).
<b>Routes of entry</b>	: Eye contact. Inhalation. Ingestion.
<b>Potential acute health effects</b>	
<b>Inhalation</b>	: Irritating to respiratory system. Exposure to decomposition products may cause a health hazard. Serious effects may be delayed following exposure.
<b>Ingestion</b>	: No known significant effects or critical hazards.
<b>Skin</b>	: Irritating to skin.
<b>Eyes</b>	: Irritating to eyes.
<b>Potential chronic health effects</b>	

12/20/2010.

1/8

Pierce® BCA Solid

## 2. Hazards identification

<b>Chronic effects</b>	: Repeated or prolonged inhalation of dust may lead to chronic respiratory irritation.
<b>Carcinogenicity</b>	: No known significant effects or critical hazards.
<b>Mutagenicity</b>	: No known significant effects or critical hazards.
<b>Teratogenicity</b>	: No known significant effects or critical hazards.
<b>Developmental effects</b>	: No known significant effects or critical hazards.
<b>Fertility effects</b>	: No known significant effects or critical hazards.
<b>Over-exposure signs/symptoms</b>	
<b>Inhalation</b>	: Adverse symptoms may include the following: respiratory tract irritation coughing
<b>Ingestion</b>	: No specific data.
<b>Skin</b>	: Adverse symptoms may include the following: irritation redness
<b>Eyes</b>	: Adverse symptoms may include the following: pain or irritation watering redness
<b>Medical conditions aggravated by over-exposure</b>	: None known.

See toxicological information (Section 11)

## 3. Composition/information on ingredients

### United States

Name	CAS number	%
[2,2'-Biquinoline]-4,4'-dicarboxylic acid, sodium salt (1:2)	979-88-4	98 - 100

### Canada

Name	CAS number	%
[2,2'-Biquinoline]-4,4'-dicarboxylic acid, sodium salt (1:2)	979-88-4	98 - 100

There are no additional ingredients present which, within the current knowledge of the supplier and in the concentrations applicable, are classified as hazardous to health or the environment and hence require reporting in this section.

## 4. First aid measures

<b>Eye contact</b>	: Check for and remove any contact lenses. Immediately flush eyes with plenty of water for at least 15 minutes, occasionally lifting the upper and lower eyelids. Get medical attention immediately.
<b>Skin contact</b>	: In case of contact, immediately flush skin with plenty of water for at least 15 minutes while removing contaminated clothing and shoes. Wash clothing before reuse. Clean shoes thoroughly before reuse. Get medical attention immediately.
<b>Inhalation</b>	: Move exposed person to fresh air. If not breathing, if breathing is irregular or if respiratory arrest occurs, provide artificial respiration or oxygen by trained personnel. Loosen tight clothing such as a collar, tie, belt or waistband. Get medical attention immediately.
<b>Ingestion</b>	: Wash out mouth with water. Do not induce vomiting unless directed to do so by medical personnel. Never give anything by mouth to an unconscious person. Get medical attention immediately.

12/20/2010.

2/8

Pierce® BCA Solid

## 4. First aid measures

**Protection of first-aiders** : No action shall be taken involving any personal risk or without suitable training. If it is suspected that fumes are still present, the rescuer should wear an appropriate mask or self-contained breathing apparatus. It may be dangerous to the person providing aid to give mouth-to-mouth resuscitation.

**Notes to physician** : In case of inhalation of decomposition products in a fire, symptoms may be delayed. The exposed person may need to be kept under medical surveillance for 48 hours.

## 5. Fire-fighting measures

**Flammability of the product** : Fine dust clouds may form explosive mixtures with air.

### Extinguishing media

**Suitable** : Use dry chemical powder.

**Not suitable** : Do not use water jet.

**Special exposure hazards** : Promptly isolate the scene by removing all persons from the vicinity of the incident if there is a fire. No action shall be taken involving any personal risk or without suitable training. Move containers from fire area if this can be done without risk. Use water spray to keep fire-exposed containers cool.

**Hazardous thermal decomposition products** : Decomposition products may include the following materials:  
carbon dioxide  
carbon monoxide  
nitrogen oxides  
metal oxide/oxides

**Special protective equipment for fire-fighters** : Fire-fighters should wear appropriate protective equipment and self-contained breathing apparatus (SCBA) with a full face-piece operated in positive pressure mode.

## 6. Accidental release measures

**Personal precautions** : No action shall be taken involving any personal risk or without suitable training. Evacuate surrounding areas. Keep unnecessary and unprotected personnel from entering. Do not touch or walk through spilled material. Shut off all ignition sources. No flares, smoking or flames in hazard area. Avoid breathing dust. Provide adequate ventilation. Wear appropriate respirator when ventilation is inadequate. Put on appropriate personal protective equipment (see Section 8).

**Environmental precautions** : Avoid dispersal of spilled material and runoff and contact with soil, waterways, drains and sewers. Inform the relevant authorities if the product has caused environmental pollution (sewers, waterways, soil or air).

### Methods for cleaning up

**Small spill** : Move containers from spill area. Vacuum or sweep up material and place in a designated, labeled waste container. Use spark-proof tools and explosion-proof equipment. Dispose of via a licensed waste disposal contractor.

**Large spill** : Move containers from spill area. Approach release from upwind. Prevent entry into sewers, water courses, basements or confined areas. Vacuum or sweep up material and place in a designated, labeled waste container. Avoid creating dusty conditions and prevent wind dispersal. Use spark-proof tools and explosion-proof equipment. Dispose of via a licensed waste disposal contractor. Note: see section 1 for emergency contact information and section 13 for waste disposal.

## 7. Handling and storage

**Handling** : Put on appropriate personal protective equipment (see Section 8). Eating, drinking and smoking should be prohibited in areas where this material is handled, stored and processed. Workers should wash hands and face before eating, drinking and smoking. Remove contaminated clothing and protective equipment before entering eating areas. Do not ingest. Avoid contact with eyes, skin and clothing. Avoid breathing dust. Avoid the creation of dust when handling and avoid all possible sources of ignition (spark or flame). Prevent dust accumulation. Use only with adequate ventilation. Wear appropriate respirator when ventilation is inadequate. Keep in the original container or

12/20/2010.

3/8

Life Science Research PO Box 117 Rockford, IL (815) 968-0747 www.thermo.com  
Pierce Biotechnology Inc. 3747 N. Meridian Road 61105 (815) 968-7316 Fax

Pierce® BCA Solid

## 7. Handling and storage

an approved alternative made from a compatible material, kept tightly closed when not in use. Electrical equipment and lighting should be protected to appropriate standards to prevent dust coming into contact with hot surfaces, sparks or other ignition sources. Empty containers retain product residue and can be hazardous. Do not reuse container.

**Storage** : Store between the following temperatures: 20 to 25°C (68 to 77°F). Store in accordance with local regulations. Store in a segregated and approved area. Store in original container protected from direct sunlight in a dry, cool and well-ventilated area, away from incompatible materials (see section 10) and food and drink. Eliminate all ignition sources. Separate from oxidizing materials. Keep container tightly closed and sealed until ready for use. Containers that have been opened must be carefully resealed and kept upright to prevent leakage. Do not store in unlabeled containers. Use appropriate containment to avoid environmental contamination.

## 8. Exposure controls/personal protection

### Canada

#### Occupational exposure limits

No exposure limit value known.

#### Consult local authorities for acceptable exposure limits.

**Recommended monitoring procedures** : If this product contains ingredients with exposure limits, personal, workplace atmosphere or biological monitoring may be required to determine the effectiveness of the ventilation or other control measures and/or the necessity to use respiratory protective equipment.

**Engineering measures** : Use only with adequate ventilation. If user operations generate dust, fumes, gas, vapor or mist, use process enclosures, local exhaust ventilation or other engineering controls to keep worker exposure to airborne contaminants below any recommended or statutory limits. The engineering controls also need to keep gas, vapor or dust concentrations below any lower explosive limits. Use explosion-proof ventilation equipment.

**Hygiene measures** : Wash hands, forearms and face thoroughly after handling chemical products, before eating, smoking and using the lavatory and at the end of the working period. Appropriate techniques should be used to remove potentially contaminated clothing. Wash contaminated clothing before reusing. Ensure that eyewash stations and safety showers are close to the workstation location.

### Personal protection

#### Respiratory

: Use a properly fitted, air-purifying or air-fed respirator complying with an approved standard if a risk assessment indicates this is necessary. Respirator selection must be based on known or anticipated exposure levels, the hazards of the product and the safe working limits of the selected respirator.

#### Hands

: Chemical-resistant, impervious gloves complying with an approved standard should be worn at all times when handling chemical products if a risk assessment indicates this is necessary.

#### Eyes

: Safety eyewear complying with an approved standard should be used when a risk assessment indicates this is necessary to avoid exposure to liquid splashes, mists or dusts. If operating conditions cause high dust concentrations to be produced, use dust goggles.

#### Skin

: Personal protective equipment for the body should be selected based on the task being performed and the risks involved and should be approved by a specialist before handling this product.

#### Environmental exposure controls

: Emissions from ventilation or work process equipment should be checked to ensure they comply with the requirements of environmental protection legislation. In some cases, fume scrubbers, filters or engineering modifications to the process equipment will be necessary to reduce emissions to acceptable levels.

12/20/2010.

4/8

Life Science Research PO Box 117 Rockford, IL (815) 968-0747 www.thermo.com  
Pierce Biotechnology Inc. 3747 N. Meridian Road 61105 (815) 968-7316 Fax

Pierce® BCA Solid

## 9. Physical and chemical properties

Physical state	: Solid. [Crystalline powder.]
Color	: White to yellowish. [Light]
Odor	: Sweetish.
Molecular weight	: 388.29 g/mole
Molecular formula	: C20H12N2O4.2Na
pH	: 6.5 to 8.5 [Conc. (% w/w): 1%]
Melting/freezing point	: >290°C (>554°F)
Solubility	: Easily soluble in the following materials: cold water and hot water. Soluble in the following materials: methanol.

## 10. Stability and reactivity

Chemical stability	: The product is stable.
Conditions to avoid	: Avoid the creation of dust when handling and avoid all possible sources of ignition (spark or flame). Take precautionary measures against electrostatic discharges. To avoid fire or explosion, dissipate static electricity during transfer by grounding and bonding containers and equipment before transferring material. Prevent dust accumulation.
Incompatible materials	: Reactive or incompatible with the following materials: oxidizing materials
Hazardous decomposition products	: Under normal conditions of storage and use, hazardous decomposition products should not be produced.
Possibility of hazardous reactions	: Under normal conditions of storage and use, hazardous reactions will not occur.

## 11. Toxicological information

### United States

#### Acute toxicity

Conclusion/Summary	: To the best of our knowledge, the toxicological properties of this substance have not been thoroughly investigated.
--------------------	---

#### Chronic toxicity

Conclusion/Summary	: Not available.
--------------------	------------------

#### Irritation/Corrosion

Conclusion/Summary	: Not available.
--------------------	------------------

#### Sensitizer

Conclusion/Summary	: Not available.
--------------------	------------------

#### Carcinogenicity

Conclusion/Summary	: Not available.
--------------------	------------------

#### Mutagenicity

Conclusion/Summary	: Not available.
--------------------	------------------

#### Teratogenicity

Conclusion/Summary	: Not available.
--------------------	------------------

#### Reproductive toxicity

Conclusion/Summary	: Not available.
--------------------	------------------

### Canada

#### Acute toxicity

Conclusion/Summary	: To the best of our knowledge, the toxicological properties of this substance have not been thoroughly investigated.
--------------------	---

#### Chronic toxicity

Conclusion/Summary	: Not available.
--------------------	------------------

12/20/2010.

Life Science Research  
Pierce Biotechnology Inc. PO Box 117  
3747 N. Meridian Road Rockford, IL 61105 (815) 968-0747  
(815) 968-7316 Fax www.thermo.com

5/8

Pierce® BCA Solid

## 11. Toxicological information

### Irritation/Corrosion

Conclusion/Summary	: Not available.
--------------------	------------------

### Sensitizer

Conclusion/Summary	: Not available.
--------------------	------------------

### Carcinogenicity

Conclusion/Summary	: Not available.
--------------------	------------------

### Mutagenicity

Conclusion/Summary	: Not available.
--------------------	------------------

### Teratogenicity

Conclusion/Summary	: Not available.
--------------------	------------------

### Reproductive toxicity

Conclusion/Summary	: Not available.
--------------------	------------------

## 12. Ecological information

Ecotoxicity	: No known significant effects or critical hazards.
-------------	---

### United States

#### Aquatic ecotoxicity

Conclusion/Summary	: Not available.
--------------------	------------------

#### Persistence/degradability

Conclusion/Summary	: Not available.
--------------------	------------------

### Canada

#### Aquatic ecotoxicity

Conclusion/Summary	: Not available.
--------------------	------------------

#### Persistence/degradability

Conclusion/Summary	: Not available.
--------------------	------------------

#### Other adverse effects

	: No known significant effects or critical hazards.
--	---

## 13. Disposal considerations

Waste disposal	: The generation of waste should be avoided or minimized wherever possible. Significant quantities of waste product residues should not be disposed of via the foul sewer but processed in a suitable effluent treatment plant. Dispose of surplus and non-recyclable products via a licensed waste disposal contractor. Disposal of this product, solutions and any by-products should at all times comply with the requirements of environmental protection and waste disposal legislation and any regional local authority requirements. Waste packaging should be recycled. Incineration or landfill should only be considered when recycling is not feasible. This material and its container must be disposed of in a safe way. Care should be taken when handling emptied containers that have not been cleaned or rinsed out. Empty containers or liners may retain some product residues. Avoid dispersal of spilled material and runoff and contact with soil, waterways, drains and sewers.
----------------	--

Disposal should be in accordance with applicable regional, national and local laws and regulations.

Refer to Section 7: HANDLING AND STORAGE and Section 8: EXPOSURE CONTROLS/PERSONAL PROTECTION for additional handling information and protection of employees.

12/20/2010.

Life Science Research  
Pierce Biotechnology Inc. PO Box 117  
3747 N. Meridian Road Rockford, IL 61105 (815) 968-0747  
(815) 968-7316 Fax www.thermo.com

6/8

Pierce® BCA Solid

## 14. Transport information

Regulatory information	UN number	Proper shipping name	Classes	PG*
DOT Classification	Not regulated.	-	-	-
IATA-DGR Class	Not regulated.	-	-	-

PG\* : Packing group

## 15. Regulatory information

### United States

- HCS Classification** : Irritating material
- U.S. Federal regulations** : **TSCA 8(a) IUR Exempt/Partial exemption**: Not determined  
**United States inventory (TSCA 8b)**: This material is listed or exempted.  
**SARA 302/304/311/312 extremely hazardous substances**: No products were found.  
**SARA 302/304 emergency planning and notification**: No products were found.  
**SARA 302/304/311/312 hazardous chemicals**: [2,2'-Biquinoline]-4,4'-dicarboxylic acid, sodium salt (1:2)  
**SARA 311/312 MSDS distribution - chemical inventory - hazard identification**: [2,2'-Biquinoline]-4,4'-dicarboxylic acid, sodium salt (1:2): Immediate (acute) health hazard

**Clean Air Act Section 112(b) Hazardous Air Pollutants (HAPs)** : Not listed

**Clean Air Act Section 602 Class I Substances** : Not listed

**Clean Air Act Section 602 Class II Substances** : Not listed

**DEA List I Chemicals (Precursor Chemicals)** : Not listed

**DEA List II Chemicals (Essential Chemicals)** : Not listed

### State regulations

**Massachusetts** : This material is not listed.

**New York** : This material is not listed.

**New Jersey** : This material is not listed.

**Pennsylvania** : This material is not listed.

**United States inventory (TSCA 8b)** : This material is listed or exempted.

### Canada

**WHMIS (Canada)** : Class D-2B: Material causing other toxic effects (Toxic).

### Canadian lists

**Canadian NPRI** : This material is not listed.

**CEPA Toxic substances** : This material is not listed.

**Canada inventory** : This material is listed or exempted.

This product has been classified in accordance with the hazard criteria of the Controlled Products Regulations and the MSDS contains all the information required by the Controlled Products Regulations.

### International regulations

12/20/2010.

7/8

Life Science Research  
Pierce Biotechnology Inc. PO Box 117  
3747 N. Meridian Road Rockford, IL 61105 (815) 968-0747  
(815) 968-7316 Fax www.thermo.com

Pierce® BCA Solid

## 15. Regulatory information

- International lists** : **Australia inventory (AICS)**: Not determined.  
**China inventory (IECSC)**: Not determined.  
**Japan inventory**: Not determined.  
**Korea inventory**: This material is listed or exempted.  
**New Zealand Inventory of Chemicals (NZIoC)**: This material is listed or exempted.  
**Philippines inventory (PICCS)**: Not determined.

## 16. Other information

**Label requirements** : CAUSES RESPIRATORY TRACT, EYE AND SKIN IRRITATION.

**Hazardous Material Information System (U.S.A.)** :

Health	2
Flammability	0
Physical hazards	0

The customer is responsible for determining the PPE code for this material.

**National Fire Protection Association (U.S.A.)** :



**Date of printing** : 12/20/2010.

**Date of issue** : 12/20/2010.

**Date of previous issue** : No previous validation.

**Version** : 1

**Prepared by** : MSDS (Regulatory Specialist)

✓ Indicates information that has changed from previously issued version.

### Notice to reader

To the best of our knowledge, the information contained herein is accurate. However, neither the above-named supplier, nor any of its subsidiaries, assumes any liability whatsoever for the accuracy or completeness of the information contained herein.

Final determination of suitability of any material is the sole responsibility of the user. All materials may present unknown hazards and should be used with caution. Although certain hazards are described herein, we cannot guarantee that these are the only hazards that exist.

12/20/2010.

8/8

Life Science Research  
Pierce Biotechnology Inc. PO Box 117  
3747 N. Meridian Road Rockford, IL 61105 (815) 968-0747  
(815) 968-7316 Fax www.thermo.com