

PDS No. 73118x	PRODUCT DATA SHEET				Page 1 of 1
Revision 02	Disposable Inoculation Needle, PS				 greiner bio-one
	Greiner Item-No. 73118x				
Valid for Item-No.:	731180 (sterile)	731181 (sterile)	731185 (sterile)		

1.	Description / Specification	
1.1	Description	PS disposable inoculation needle
1.2	Dimensions	Length: 200 mm
1.3	Volume	-
1.4	Material / Resin	PS (Polystyrene), free from heavy metal
1.5	Colour	Yellow
1.6	Sterilization	SAL 10 ⁻³
1.7	Quality Control	- <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	Hydrophobic
2.2	Temperature range	-20°C to +60°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: www.gbo.com/bioscience →Products →Literature →Technical Information →Chemical Resistance of Resins
2.6	Shelf life	5 years after month of production
2.7	Other Information	-

3.	Packaging	731180	731181	731185
3.1	Pieces / Bag	50	1	10
3.2	Pieces / Box	2000	600	3000
3.3	Lot-No.	PPYWWX (working place, year, week, index)	PPYWWX (working place, year, week, index)	PPYWWX (working place, year, week, index)
3.4	Other Information	-	-	

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 01	Date 14 November 2006	Date 15 November 2006	Date 20 November 2006	
Date 15.03.05	Name S. Kaelberer	Name Dr. L. Marchetti	Name A. Schulz	