

KIMTECH™


EXAM GLOVES QUICK REFERENCE GUIDE



Kimtech™ Exam Gloves provide protection from known and unknown risks while delivering comfort, dexterity and the assurance to do your best science. They exceed industry standards and work hand-in-hand with other head-to-toe PPE solutions from Kimberly-Clark Professional. After use, give your gloves new life through The RightCycle™ Program, the first large-scale recycling program for non-hazardous PPE waste.

▶ PRODUCT FEATURES

- FDA cleared as examination gloves
- Powder-free
- Nitrile gloves not made with natural rubber latex
- Beaded cuff for donning
- Ambidextrous

Description	VWR Code	Color	Size	Length	Chemo Tested	Food Handling	Thickness (Fingertip)	Texture Area	Gloves/Box	Total Gloves/Case
Polaris™ Nitrile-Xtra™  	77308-526, 77308-528, 77308-530, 77308-532, 77308-534		XS-XL	12"	▲*	▲	7.5 mil	Fingertips	50/10	500 each
Polaris™ Nitrile  	77308-536, 77308-538, 77308-540, 77308-542, 77308-544		XS-XL	9.5"	▲*	▲	5.9 mil	Fingertips	100/10 (XL = 90/10)	1,000 each (XL=900)
Purple Nitrile-Xtra™ 	32933-968, 32933-970, 32933-972, 32933-974, 32933-976		XS-XL	12"	▲		5.9 mil	Fingertips	50/10	500 each
Purple Nitrile™ 	32934-076, 32934-078, 32934-080, 32934-082, 32934-084		XS-XL	9.5"	▲		5.9 mil	Fingertips	100/10 (XL=90/10)	1,000 each (XL=900)
Vista™ Nitrile 	76520-546, 76520-548, 76520-550, 76520-552, 76520-554		XS-XL	9.5"	▲	▲	4.3 mil	Fingertips	200/10 (XL=180/10)	2,000 each (XL=1,800)
Sterling Nitrile-Xtra™ 	94023-906, 94023-908, 94023-910, 94023-912, 94023-914		XS-XL	12"	▲	▲	3.5 mil	Fingertips	100/10	1,000 each
Sterling™ Nitrile 	89140-816, 89140-818, 89140-820, 89140-822, 89140-824		XS-XL	9.5"	▲	▲	3.5 mil	Fingertips	200/10 (XL=170/10)	2,000 each (XL=1,700)
Element™ Nitrile 	76520-536, 76520-538, 76520-540, 76520-542, 76520-544		XS-XL	9.0"		▲	3.2 mil	Fingertips	250/10 (XL=230/10)	2,500 each (XL=2,300)



Apex™ Make your laboratory exceptional

APEX is a unique employee training program that helps facilitate the changeover to the Kimtech™ Brand of products. By effectively engaging employees, it ensures protection of the scientific process and enhances efficiency and performance.



United States
Ergonomics Certified



Acceptable quality level.
Test Method: ASTM D5151



All Kimtech™ Exam Gloves are recyclable
through The RightCycle™ Program

* "Opioid" fentanyl and gastric acid tested.
Test Method: ASTM D6978-05

©/™ Trademarks of Kimberly-Clark Worldwide, Inc. or its affiliates. Marques de commerce de Kimberly-Clark Worldwide, Inc. ou de ses sociétés affiliées. © KCWW. T23H53-01 VWR 2/24



VWR.COM

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada unless otherwise noted, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit [vwr.com](https://www.vwr.com) to view our privacy policy, trademark owners, and additional disclaimers. © 2024 Avantor, Inc. All rights reserved.