



CHI Scientific
Accelerating Primary Cell Research



The fibroblast growth control system, FibrOut® is a combination of several biochemical compounds and reagents, which prevent overgrowth of fibroblastic and contaminated cells and increase targeted cell growth during primary cell culture. Each fibroblast growth inhibitory system is customized to be tissue-specific and cell-specific, and proven to work most effectively with our corresponding PrimaCell™ system, the complete tissue-specific and cell-type-specific primary cell culture systems.

Exclusively available from CHI Scientific through VWR!

Each FibrOut® system may contain one or more of the following chemicals:

Trypsin	Collagenase	D-valine
Cis-OH-proline	Ethylmercurithiosalicylate	Phenobarbitone
Formulated serum substitutes	Antimesodermal antibodies	

	Liver	Intestines	Kidney	Adipose	Neural	Lung	Reproductive Tissues	Thymus	Mammary	Pituitary	Fat	Blood Vessels	Pancreas	Brain	Skin	Colon	Bone	Endothelial Tissues	Heart Muscle	Cartilage	Adrenal
Human	+	-	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+
Avian	+	-	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+
Mouse	+	-	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+
Rat	+	-	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+

Selected FibrOut® (more than 80 selections available on vwr.com)

VWR Part No.	Description
10787-032	Human FibrOut® 1, for adipose, fat, eye, lung, mammary
10786-022	Human FibrOut® 9, for tumors
10785-998	Mouse FibrOut® 3, for bone, heart, muscle
10786-058	Rat FibrOut® 7, for reproductive tissues



1.800.932.5000 | vwr.com

Prices and product details are current when published; subject to change without notice. Certain products may be limited by federal, state, provincial, or local regulations. VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US, void where prohibited by law or company policy, while supplies last. VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks, of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). Visit vwr.com to view our privacy policy, trademark owners and additional disclaimers. ©2016 VWR International, LLC. All rights reserved.