

Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® Grade 415 Filter Paper, Qualitative, Crepe
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2021-12-03

VWR Catalogue Number	Diameter [cm]
28320-019	5.5
28320-020	7.5
28320-041	9.0
28320-085	11.0
28320-100	12.5
28320-121	15.0
28320-143	18.5
28320-187	20.5
28320-201	25.0
28320-223	33.0
28297-846	38.5

Quality System Compliance

Products are manufactured under the [ISO 9001:2015 and ISO 13485: 2016](#) standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Test	Range	Result
Grammage, gsm	68-84	Pass
Thickness, µm@53Kpa	195-325	Pass
Porosity, s/300ml/0.1sq.in	40-80	Pass

Product Specifications

Material	Grade 415 Filter Paper. Particle Retention: 25 µm Porosity: Coarse
Sterilization	Not applicable.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	Not applicable.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	Not applicable.
REACH Statement	Not applicable.
Storage Conditions	Store at 22-24°C and Humidity 48-52%.
Shelf Life	5 years.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.