



Certificate of Analysis

Product:
Ficoll-Paque™ PLUS
Code Number:

17144003

Lot No: 10298772

Test/Characteristic:	Limit:	Result:
Sterility, F₀-value minutes	min. 15	Passed
Density g/ml, 20°C	1.076 - 1.078	1.077
pH	5.5 - 7.5	6.7
UV / Visible spectra	to pass test	Passed
Sodium Diatrizoate g/dl	8.1 - 9.9	9.2
Colour; A = 405 nm AU	max. 0.20	0.11
Endotoxin activity EU/ml	less than 0.12	<0.12

Note: Sterility is ensured by in-process control.

Manufactured in compliance with our ISO 9001 certified quality management system.

No animal derived material is used in raw materials or in the manufacture of this product.

Approval Date (Year-Month-Day): 2020-10-05

Expiry Date (Year-Month): 2023-09

Manufacturing Date (Year-Month): 2020-09

Cytiva and the Drop logo are trademarks of Global Life Sciences
IP Holdco LLC or an affiliate

Cytiva Sweden AB
Björkgatan 30
SE-751 84 Uppsala, Sweden
T + 46 (0)18 612 00 00
Cytiva.com
Reg.No. SE 55 61 08 1919 01

Tests and limits according to AS 45-6006-81 Ed. AE
Issued (Year-Month-day): 2020-10-05 by Quality Assurance
Medin, Britta

This document has been electronically produced and is valid without a signature.

DOC1103901 / 7