



BioPharma Protection Solutions

Pharmaceutical

Applications:

- Loading and unloading of process equipment
- Transferring liquids and solids
- Loading wetcakes and unloading dry products
- Packaging and labeling

User needs:

- Chemical protection
- Dexterity and grip
- Abrasion and puncture resistance
- Sterile and/or clean

Biotechnology

Applications:

- Transferring liquids
- Loading centrifuges and chromatography columns
- Sterilizing vials and ampules
- Aseptic filling and labeling

User needs:


- Wet grip
- Sterile and/or clean
- Dexterity
- Double donning

BioClean™ N-Plus (BNPS)



Benefits & Features:

- 400mm/ 16" length for extra protection
- Hand specific
- Sterile nitrile
- Approved for use with chemotherapy drugs (ASTM D 6978)



Size	VWR Cat. No.
6	75834-048
6.5	75834-050
7	75834-052
7.5	75834-054
8	75834-056
8.5	75834-058
9	75834-060
10	75834-118

BioClean-D™ Drop-down Garment (S-BDSH)



Benefits & Features:

- Innovative quick & easy donning design
- Lightweight for increased user comfort/ wearability
- Sterile & Class 10 (ISO 4) compatible
- Low-linting antistatic CleanTough material



Size	VWR Cat. No.
XS	75845-800
S	75845-796
M	75845-794
L	75845-792
XL	75845-798
2XL	75845-802
3XL	75845-804
4XL	75845-806
5XL	76271-424
6XL	76271-426

BioClean™ Clearview Autoclavable Goggles (BCAP)



Benefits & Features:

- Indirect ventilation system to reduce risk of contamination into the controlled environment
- Toughened, anti-scratch, anti-fog lens for clear vision
- Anti-fog performance maintained for up to 25 autoclave cycles (under laboratory conditions)
- Panoramic lens for increased field of vision



VWR Cat. No.
75845-813

