

Boreal2 AP Stereo	
	470014-538
Eyepiece	Paired 10x, 20mm
Head	Fixed 45° inclined head
Magnification	2X 4X
Stage	Frosted glass and reversible black & white
Working Distance	89mm
Focusing	Coarse focus with tension adjustment and slip clutch
Top Illumination	Corded LED
Bottom Illumination	Corded LED
Replacement Bulb	(T)Ward's Item #470330-980 (B) Ward's #470330-978
Product Includes	Removable soft eyecups and vinyl dust cover
Dimensions	12" (H) x 7" (W) x 9" (D)
Warranty	Boreal Lifetime (includes shipping, excludes bulb & fuse)