



Quality Certificate

VWR® POWDER FUNNELS, PARALLEL STEMS, POLYPROPYLENE

As per the manufacturer, the product identified below meets the following criteria:

North American Catalog No:	30252-957
Lot Number:	1288991
Description:	VWR Funnel Pwdr 150mm Od Pk4.
Manufactured Date:	June 12, 2020

The product is inspected and released according to the requirements of product and quality specifications as maintained in the manufacturer's quality management system which is compliant to ISO 13485:2016 in Monterrey Mexico.

Description	Common Name	DMF #	Cytotoxicity	USP Class VI	FDA Compliance - 21 CFR
FUN, POWDER, 150MM, PP	COMPONENT PART				
RESIN, PP, COPOLYMER, INJ	POLYPROPYLENE COPOLYMER	7478	N/A	PASSED	177.1520(a)(3)(i), (c)3.1a, except for cooking, (use conditions C-H)

If N/A appears in any of the columns above, it means the information is not available. Any item listed as "COMPONENT PART" will show blank in the DMF#, Cytotoxicity, USP Class VI, and FDA Compliance Information columns.

If the word "PASSED" appears in the USP Class VI column next to the resin listing, this material has passed USP Class VI requirements, latest Volume, as part of the manufacturer's initial test approval protocol.

If the word "PASSED" appears in the Cytotoxicity column next to the resin listing, this material was tested and shown to be non-cytotoxic as part of the manufacturer's initial test approval protocol, using either mouse fibroblast L929 cells or the more sensitive human diploid lung cell lines WI-38 or MRC-5.

Signed:

Jamie Ethier
VP Global Quality
VWR, Part of Avantor

Date: October 12, 2020