



Interior

Front View



Side View



Cat# 10819-768/FRCV-30

115VAC/60Hz/1PH
NEMA-5-15P

Compartment	Cu.Ft.	Defrost	Ft Door	Int Door	Shelves	Interior Dimensions			Refrigerant	H.P.	Amps	Weight
						W"	D"	H"				
Freezer	30	Manual	1 Solid	1	0	29	22 1/4	12	R404a	1/3	6	336 lbs.
Refrigerator				0	3 Adj	32	22 1/2	50 1/4				

Note: Interior dimensions are estimates and subject to change without notice