

Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® Premium Plus Slide Boxes, 100-Place
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2022-02-23

VWR Catalogue Number	Colour
76265-652	Blue
76265-656	Green
76265-654	Orange
76265-622	Purple
76265-620	Red
76265-640	Assorted

Quality System Compliance

Products are manufactured in accordance with the manufacture's quality system. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Product Specifications

Material	Base injection-molded of acrylonitrile-butadiene-styrene homopolymer with non-toxic blue pigment, which are inert or minimally reactive with most common laboratory reagents. Lid injection-molded of polycarbonate with no pigment, which is inert or minimally reactive with most common laboratory reagents. Inspected in the molding facility for cracks, incomplete or distorted shape, or other damage or defects which would compromise proper function.
Sterilization	Not applicable.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.

BSE/TSE	No use of any raw material produced from, or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	The product is Latex free.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Please see RoHS information below.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr ⁶⁺), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP), DEHP (di(2- ethylhexyl) phthalate)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2021-07-08 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
Storage Conditions	Store down to -80°C. Must use tubes capable of being stored at the required temperature of the user.
Shelf Life	Not applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.