


Component Schematic

96-well Microplate (PP) Deep Well (Round) Rimless



L	128.00 ± 0.10	Hw	39.00 ± 0.05
W	86.00 ± 0.10	P	9.00 ± 0.05
H	41.60 ± 0.05	D	7.05 ± 0.05
Hs	3.35 ± 0.10	C	12
X1	14.40 ± 0.05	R	8
Xn	113.40 ± 0.05		
Y1	11.50 ± 0.05		
Yn	74.50 ± 0.05		

Volume 1,300μl

Well form
Material

U Bottom
PP