



SAFEGUARD YOUR SAMPLE TRACKING AND STORAGE

CRYO.S WITH DATAMATRIX AND LINEAR BARCODING

With the option for barcoding, Greiner Bio-One adds a beneficial feature to its well-established range of Cryo.s.

Barcoded Cryo.s carry a 1D barcode on tube side and a 2D code ECC at starfoot bottom.

Qualified for use for use within automated sample processing and storage systems. Additionally, each tube features a human-readable representation of code contents, enabling user-friendly readout of tube ID without barcode readers.

Barcodes are applied using state-of-the-art laser labelling technologies, resulting in extremely robust and temperature-resistant labels.

Unique ID, pre-coded tubes available off-the-shelf, with customized barcoding services also available.

KEY FACTS :

- / Laser-written 1D barcode type 128 and human readable code on tube side
- / Starfoot base tubes contain 2D Datamatrix code ECC 200 with 14x14 dots
- / Reed-Solomon error correction for high redundancy of encrypted information
- / Strong resistance to mechanical and chemical damage
- / High data capacity
- / Off-the-shelf pre-coded tubes, as well as customized barcoding, with unique, sequential barcodes

PRODUCT OVERVIEW

Ordering information:-

VWR PN	Product Description	Pack QTY	Case QTY
89176-674	2 mL Cryo.s w/Datamatrix (generic 2D sequence), PP, Sterile, Internal Thread Natural Cap, ID Field, Starfoot Base, Round Bottom, 150 White Cap Inserts per CS	100	500
89176-676	1 mL Cryo.s w/Datamatrix (generic 2D sequence), PP, Sterile, Internal Thread Natural Cap, ID Field, Starfoot Base, Conical Bottom, 150 White Cap Inserts per CS	100	500
89176-678	2 mL Cryo.s w/Datamatrix (generic 2D sequence), PP, Sterile, External Thread Natural Cap, ID Field, Starfoot Base, Round Bottom, 150 White Cap Inserts per CS	100	500
89176-680	4 mL w/Datamatrix (generic 2D sequence), PP, Sterile, External Thread Natural Cap, ID Field, Starfoot Base, Round Bottom, 100 White Cap Inserts per CS	50	300



VWR.COM

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada unless otherwise noted, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit [vwr.com](https://www.vwr.com) to view our privacy policy, trademark owners, and additional disclaimers. © 2023 Avantor, Inc. All rights reserved.