

# Fall Discounts on VWR® ULT Freezers Eco Premium

**SAVE 15% ON VWR® ULT FREEZERS ECO PREMIUM THROUGH DECEMBER 31, 2022!**

The VWR® ULT Freezers Eco Premium are designed for safe and reliable storage of sensitive samples. They feature ergonomic control and low energy consumption.

- Maximum 2" box capacity: 352/528
- Energy star certified
- Temperatures between –86 and –50 °C
- Long warm-up time (Warm up time from –86 and –50 °C: 360 min)
- Highly efficient VIP insulation and long-lasting rust-proof sample chamber
- Green refrigerants
- USB and Ethernet interface for data transfer

The VWR® ULT Freezers Eco Premium have a stainless steel inner chamber, flexible rack configuration and are further characterized by low energy consumption and electricity costs. They have individually adaptable chamber compartments and samples remain frozen for >30 hours during power outage or system failure. The climate-friendly refrigerants used help to protect the environment.



Description	Capacity	Storage capacity	Temperature range	Ext. WxDxH	Interior WxDxH	Cat. No.	NA Price List
<b>120V models</b>							
VWR ULT freezer 352 Eco Premium	500 L (17.7 cu. ft.)	352 (2") cryoboxes	–86...–50 °C	826×938×1966 mm (33×37×77")	606×605×1300 mm (24×24×51") mm	<b>76514-174</b>	<b>\$10,999.99</b>
VWR ULT freezer 528 Eco Premium	700 L (24.7 cu. ft.)	528 (2") cryoboxes	–86...–50 °C	1110×938×1966 mm (44×37×77")	890×605×1300 mm (35×24×51") mm	<b>76514-176</b>	<b>\$11,999.99</b>
<b>240V models</b>							
VWR ULT freezer 352 Eco Premium, 240 V	500 L (17.7 cu. ft.)	352 (2") cryoboxes	–86...–50 °C	826×938×1966 mm (33×37×77")	606×605×1300 mm (24×24×51") mm	<b>76533-504</b>	<b>\$10,999.99</b>
VWR ULT freezer 528 Eco Premium, 240 V	700 L (24.7 cu. ft.)	528 (2") cryoboxes	–86...–50 °C	1110×938×1966 mm (44×37×77")	890×605×1300 mm (35×24×51") mm	<b>76533-506</b>	<b>\$11,999.99</b>