



Moticam X5 PLUS



SAP Code	1100600101322
Sensor	1/3" CMOS Sensor
Imaging area	5.44 x 3.07mm
Capture resolution	4MP
Pixel size (um)	2.0umx2.0um
Scan mode	Progressive
Shutter mode	Rolling Shutter
Operating Temperature	-10 ~ +50°C
Operating Humidity	30~80%RH
Lens mount	CS-Mount
Software	MotiConnect for Android/iOS/Windows/Mac OS Motic Images Plus 3.1 Windows/Mac OS/Linux
Package Includes	12V Power Adapter, CS Ring, Motic 4-dot calibration slide, Dust Cap, Manual, 12mm lens, Adapter 30/38, Macro Tube, Macro Calibration card

INTERFACE & BUTTON FUNCTIONS

USB	USB 5V Output
RJ45	LAN Port, Default IP: 192.168.2.151
WiFi	2.4G/5G (Support 802.11ac), Default IP: 192.168.1.151
DC12V	DC12V/2.5A
LED	Orange LED, Blue LED
RST	Reset, 10s

VIDEO

Image Compression Engine	H264 or MJPEG
Major Stream	1920x1080, 1280x720
Minor Stream	640x360
MJPG	1280x720
Still(Capture)	2560x1440
Frame rate	30fps



*Frames per second under optimal lighting conditions and in compliance with computer technical requirements.