

VWR disposable glove guide

01. CLEANROOM GLOVES

02. EXAM OR SAFETY
GLOVES







VWR disposable glove guide

Visit
**VWR.COM/
VWRGLOVES**
to search
online!

The VWR Collection line of disposable gloves offers a variety of choices depending on your application. Whether your needs are production, laboratory, or industrial focused, we have you covered. Use the guide below to understand the features of our VWR® Disposable Gloves and choose the right level of protection from various material, thickness, and length options.

CLEANROOM GLOVES, POWDER-FREE

	Product Line	Packed	Static Dissipative	Thickness, Finger	Length	Texture	Cuff	Color	Ambidextrous	Packaging	Cat. No. Series
Nitrile Gloves											
	VWR® CERTICLEAN® Class 10 Nitrile Gloves	ISO Class 4 (FED-STD-209E Class 10/M2.5)	Yes	6.3 mil	12 in.	Fingertips	Beaded	White	Yes	100/Bag, 10 Bags/Case	40101-352
	VWR® CERTICLEAN® Class 100 Nitrile Gloves	ISO Class 5 (FED-STD-209E Class 100/M3.5)	Yes	6.3 mil	10 in.	Fingertips	Beaded	White	Yes	100/Bag, 10 Bags/Case	40101-222
	VWR® Maximum Protection Class 10 Nitrile Gloves	ISO Class 4 (FED-S TD-209E Class 10/M2.5)	Yes	5.9 mil	12 in.	Fingers	Tapered, Beaded	White	Yes	100/Bag, 10 Bags/Case	89130-846
	VWR® Advanced Protection Class 100 Nitrile Gloves	ISO Class 5 (FED-STD-209E Class 100/M3.5)	Yes	5.9 mil	12 in.	Fingers	Tapered, Beaded	White	Yes	100/Bag, 10 Bags/Case	89130-836
Latex Gloves											
	VWR® CERTICLEAN® Class 10 Latex Gloves	ISO Class 4 (FED-S TD-209E Class 10/M2.5)	No	8.7 mil	12 in.	Fully Textured	Tapered, Beaded	Natural	Yes	100/Bag, 10 Bags/Case	40101-146
			No	8.7 mil	12 in.	Fully Textured	Tapered, Beaded	Natural	Hand-Specific	50 Pair/Bag, 4 Bags/Case	40101-044
	VWR® CERTICLEAN® Class 100 Latex Gloves	ISO Class 5 (FED-STD-209E Class 100/M3.5)	No	8.7 mil	10 in.	Fully Textured	Beaded	Natural	Yes	100/Bag, 10 Bags/Case	40101-142
			No	8.7 mil	12 in.	Fully Textured	Beaded	Natural	Yes	100/Bag, 10 Bags/Case	40101-158

GENERAL USE GLOVES, POWDER-FREE

	Product Line	Feature	Material	Chemo Rated	Thickness, Finger	Length	Texture	Cuff	Color	Packaging	Cat. No. Series
	VWR® MICROGRIP® PURPLE NITRILE® Poly-Coated Powder-Free Gloves	Polymer Coating for Easy Donning	Nitrile	No	5.5 mil	10 in.	Fully Textured	Beaded	Purple	100/Pack, 10 Packs/Case	40101-346
	VWR® Powder-Free Nitrile Examination Gloves	Everyday Glove, Ideal for Many Applications	Nitrile	Yes	4.7 mil	9 in.	Fully Textured	Rolled, Beaded	Light Blue	100/Pack, 10 Packs/Case	82026-424
	VWR® Finger Textured Nitrile Examination Gloves	Everyday Glove, Ideal for Many Applications	Nitrile	Yes	4.7 mil	9 in.	Finger	Rolled, Beaded	Light Blue	100/Pack, 10 Packs/Case	76411-506
	VWR® Black Nitrile Examination Gloves	Exceptional Barrier Protection, Sleek Look, Fentanyl Approved	Nitrile	Yes	5.5 mil	9.5 in.	Fully Textured	Beaded	Black	100/Pack, 10 Packs/Case	76307-464
	VWR® Ultra Light Black Nitrile Glove	Bulk Packed for Everyday Use	Nitrile	Yes	2.5 mil	9.5 in	Finger	Rolled, Beaded	Black	250/pack (XL=225), 10 packs/case	76457-130
	VWR® Soft Nitrile Examination Gloves	Exceptional Tactile Sensitivity and Comfort	Nitrile	Yes	3.9 mil	9 in.	Finger	Beaded	Blue	100/Pack, 10 Packs/Case	89038-268
	VWR® Nitrile 200 Examination Gloves	Space-Saving and Waste-Reducing Packaging	Nitrile	Yes	3.1 mil	9 in.	Finger	Rolled, Beaded	Blue	200/Pack, 10 Packs/Case	89428-748
	VWR® Nitrile Extended Cuff Exam Gloves	Extended Cuff and Thicker High Risk Protection	Nitrile	Yes	6.1 mil	12 in.	None	Rolled, Beaded	Blue	50/Pack, 10 Packs/Case	89521-896
	VWR Vinyl Nitrile Hybrid Exam Glove	Durable, Cost Effective Blended Nitrile Solution	Vinyl/Nitrile	No	4.0 mil	9.5 in	None	Rolled	Blue	100/pack, 10 packs/case	76456-892
	VWR® Powder-Free Latex Examination Gloves	Natural Rubber Latex For Everyday Use, High Barrier Protection	Latex	Yes	5.5 mil	9.5 in	Fully Textured	Beaded	Natural	100/Pack, 10 Packs/Case	76319-668
	VWR Vinyl Examination Glove	Great for Everyday Use in Low Risk Environments	Vinyl	No	4.5 mil	9.5 in	None	Rolled	Transparent	100/pack, 10 packs/case	76459-214
	VWR® Polyethylene Disposable Gloves	Durable Food and Service Gloves	Polyethylene	No	1.2 mil	12 in.	Fully Textured	N/A	Transparent	100/Pack, 10 Packs/Case	32915-224
						21.5 in. (Elbow Length)				100/Pack, 5 Packs/Case	32915-272



SCIENCE *DELIVERED*

From scientific discovery to scale-up and commercial delivery, you need mission-critical products, services and solutions on a global scale.

Avantor® offers global access to our own portfolio of trusted, quality brands and critical products through our premier delivery platform, VWR®. All of this, combined with infrastructure strategically located to help serve your specific needs, helps move science forward – fast. That's **science delivered**.

Get the mission-critical scientific products and solutions you need from Avantor, delivered through VWR.

➔ vwr.com/science-delivered

 **avantor**™
delivered by **VWR**™

 **avantor**™
delivered by **VWR**™

[VWR.COM](https://vwr.com)

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada unless otherwise noted, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2021 Avantor, Inc. All rights reserved.