


Features:

- Easy to read labeling, including specific chemical information and safety warnings
- Color-coded screw caps for easy identification
- Made with high-quality, low-density polyethylene (LDPE)
- Performance Plastic - thin wall design provides easy to control dispensing with smooth bounce back, won't get stiff over time
- Innovative Flex-Ink - labels won't fade, crack or peel
- Each bottle features a two-piece construction, consisting of the bottle and leak-proof cap with built-in nozzle.
- Wide mouth for easy refilling
- Self-venting design to eliminate leaks and dripping





Chemical Specific NFPA
Hazard Rating Diamond

Dedicated Write in
Space For chemical
Manufacturer Information

eisco Part #CHWB1057

Methanol
CH₃OH

USER FILL IN :

Chemical Manufacturer Name :

Manufacturer address :

Manufacturer website :

Manufacturer phone number : Date :

Chemical Name
Chemical Formula

Chemical Specific
CAS Registry Number

Chemical Specific
GHS Hazard Statement

Chemical Specific
First - Aid Measured

Methanol - CH₃OH

Cas No. 67-56-1 Max Fill Level

DANGER

HIGHLY FLAMMABLE LIQUID AND VAPOUR. TOXIC IF SWALLOWED, IN CONTACT WITH SKIN OR IF INHALED. CAUSES DAMAGE TO EYES BY INGESTION.

- Prevention & Precaution**
Keep away from heat, sparks, and open flames. No smoking, eating or drinking while using. Keep container tightly closed.
Avoid breathing vapours. Use only outdoors or in a well-ventilated area. Wear eye protection. Wear protective gloves and clothing. Wash hand thoroughly after handling.
- Response**
If inhaled : Remove person to fresh air, keep comfortable and breathing. Call a poison center immediately.
If swallowed : Rinse mouth and call a poison center immediately.
If on skin : Rinse with soap and water. Call a poison center if you feel unwell.
In case of fire : Use water spray, alcohol-resistant foam, dry chemical, or carbon dioxide for extinction.

www.eiscolabs.com

Chemical Name
and Formula

Chemical Specific GHS
Signal Word and Pictograms

Chemical Specific
Precautionary statements



Acetone Wash Bottles (Wide Mouth)

Product No.	Capacity	Color	Pack Size
CHWB1010	500 ml	1 Color	Pack of 1
CHWB1010PK2	500 ml	1 Color	Pack of 2
CHWB1010PK6	500 ml	1 Color	Pack of 6
CHWB1011	500 ml	2 Color	Pack of 1
CHWB1011PK2	500 ml	2 Color	Pack of 2
CHWB1011PK6	500 ml	2 Color	Pack of 6
CHWB1012	500 ml	4 Color	Pack of 1
CHWB1012PK2	500 ml	4 Color	Pack of 2
CHWB1012PK6	500 ml	4 Color	Pack of 6
CHWB1013	1000 ml	1 Color	Pack of 1
CHWB1013PK2	1000 ml	1 Color	Pack of 2
CHWB1013PK6	1000 ml	1 Color	Pack of 6
CHWB1014	1000 ml	2 Color	Pack of 1
CHWB1014PK2	1000 ml	2 Color	Pack of 2
CHWB1014PK6	1000 ml	2 Color	Pack of 6
CHWB1015	1000 ml	4 Color	Pack of 1
CHWB1015PK2	1000 ml	4 Color	Pack of 2
CHWB1015PK6	1000 ml	4 Color	Pack of 6



Distilled Water Wash Bottles (Wide Mouth)

Product No.	Capacity	Color	Pack Size
CHWB1019	500 ml	1 Color	Pack of 1
CHWB1019PK2	500 ml	1 Color	Pack of 2
CHWB1019PK6	500 ml	1 Color	Pack of 6
CHWB1020	500 ml	2 Color	Pack of 1
CHWB1020PK2	500 ml	2 Color	Pack of 2
CHWB1020PK6	500 ml	2 Color	Pack of 6
CHWB1021	500 ml	4 Color	Pack of 1
CHWB1021PK2	500 ml	4 Color	Pack of 2
CHWB1021PK6	500 ml	4 Color	Pack of 6
CHWB1022	1000 ml	1 Color	Pack of 1
CHWB1022PK2	1000 ml	1 Color	Pack of 2
CHWB1022PK6	1000 ml	1 Color	Pack of 6
CHWB1023	1000 ml	2 Color	Pack of 1
CHWB1023PK2	1000 ml	2 Color	Pack of 2
CHWB1023PK6	1000 ml	2 Color	Pack of 6
CHWB1024	1000 ml	4 Color	Pack of 1
CHWB1024PK2	1000 ml	4 Color	Pack of 2
CHWB1024PK6	1000 ml	4 Color	Pack of 6



Isopropanol Wash Bottles (Wide Mouth)

Product No.	Capacity	Color	Pack Size
CHWB1028	500 ml	1 Color	Pack of 1
CHWB1028PK2	500 ml	1 Color	Pack of 2
CHWB1028PK6	500 ml	1 Color	Pack of 6
CHWB1029	500 ml	2 Color	Pack of 1
CHWB1029PK2	500 ml	2 Color	Pack of 2
CHWB1029PK6	500 ml	2 Color	Pack of 6
CHWB1030	500 ml	4 Color	Pack of 1
CHWB1030PK2	500 ml	4 Color	Pack of 2
CHWB1030PK6	500 ml	4 Color	Pack of 6
CHWB1031	1000 ml	1 Color	Pack of 1
CHWB1031PK2	1000 ml	1 Color	Pack of 2
CHWB1031PK6	1000 ml	1 Color	Pack of 6
CHWB1032	1000 ml	2 Color	Pack of 1
CHWB1032PK2	1000 ml	2 Color	Pack of 2
CHWB1032PK6	1000 ml	2 Color	Pack of 6
CHWB1033	1000 ml	4 Color	Pack of 1
CHWB1033PK2	1000 ml	4 Color	Pack of 2
CHWB1033PK6	1000 ml	4 Color	Pack of 6



Ethanol Wash Bottles (Wide Mouth)

Product No.	Capacity	Color	Pack Size
CHWB1037	500 ml	1 Color	Pack of 1
CHWB1037PK2	500 ml	1 Color	Pack of 2
CHWB1037PK6	500 ml	1 Color	Pack of 6
CHWB1038	500 ml	2 Color	Pack of 1
CHWB1038PK2	500 ml	2 Color	Pack of 2
CHWB1038PK6	500 ml	2 Color	Pack of 6
CHWB1039	500 ml	4 Color	Pack of 1
CHWB1039PK2	500 ml	4 Color	Pack of 2
CHWB1039PK6	500 ml	4 Color	Pack of 6
CHWB1040	1000 ml	1 Color	Pack of 1
CHWB1040PK2	1000 ml	1 Color	Pack of 2
CHWB1040PK6	1000 ml	1 Color	Pack of 6
CHWB1041	1000 ml	2 Color	Pack of 1
CHWB1041PK2	1000 ml	2 Color	Pack of 2
CHWB1041PK6	1000 ml	2 Color	Pack of 6
CHWB1042	1000 ml	4 Color	Pack of 1
CHWB1042PK2	1000 ml	4 Color	Pack of 2
CHWB1042PK6	1000 ml	4 Color	Pack of 6



Sodium Hypochlorite Wash Bottles (Wide Mouth)

Product No.	Capacity	Color	Pack Size
CHWB1046	500 ml	1 Color	Pack of 1
CHWB1046PK2	500 ml	1 Color	Pack of 2
CHWB1046PK6	500 ml	1 Color	Pack of 6
CHWB1047	500 ml	2 Color	Pack of 1
CHWB1047PK2	500 ml	2 Color	Pack of 2
CHWB1047PK6	500 ml	2 Color	Pack of 6
CHWB1048	500 ml	4 Color	Pack of 1
CHWB1048PK2	500 ml	4 Color	Pack of 2
CHWB1048PK6	500 ml	4 Color	Pack of 6
CHWB1049	1000 ml	1 Color	Pack of 1
CHWB1049PK2	1000 ml	1 Color	Pack of 2
CHWB1049PK6	1000 ml	1 Color	Pack of 6
CHWB1050	1000 ml	2 Color	Pack of 1
CHWB1050PK2	1000 ml	2 Color	Pack of 2
CHWB1050PK6	1000 ml	2 Color	Pack of 6
CHWB1051	1000 ml	4 Color	Pack of 1
CHWB1051PK2	1000 ml	4 Color	Pack of 2
CHWB1051PK6	1000 ml	4 Color	Pack of 6



Methanol Wash Bottles (Wide Mouth)

Product No.	Capacity	Color	Pack Size
CHWB1061	500 ml	1 Color	Pack of 1
CHWB1061PK2	500 ml	1 Color	Pack of 2
CHWB1061PK6	500 ml	1 Color	Pack of 6
CHWB1062	1000 ml	1 Color	Pack of 1
CHWB1062PK2	1000 ml	1 Color	Pack of 2
CHWB1062PK6	1000 ml	1 Color	Pack of 6
CHWB1056	500 ml	2 Color	Pack of 1
CHWB1056PK2	500 ml	2 Color	Pack of 2
CHWB1056PK6	500 ml	2 Color	Pack of 6
CHWB1057	500 ml	4 Color	Pack of 1
CHWB1057PK2	500 ml	4 Color	Pack of 2
CHWB1057PK6	500 ml	4 Color	Pack of 6
CHWB1059	1000 ml	2 Color	Pack of 1
CHWB1059PK2	1000 ml	2 Color	Pack of 2
CHWB1059PK6	1000 ml	2 Color	Pack of 6
CHWB1060	1000 ml	4 Color	Pack of 1
CHWB1060PK2	1000 ml	4 Color	Pack of 2
CHWB1060PK6	1000 ml	4 Color	Pack of 6



Chemical Wash Bottles, No Labels (Wide Mouth)

Product No.	Capacity	Pack Size
CHWB1001	250 ml	Pack of 1
CHWB1001PK2	250 ml	Pack of 2
CHWB1001PK6	250 ml	Pack of 6
CHWB1002	500 ml	Pack of 1
CHWB1002PK2	500 ml	Pack of 2
CHWB1002PK6	500 ml	Pack of 6
CHWB1003	1000 ml	Pack of 1
CHWB1003PK2	1000 ml	Pack of 2
CHWB1003PK6	1000 ml	Pack of 6



Chemical Wash Bottles, Blank Labels (Wide Mouth)

Product No.	Capacity	Pack Size
CHWB1004	250 ml	Pack of 1
CHWB1004PK2	250 ml	Pack of 2
CHWB1004PK6	250 ml	Pack of 6
CHWB1005	500 ml	Pack of 1
CHWB1005PK2	500 ml	Pack of 2
CHWB1005PK6	500 ml	Pack of 6
CHWB1006	1000 ml	Pack of 1
CHWB1006PK2	1000 ml	Pack of 2
CHWB1006PK6	1000 ml	Pack of 6