






REVISION NO	REVISION DATE	DESCRIPTION	REVISED BY	APPROVED BY

ROUND WIDE MOUTH BOTTLE SPECIFICATION

VWR-US CAT NO	VWR EUROPE CAT NO	MATERIAL	COLOR	CAPACITY	APROX BRIM CAPACITY (ml)	CAP SIZE	HEIGHT WITH OUT CAP (mm)	BODY DIA (mm)
414004-186	215-1889	HDPE	TRANSLUCENT	4L	4.5L	100-415	298	163
414004-187	215-1898	PP	TRANSLUCENT	2L	2.28L	100-415	240	119
414004-188	215-1899	PP	TRANSLUCENT	4L	4.5L	100-415	298	163
414004-205	215-1918	HDPE	TRANSLUCENT	2L	2.28L	100-415	240	119

NOTE :
1) PP MATERIAL BOTTLES ARE AUTOCLABLE.
2) 4L BOTTLE HAS INTEGRATED HANDLE DESIGN.

Closure Application Torque		
Closure Size	Minimum Torque Kg cm	Maximum Torque Kg cm
13-415	6	8
20-415	11	16
24-415	13	19
28-415	18	25
38-415	25	35
38-430	31	38
43-415	32	44
48-415	34	48
53-415	38	52
53T	43	60
63-415	46	64
83T	69	96
100-415	54	74



TITLE : WIDE MOUTH ROUND BOTTLE

Material : AS PER TABLE

Rev : 0 SCALE : NOT TO SCALE



DRAWN BY	
DRAWN DATE	2/2/2009
DRAWING NO	

