

# VWR<sup>®</sup> MOBILE LABORATORY FURNITURE SOLUTIONS



# MOBILE LABORATORY FURNITURE SOLUTIONS



## LEARN MORE ABOUT VWR FURNITURE PRODUCTS, SERVICES, AND CUSTOMER SUPPORT

### VWR REDISHIP

VWR REDISHIP offers you many laboratory furniture items in a matter of days\*. If you're ready to renovate your lab at a moment's notice, or if you need to quickly replace lab chairs, tables, carts, shelving, steel casework, or fume hoods, VWR REDISHIP is ready when you are.

Now, to help you easily identify all REDISHIP products in this catalog, a blue truck icon appears next to the catalog number!



\*Most REDISHIP items ship within 72 hours of order. All items are subject to availability at time of order. Large-quantity orders may require longer lead times. 72-hour shipping time is based on immediate credit approval.

The VWR mobile laboratory makes changing or expanding a work space easier than ever.

Need to change a single-sided station into a double-sided one? How about deeper shelves or electrical inserts? It's a snap to rearrange an existing workspace or start a new one from scratch. Everything from storage to lighting to keyboard trays can be placed wherever you want, whenever you want.

VWR\* Mobile Laboratory Furniture adapts to your growing needs and helps increase productivity, worker comfort, and workspace efficiency. And the modular, intelligent design helps take the hassle out of change. All parts and accessories will attach securely, so no parts are wasted. You get the flexibility to create the system you need and the ability to modify your work space as necessary.

With a wide variety of benches, configurations, and add-on accessories, our furniture conforms to your exact specifications.

Simply visit [vwr.com](http://vwr.com) for details, or contact your VWR Furniture Specialist at **1.800.932.5000**.

## CONTENTS

Atlas Mobile Laboratory Benches.....	4-9
VWR Laboratory Benches.....	33-43
Worksurfaces .....	10-12
Suspended Cabinets.....	13
Shelves.....	14-15
Overhead Cabinets.....	16-17
Fixtures/Power Bars.....	18-19, 44, 51-53
Technical Sheets.....	48-50
Accessories.....	20-24, 45, 52-54

# VWR® Atlas Mobile Carts




**PLEASE NOTE:**  
Atlas Mobile Carts accommodate loads up to 1000 pounds. Perfect for mobilizing heavy equipment.

## Parts and Accessories

A. 72" high column.....14	F. Tack panel .....20	K. VWR Contour suspended cabinet...13
B. Standard shelf.....14	G. Worksurface ..... 10-12	L. Standard caster ..... 5-9
C. Overhead storage cabinet .....16	H. Front support leg ..... 5-9	M. Data outlets.....19
D. Powered backplash.....19	I. Modesty panel.....5-9	N. Service fixtures.....18
E. Center column.....14	J. Standard shelf .....14	

# How to Select & Specify Your Atlas Mobile Cart



# Preconfigured Atlas Mobile Carts with Casters

Preconfigured Atlas mobile carts include two (2) columns with bases, modesty panel, front apron rail and Kemresin® epoxy resin worksurface. Carts are 36" high x 31" deep and are available in 36", 48" and 72" widths. Dimensions include caster height. 72"-wide carts are equipped with a center leg for added support. Weight Load: 1000 lbs. Size: As selected. Color: Light Neutral.

48"-wide carts include one (1) Alpha storage unit. 72"-wide carts include two (2) Alpha storage units. Additional catalog accessories may be added.

Glides are included with carts and are interchangeable with standard casters. In the event casters are replaced with glides, cart height will decrease by 3".



**VWRAMC001** 

Atlas Mobile Cart  
36" Wide x 36" High x 31" Deep  
Kemresin Epoxy Resin Worksurface



**VWRAMC002** 

Atlas Mobile Cart  
48" Wide x 36" High x 31" Deep  
Kemresin Epoxy Resin Worksurface  
Alpha Contour Drawers -19"H x 18"W



**VWRAMC003** 

Atlas Mobile Cart  
72" Wide x 36" High x 31" Deep  
Kemresin Epoxy Resin Worksurface  
Alpha Contour Drawers -19"H x 22"D x 18"W  
Alpha Contour Drawer/Cabinet -25"H x 22"D x 18"W

# Step 1 - Specify Atlas Mobile Cart(s) with Casters

Atlas mobile carts with standard casters include two (2) columns with bases, modesty panel and front apron rail. Carts are available in 31" and 34.5" depths and 30", 36" and 76" heights. Dimensions include caster height. 72"- wide carts are equipped with a center leg for added support. Weight Capacity: 1000 lbs. See steps 3 - 8 for worksurface options. Size: As selected. Color: Light Neutral.

Glides are included with carts and are interchangeable with standard casters. In the event casters are replaced with glides, cart height will decrease by 3".



Worktop Not Included

## Sitting Height w/Std Casters

	31" Depth	34.5" Depth
33" W x 30" H	AMC313030C	AMC353030C
39" W x 30" H	AMC313630C	AMC353630C
51" W x 30" H	AMC314830C	AMC354830C
63" W x 30" H	AMC316030C	AMC356030C
75" W x 30" H	AMC317230C	AMC357230C

Note: 72" Cart Provided with Center Column

## Standing/Bench Height w/Std Casters

	31" Depth	34.5" Depth
33" W x 36" H	AMC313036C	AMC353036C
39" W x 36" H	AMC313636C	AMC353636C
51" W x 36" H	AMC314836C	AMC354836C
63" W x 36" H	AMC316036C	AMC356036C
75" W x 36" H	AMC317236C	AMC357236C

Note: 72" Cart Provided with Center Column



Worktop Not Included

## High-Column w/ Std Casters

	31" Depth	34.5" Depth
33" W x 76" H	AMC313076C	AMC353076C
39" W x 76" H	AMC313676C	AMC353676C
51" W x 76" H	AMC314876C	AMC354876C
63" W x 76" H	AMC316076C	AMC356076C
75" W x 76" H	AMC317276C	AMC357276C

Note: 72" Cart Provided with Center Column

# Step 2 - Specify Atlas Mobile Cart(s) with Heavy-Duty Leveling Casters

Atlas mobile carts with heavy-duty leveling casters include two (2) columns with bases, modesty panel and front apron rail. Carts are available in 31" and 34.5" depths and 30", 36" and 76" heights. Dimensions include caster height. 72"-wide carts are equipped with a center leg for added support. Weight Capacity: 1000 lbs. See steps 3 - 8 for worksurface options. Size: As selected. Color: Light Neutral.



Worktop Not Included

## Sitting Height w/ HD Leveling Casters

	31" Depth	34.5" Depth
33" W x 30" H	AMC313030H	AMC353030H
39" W x 30" H	AMC313630H	AMC353630H
51" W x 30" H	AMC314830H	AMC354830H
63" W x 30" H	AMC316030H	AMC356030H
75" W x 30" H	AMC317230H	AMC357230H

Note: 72" Cart Provided with Center Column

## Standing/Bench Height w/ HD Leveling

	31" Depth	34.5" Depth
33" W x 36" H	AMC313036H	AMC353036H
39" W x 36" H	AMC313636H	AMC353636H
51" W x 36" H	AMC314836H	AMC354836H
63" W x 36" H	AMC316036H	AMC356036H
75" W x 36" H	AMC317236H	AMC357236H

Note: 72" Cart Provided with Center Column



Worktop Not Included

## High-Column Carts w/ HD Leveling

	31" Depth	34.5" Depth
33" W x 76" H	AMC313076H	AMC353076H
39" W x 76" H	AMC313676H	AMC353676H
51" W x 76" H	AMC314876H	AMC354876H
63" W x 76" H	AMC316076H	AMC356076H
75" W x 76" H	AMC317276H	AMC357276H

Note: 72" Cart Provided with Center Column

# Step 1 - Specify Atlas Mobile Cart(s) w/o Front Support - Standard Casters



Atlas mobile carts with standard casters include two (2) columns with bases, modesty panel and front apron rail without front supports. Carts are available in 31" and 34.5" depths and 30", 36" and 76" heights. Dimensions include caster height. Weight Capacity: 600 lbs. See steps 3 - 8 for worksurface options. Size: As selected. Color: Light Neutral.

Glides are included with carts and are interchangeable with standard casters. In the event casters are replaced with glides, cart height will decrease by 3".



Worktop Not Included

## Sitting Height w/ Std Casters

	31" Depth	34.5" Depth
33" W x 30" H	AMC313030CN	AMC353030CN
39" W x 30" H	AMC313630CN	AMC353630CN
51" W x 30" H	AMC314830CN	AMC354830CN
63" W x 30" H	AMC316030CN	AMC356030CN
75" W x 30" H	AMC317230CN	AMC357230CN

Note: 72" Cart Provided with Center Column

## Standing/Bench Height w/ Std Casters

	31" Depth	34.5" Depth
33" W x 36" H	AMC313036CN	AMC353036CN
39" W x 36" H	AMC313636CN	AMC353636CN
51" W x 36" H	AMC314836CN	AMC354836CN
63" W x 36" H	AMC316036CN	AMC356036CN
75" W x 36" H	AMC317236CN	AMC357236CN

Note: 72" Cart Provided with Center Column



Worktop Not Included

## High-Column w/ Std Casters

	31" Depth	34.5" Depth
33" W x 76" H	AMC313076CN	AMC353076CN
39" W x 76" H	AMC313676CN	AMC353676CN
51" W x 76" H	AMC314876CN	AMC354876CN
63" W x 76" H	AMC316076CN	AMC356076CN
75" W x 76" H	AMC317276CN	 AMC357276CN

Note: 72" Cart Provided with Center Column

# Step 2 - Specify Atlas Mobile Cart(s) w/o Front Support - HD Leveling Casters

Atlas mobile carts with heavy-duty leveling casters include two (2) columns with bases, modesty panel and front apron rail without front supports. Carts are available in 31" and 34.5" depths and 30", 36" and 76" heights. Dimensions include caster height. 72"-wide carts are equipped with a center leg for added support. Weight Capacity: 600 lbs. See steps 3 - 8 for worksurface options. Size: As selected. Color: Light Neutral.



Worktop Not Included

## Sitting Height w/ HD Leveling Casters

	31" Depth	34.5" Depth
33" W x 30" H	AMC313030HN	AMC353030HN
39" W x 30" H	AMC313630HN	AMC353630HN
51" W x 30" H	AMC314830HN	AMC354830HN
63" W x 30" H	AMC316030HN	AMC356030HN
75" W x 30" H	AMC317230HN	AMC357230HN

Note: 72" Cart Provided with Center Column

## Standing/Bench Height w/ HD Leveling Casters

	31" Depth	34.5" Depth
33" W x 36" H	AMC313036HN	AMC353036HN
39" W x 36" H	AMC313636HN	AMC353636HN
51" W x 36" H	AMC314836HN	AMC354836HN
63" W x 36" H	AMC316036HN	AMC356036HN
75" W x 36" H	AMC317236HN	AMC357236HN

Note: 72" Cart Provided with Center Column



Worktop Not Included

## High-Column Carts w/ HD Leveling Casters


	31" Depth	34.5" Depth
33" W x 76" H	AMC313076HN	AMC353076HN
39" W x 76" H	AMC313676HN	AMC353676HN
51" W x 76" H	AMC314876HN	AMC354876HN
63" W x 76" H	AMC316076HN	AMC356076HN
75" W x 76" H	AMC317276HN	AMC357276HN

Note: 72" Cart Provided with Center Column

## Step 3 - Specify Kemresin® Epoxy Resin Worksurface for Sitting/Bench Height



Kemresin® worksurfaces are composed of molded, modified epoxy resin. Each worksurface has three square edges and one edge with a 1/16" chamfer. Size: 1"H x Selected Width x Selected Depth. Color: Black.




### Black Kemresin® Epoxy Resin Worksurfaces

	31" Depth	34.5" Depth
33" W	AMCRWB3130M	AMCRWB3530M
39" W	AMCRWB3136M	AMCRWB3536M
51" W	 AMCRWB3148M	 AMCRWB3548M
63" W	 AMCRWB3160M	 AMCRWB3560M
75" W	AMCRWB3172M	AMCRWB3572M

## Step 4 - Specify Kemresin® Epoxy Resin Worksurface for High-Column

Kemresin® worksurfaces are composed of molded, modified epoxy resin. Each worksurface has three square edges and one edge with a 1/16" chamfer. Two or three cut-outs accommodate 76" high columns, depending on bench width. High-Column worksurfaces may be used with sitting/bench height carts when height-adjustment is required. Size: 1"H x Selected Width x Selected Depth. Color: Black.



### Black Kemresin® Epoxy Resin Worksurfaces

	31" Depth	34.5" Depth
33" W	AMCRWB3130H	AMCRWB3530H
39" W	 AMCRWB3136H	AMCRWB3536H
51" W	 AMCRWB3148H	 AMCRWB3548H
63" W	 AMCRWB3160H	AMCRWB3560H
75" W	 AMCRWB3172H	 AMCRWB3572H



## Step 5 - Specify Plastic Laminate Worksurface for Sitting/Bench Height

Laminate worksurfaces are composed of high-pressure plastic laminate. Each worksurface has three square edges and one full 180-degree wrapped front edge. Size: 11/16"H x Selected Width x Selected Depth. Color: Black.



### Black Plastic Laminate Worksurfaces

	31" Depth	34.5" Depth
33" W	AMCLWB3130M	AMCLWB3530M
39" W	AMCLWB3136M	AMCLWB3536M
51" W	 AMCLWB3148M	AMCLWB3548M
63" W	 AMCLWB3160M	AMCLWB3560M
75" W	AMCLWB3172M	AMCLWB3572M

## Step 6 - Specify Plastic Laminate Worksurface for High-Column

Laminate worksurfaces are composed of high-pressure plastic laminate. Each worksurface has three square edges and one full 180-degree wrapped front edge. Two or three cut-outs accommodate 76" high columns, depending on bench width. High-Column worksurfaces may be used with sitting/bench height carts when height-adjustment is required. Size: 11/16"H x Selected Width x Selected Depth. Color: Black.




### Black Plastic Laminate Worksurfaces

	31" Depth	34.5" Depth
33" W	AMCLWB3130H	AMCLWB3530H
39" W	AMCLWB3136H	AMCLWB3536H
51" W	AMCLWB3148H	AMCLWB3548H
63" W	AMCLWB3160H	AMCLWB3560H
75" W	AMCLWB3172H	AMCLWB3572H

## Step 7 - Specify Black Phenolic Resin Worksurface for Sitting/Bench Height

Phenolic resin worksurfaces are composed of cellulose fiber-reinforced, phenolic resin core with a highly-crosslinked, polyurethane, copolymer surface. Each worksurface has three square edges and one edge with a 1/16" chamfer. Worksurfaces for High-Column Carts may be used with Mid-Column carts when height adjustment is required. Size: 1"H x Selected Width x Selected Depth. Color: Black.




### Phenolic Resin Worksurfaces

	31" Depth	34.5" Depth
33" W	AMCTWB3130M	AMCTWB3530M
39" W	AMCTWB3136M	AMCTWB3536M
51" W	AMCTWB3148M	AMCTWB3548M
63" W	AMCTWB3160M	AMCTWB3560M
75" W	AMCTWB3172M	AMCTWB3572M

## Step 8 - Specify Black Phenolic Resin Worksurface for High-Column

Phenolic resin worksurfaces are composed of cellulose fiber-reinforced, phenolic resin core with a highly-crosslinked, polyurethane, copolymer surface. Two or three cut-outs accommodate 76" high columns, depending on bench width. Each worksurface has three square edges and one edge with a 1/16" chamfer. High-Column worksurfaces may be used with sitting/bench height carts when height-adjustment is required. Size: 1"H x Selected Width x Selected Depth. Color: Black.




### Phenolic Resin Worksurfaces

	31" Depth	34.5" Depth
33" W	AMCTWB3130H	AMCTWB3530H
39" W	AMCTWB3136H	AMCTWB3536H
51" W	AMCTWB3148H	AMCTWB3548H
63" W	AMCTWB3160H	AMCTWB3560H
75" W	AMCTWB3172H	AMCTWB3572H


# VWR® Contour™ Suspended Cabinets




All suspended height cabinets are 19" or 25" high, 22" deep and 18" wide. Cabinets are designed to hang beneath table frames. Cabinets are phosphate-coated with a baked, chemical-resistant, powder-coat finish. Color: light neutral.



## 19" Drawer Cabinets



Three Drawer w/ Pull Board  
18" W CFD-5002-B




One Drawer w/ Pull Board  
One File Drawer w/ Hangers  
18" W CFD-5852-B




Swings Right  
18" W CFC-5032-A  
Swing Right w/ Pull Board  
18" W CFC-5032-B

Swings Left  
18" W CFC-5032-L  
Swing Left w/ Pull Board  
18" W CFC-5032-BL


## 25" Drawer Cabinets



Three Drawers  
18" W CFD-4002-A  
Three Drawers w/ Pull Board  
18" W CFD-4002-B




Two File Drawers with Hangers  
18" W CFD-4812-A




Three Drawers w/ Pull Board  
One File Drawer w/ Hangers  
18" W CFD-4852-B

## 25" Cupboard Cabinets



Swings Right  
18" W CFC-4032-A  
Swing Left  
18" W CFC-4032-L




Swings Right w/ Pull Board  
18" W CFC-4032-B  
Swings Left w/ Pull Board  
18" W CFC-4032-BL

## 25" Cupboard and Drawer Cabinets



Swings Right  
18" W CFE-4402-A  
Swings Left  
18" W CFE-4402-L



Swings Right w/ Pull Board  
18" W CFE-4402-B  
Swings Left w/ Pull Board  
18" W CFE-4402-BL

# VWR® Atlas Modular Shelving



## Single Atlas Uprights


Individual Atlas uprights bolt directly to a new or existing counter top. A ceiling trim ring is provided with uprights 60"-96" high. Mounting hardware included. Color: Light Neutral. Column Inserts: Black.



24"H x 8"W x 1¾"D		VWRC24
30"H x 8"W x 1¾"D		VWRC30
36"H x 8"W x 1¾"D		VWRC36
48"H x 8"W x 1¾"D		VWRC48
54"H x 8"W x 1¾"D		VWRC54
60"H x 8"W x 1¾"D		VWRC60
66"H x 8"W x 1¾"D		VWRC66
72"H x 8"W x 1¾"D		VWRC72
84"H x 8"W x 1¾"D		VWRC84
96"H x 8"W x 1¾"D		VWRC96

## Single Atlas Cross Braces


Atlas cross braces support Atlas uprights and stabilize an Atlas shelf assembly. One cross brace is required for all single- and double-sided assemblies. Mounting hardware included. Color: Light Neutral.



36"L x 3"H x 6⅞"W		VWRCB36
42"L x 3"H x 6⅞"W		VWRCB42
48"L x 3"H x 6⅞"W		VWRCB48
54"L x 3"H x 6⅞"W		VWRCB54
60"L x 3"H x 6⅞"W		VWRCB60
66"L x 3"H x 6⅞"W		VWRCB66
72"L x 3"H x 6⅞"W		VWRCB72

## Standard Shelves (No Lip)


Standard Atlas shelves are designed to be infinitely adjustable along vertical Atlas uprights. Color: Light Neutral. Size: As selected.



	12" Depth	15" Depth	18" Depth	24" Depth
36" W	VWRNL1236	VWRNL1536	VWRNL1836	VWRNL2436
42" W	VWRNL1242	VWRNL1542	VWRNL1842	VWRNL2442
48" W	VWRNL1248	VWRNL1548	VWRNL1848	VWRNL2448
54" W	VWRNL1254	VWRNL1554	VWRNL1854	VWRNL2454
60" W	VWRNL1260	VWRNL1560	VWRNL1860	VWRNL2460
66" W	VWRNL1266	VWRNL1566	VWRNL1866	VWRNL2466
72" W	VWRNL1272	VWRNL1572	VWRNL1872	VWRNL2472

## Single-Lipped Shelves


Single-lipped Atlas shelves are designed to be infinitely adjustable along vertical Atlas uprights. Lipped shelves can be used with the lip in the front or the back and can be switched in the field. Color: Light Neutral. Size: As selected.



	12" Depth	15" Depth	18" Depth	24" Depth
36" W	VWR1L1236	VWR1L1536	VWR1L1836	VWR1L2436
42" W	VWR1L1242	VWR1L1542	VWR1L1842	VWR1L2442
48" W	VWR1L1248	VWR1L1548	VWR1L1848	VWR1L2448
54" W	VWR1L1254	VWR1L1554	VWR1L1854	VWR1L2454
60" W	VWR1L1260	VWR1L1560	VWR1L1860	VWR1L2460
66" W	VWR1L1266	VWR1L1566	VWR1L1866	VWR1L2466
72" W	VWR1L1272	VWR1L1572	VWR1L1872	VWR1L2472

## Double-Lipped Shelves


Double-lipped Atlas shelves are designed to be infinitely adjustable along vertical Atlas uprights. Color: Light Neutral. Size: As selected.



	12" Depth	15" Depth	18" Depth	24" Depth
36" W	VWR2L1236	VWR2L1536	VWR2L1836	VWR2L2436
42" W	VWR2L1242	VWR2L1542	VWR2L1842	VWR2L2442
48" W	VWR2L1248	VWR2L1548	VWR2L1848	VWR2L2448
54" W	VWR2L1254	VWR2L1554	VWR2L1854	VWR2L2454
60" W	VWR2L1260	VWR2L1560	VWR2L1860	VWR2L2460
66" W	VWR2L1266	VWR2L1566	VWR2L1866	VWR2L2466
72" W	VWR2L1272	VWR2L1572	VWR2L1872	VWR2L2472

## Specify End-Mounted Shelves

End-mounted Atlas shelves are designed to be infinitely adjustable along vertical Atlas uprights. Mounting hardware included. Color: Light Neutral. Size: As selected.



12"D x 24"W	VWRES1224
12"D x 30"W	VWRES1230




### Shelf Brackets

Atlas shelf brackets are designed to be infinitely adjustable along vertical Atlas upright. Select style of bracket as required. Select depth of bracket according to shelf depth. Color: Light Neutral. Size: As selected.

#### High-Side Shelf Brackets for Single-Shelf Applications

12" High-side bracket, pair		VWRBH12
15" High-side bracket, pair		VWRBH15
18" High-side bracket, pair		VWRBH18
24" High-side bracket, pair		VWRBH24



These brackets may be used upright or inverted



#### Low Profile Shelf Brackets for Single-Shelf Applications

12" Low-profile bracket, pair		VWRBL12
15" Low-profile bracket, pair		VWRBL15
18" Low-profile bracket, pair		VWRBL18
24" Low-profile bracket, pair		VWRBL24

These brackets for under mount only



### Umbrella Shelves

Center-mounted Atlas shelves are designed to be infinitely adjustable, mounted evenly between Atlas uprights. Umbrella shelves are mounted at the top of uprights and replace the cross brace when used. Color: Light Neutral. Size: As selected.

#### Parts and Accessories


- A. Umbrella shelf
- B. Umbrella Bracket
- C. Center-mounted standard shelf
- D. High-side center-mounting bracket
- E. Center-mounted standard shelf
- F. Center-mounted standard shelf
- G. Low-profile center-mounting bracket

Standard shelves may be mounted in the center of the column using center-mounted shelf brackets.

	12" Depth	15" Depth	18" Depth	24" Depth	36" Depth
36" W	VWRNLU1236	VWRNLU1536	VWRNLU1836	VWRNLU2436	VWRNLU3636
42" W	VWRNLU1242	VWRNLU1542	VWRNLU1842	VWRNLU2442	VWRNLU3642
48" W	VWRNLU1248	VWRNLU1548	VWRNLU1848	VWRNLU2448	VWRNLU3648
54" W	VWRNLU1254	VWRNLU1554	VWRNLU1854	VWRNLU2454	VWRNLU3654
60" W	VWRNLU1260	VWRNLU1560	VWRNLU1860	VWRNLU2460	VWRNLU3660
66" W	VWRNLU1266	VWRNLU1566	VWRNLU1866	VWRNLU2466	VWRNLU3666
72" W	VWRNLU1272	VWRNLU1572	VWRNLU1872	VWRNLU2472	VWRNLU3672

### Center-Mounted Shelf Brackets

Center-mounted Atlas shelf brackets are designed to be infinitely adjustable along vertical Atlas uprights. Select style of bracket as required. Select depth of bracket according to shelf depth. Color: Light Neutral. Size: As selected.



#### High-Side, Center-Mounted Shelf Brackets

12" High-side bracket, pair		VWRBHC12
15" High-side bracket, pair		VWRBHC15
18" High-side bracket, pair		VWRBHC18
24" High-side bracket, pair		VWRBHC24
36" High-side bracket, pair		VWRBHC36

These brackets may be used upright or inverted

#### Low-Profile, Center-Mounted Shelf Brackets


12" Low-profile bracket, pair		VWRBLC12
15" Low-profile bracket, pair		VWRBLC15
18" Low-profile bracket, pair		VWRBLC18
24" Low-profile bracket, pair		VWRBLC24
36" Low-profile bracket, pair		VWRBLC36

These brackets for under mount only

#### Umbrella Shelf Brackets

12" Umbrella bracket, pair		VWRBU12
15" Umbrella bracket, pair		VWRBU15
18" Umbrella bracket, pair		VWRBU18
24" Umbrella bracket, pair		VWRBU24
36" Umbrella bracket, pair		VWRBU36

Designed to stabilize top shelf and replace cross brace




# VWR® Atlas Overhead Storage Cabinets



## Overhead Storage Cabinets

Center-mounted Atlas overhead storage cabinets are used in place of open shelving when enclosed storage and/or resistance to dust is required. For double-sided overhead storage without pass thru, order two overhead storage units in the size required and mount back to back. Includes hardware for single- and double-sided mount. Weight capacity: 40 lbs./ft.2 per bottom and top surfaces, with a maximum of 600 lbs. Color: Light Neutral. Size: 16<sup>3</sup>/<sub>8</sub>"H x 15"D x Selected Width.



### 16<sup>3</sup>/<sub>8</sub>" High Overhead Storage with Sliding Acrylic Doors

36" Overhead Storage	34"W x 15"D	VWR407520A
42" Overhead Storage	40"W x 15"D	VWR407524A
48" Overhead Storage	46"W x 15"D	VWR407521A
54" Overhead Storage	52"W x 15"D	VWR407525A
60" Overhead Storage	58"W x 15"D	VWR407522A
66" Overhead Storage	64"W x 15"D	VWR407526A
72" Overhead Storage	70"W x 15"D	VWR407523A

### 16<sup>3</sup>/<sub>8</sub>" High Overhead Storage with Sliding Safety Glass Doors


36" Overhead Storage	34"W x 15"D	VWR407520G
42" Overhead Storage	40"W x 15"D	VWR407524G
48" Overhead Storage	46"W x 15"D	VWR407521G
54" Overhead Storage	52"W x 15"D	VWR407525G
60" Overhead Storage	58"W x 15"D	VWR407522G
66" Overhead Storage	64"W x 15"D	VWR407526G
72" Overhead Storage	70"W x 15"D	VWR407523G

### 16<sup>3</sup>/<sub>8</sub>" High Overhead Storage with Sliding Safety Glass Doors

36" Overhead Storage	34"W x 15"D	VWR407520M
42" Overhead Storage	40"W x 15"D	VWR407524M
48" Overhead Storage	46"W x 15"D	VWR407521M
54" Overhead Storage	52"W x 15"D	VWR407525M
60" Overhead Storage	58"W x 15"D	VWR407522M
66" Overhead Storage	64"W x 15"D	VWR407526M
72" Overhead Storage	70"W x 15"D	VWR407523M

## Overhead Storage Cabinets

Atlas overhead storage cabinets are used in place of open shelving when enclosed storage and/or resistance to dust is required. Shelves are adjustable and mounting hardware is included. For double-sided overhead storage without pass thru, order two overhead storage units in the size required and mount back to back. Includes hardware for single- and double-sided mount. Weight capacity: 40 lbs./ft.2 per bottom, shelf and top surfaces, with a maximum of 600 lbs. Color: Light Neutral. Size: 30"H x 16"D x Selected Width.



### 30" High Overhead Storage with Hinged Metal Doors

36" Overhead Storage	34"W x 16"D	VWR407520HM
42" Overhead Storage	40"W x 16"D	VWR407524HM
48" Overhead Storage	46"W x 16"D	VWR407521HM
54" Overhead Storage	52"W x 16"D	VWR407525HM
60" Overhead Storage	58"W x 16"D	VWR407522HM
66" Overhead Storage	64"W x 16"D	VWR407526HM
72" Overhead Storage	70"W x 16"D	VWR407523HM



### 30" High Overhead Storage with Hinged, Framed, Safety Glass Doors

36" Overhead Storage	34"W x 16"D	VWR407520HG
42" Overhead Storage	40"W x 16"D	VWR407524HG
48" Overhead Storage	46"W x 16"D	VWR407521HG
54" Overhead Storage	52"W x 16"D	VWR407525HG
60" Overhead Storage	58"W x 16"D	VWR407522HG
66" Overhead Storage	64"W x 16"D	VWR407526HG
72" Overhead Storage	70"W x 16"D	VWR407523HG



## Pass-Thru Overhead Storage Cabinets

Atlas overhead storage cabinets are used in place of open shelving when enclosed storage and/or resistance to dust is required. Mounting hardware included. Weight capacity: 40 lbs./ft.2 per bottom and top surfaces, with a maximum of 600 lbs. Color: Light Neutral. Size: 16<sup>3</sup>/<sub>8</sub>"H x 22<sup>1</sup>/<sub>8</sub>"D x Selected Width or 30"H x 22"D x Selected Width.





### 16<sup>3</sup>/<sub>8</sub>" High Pass-Thru Overhead Storage with Sliding Acrylic Doors

36" Pass-thru OH Storage	34"W x 22 <sup>1</sup> / <sub>8</sub> "D	VWR407530A
42" Pass-thru OH Storage	40"W x 22 <sup>1</sup> / <sub>8</sub> "D	VWR407534A
48" Pass-thru OH Storage	46"W x 22 <sup>1</sup> / <sub>8</sub> "D	VWR407531A
54" Pass-thru OH Storage	52"W x 22 <sup>1</sup> / <sub>8</sub> "D	VWR407535A
60" Pass-thru OH Storage	58"W x 22 <sup>1</sup> / <sub>8</sub> "D	VWR407532A
66" Pass-thru OH Storage	64"W x 22 <sup>1</sup> / <sub>8</sub> "D	VWR407536A
72" Pass-thru OH Storage	70"W x 22 <sup>1</sup> / <sub>8</sub> "D	VWR407533A

### 16<sup>3</sup>/<sub>8</sub>" High Pass-Thru Overhead Storage with Sliding Safety Glass Doors

36" Pass-thru OH Storage	34"W x 22 <sup>1</sup> / <sub>8</sub> "D	VWR407530G
42" Pass-thru OH Storage	40"W x 22 <sup>1</sup> / <sub>8</sub> "D	VWR407534G
48" Pass-thru OH Storage	46"W x 22 <sup>1</sup> / <sub>8</sub> "D	VWR407531G
54" Pass-thru OH Storage	52"W x 22 <sup>1</sup> / <sub>8</sub> "D	VWR407535G
60" Pass-thru OH Storage	58"W x 22 <sup>1</sup> / <sub>8</sub> "D	VWR407532G
66" Pass-thru OH Storage	64"W x 22 <sup>1</sup> / <sub>8</sub> "D	VWR407536G
72" Pass-thru OH Storage	70"W x 22 <sup>1</sup> / <sub>8</sub> "D	VWR407533G




### 16<sup>3</sup>/<sub>8</sub>" High Pass-Thru Overhead Storage with Sliding Metal Doors

36" Pass-thru OH Storage	34"W x 22 <sup>1</sup> / <sub>8</sub> "D	VWR407530M
42" Pass-thru OH Storage	40"W x 22 <sup>1</sup> / <sub>8</sub> "D	VWR407534M
48" Pass-thru OH Storage	46"W x 22 <sup>1</sup> / <sub>8</sub> "D	VWR407531M
54" Pass-thru OH Storage	52"W x 22 <sup>1</sup> / <sub>8</sub> "D	VWR407535M
60" Pass-thru OH Storage	58"W x 22 <sup>1</sup> / <sub>8</sub> "D	VWR407532M
66" Pass-thru OH Storage	64"W x 22 <sup>1</sup> / <sub>8</sub> "D	VWR407536M
72" Pass-thru OH Storage	70"W x 22 <sup>1</sup> / <sub>8</sub> "D	VWR407533M

### 30" High Pass-Thru Overhead Storage with Hinged Metal Doors

36" Pass-thru OH Storage	34"W x 22"D	VWR407530HM
42" Pass-thru OH Storage	40"W x 22"D	VWR407534HM
48" Pass-thru OH Storage	46"W x 22"D	VWR407531HM
54" Pass-thru OH Storage	52"W x 22"D	VWR407535HM
60" Pass-thru OH Storage	58"W x 22"D	VWR407532HM
66" Pass-thru OH Storage	64"W x 22"D	VWR407536HM
72" Pass-thru OH Storage	70"W x 22"D	VWR407533HM



### 30" High Pass-Thru Overhead Storage with Hinged, Framed, Safety Glass Doors


36" Pass-thru OH Storage	34"W x 22"D	VWR407530HG
42" Pass-thru OH Storage	40"W x 22"D	VWR407534HG
48" Pass-thru OH Storage	46"W x 22"D	VWR407531HG
54" Pass-thru OH Storage	52"W x 22"D	VWR407535HG
60" Pass-thru OH Storage	58"W x 22"D	VWR407532HG
66" Pass-thru OH Storage	64"W x 22"D	VWR407536HG
72" Pass-thru OH Storage	70"W x 22"D	VWR407533HG

# VWR® Slotted Atlas Shelving



## Field-Applied Shelf Lips


Shelf lips may be used on the front or back of Atlas shelves as needed. Select size based on shelf size. Mounting hardware not included. Color: Light Neutral.



36" Removable Shelf Lip	<b>VWRLIP36</b>
42" Removable Shelf Lip	<b>VWRLIP42</b>
48" Removable Shelf Lip	<b>VWRLIP48</b>
54" Removable Shelf Lip	<b>VWRLIP54</b>
60" Removable Shelf Lip	<b>VWRLIP60</b>
66" Removable Shelf Lip	<b>VWRLIP66</b>
72" Removable Shelf Lip	<b>VWRLIP72</b>

## Upright Extensions


Upright extensions are used to extend Atlas uprights to or through ceilings to provide a chase for electric, data or gas service lines. Each column extension includes mounting hardware and finish plate. Color: Light Neutral Column insert: Black. Size: 8"W x 1¾"D x Selected Height.



24" Upright Extension	<b>VWRCE24</b>
30" Upright Extension	<b>VWRCE30</b>
36" Upright Extension	<b>VWRCE36</b>
48" Upright Extension	<b>VWRCE48</b>
54" Upright Extension	<b>VWRCE54</b>
60" Upright Extension	<b>VWRCE60</b>
66" Upright Extension	<b>VWRCE66</b>
72" Upright Extension	<b>VWRCE72</b>


## Pipe Chase Assembly

Atlas pipe chase assemblies are used when service lines are too numerous to fit into a standard Atlas upright. Includes metal panels and brackets for two-sided assembly. Panels are removable for access. Order (1) upright when using with existing upright. Color: Light Neutral. Size: As selected.



## Hat Channel

Hat channels are used to extend the depth of standard Atlas uprights to accommodate gas service lines and surface mount fixtures by adding an additional 2" to the depth of a standard Atlas upright to total 3½". These hat channels are easily removed for access. All holes required must be field drilled. Mounting hardware included. Color: Black.



14" Hat Channel	<b>VWRHC14</b>
18" Hat Channel	<b>VWRHC18</b>

\*Turretts not included

## Counter-Mounted Pipe Chase Assemblies

	6" Wide	8" Wide	10" Wide	12" Wide
24" L	<b>VWRPC246</b>	<b>VWRPC248</b>	<b>VWRPC2410</b>	<b>VWRPC2412</b>
36" L	<b>VWRPC366</b>	<b>VWRPC368</b>	<b>VWRPC3610</b>	<b>VWRPC3612</b>
42" L	<b>VWRPC426</b>	<b>VWRPC428</b>	<b>VWRPC4210</b>	<b>VWRPC4212</b>
48" L	<b>VWRPC486</b>	<b>VWRPC488</b>	<b>VWRPC4810</b>	<b>VWRPC4812</b>
54" L	<b>VWRPC546</b>	<b>VWRPC548</b>	<b>VWRPC5410</b>	<b>VWRPC5412</b>
60" L	<b>VWRPC606</b>	<b>VWRPC608</b>	<b>VWRPC6010</b>	<b>VWRPC6012</b>
66" L	<b>VWRPC666</b>	<b>VWRPC668</b>	<b>VWRPC6610</b>	<b>VWRPC6612</b>
72" L	<b>VWRPC726</b>	<b>VWRPC728</b>	<b>VWRPC7210</b>	<b>VWRPC7212</b>

## Atlas Pipe Chase Ceiling Trim Kit

6" Deep Pipe Chase	<b>CFPC-VWRPC6TK</b>
8" Deep Pipe Chase	<b>CFPC-VWRPC8TK</b>
10" Deep Pipe Chase	<b>CFPC-VWRPC10TK</b>
12" Deep Pipe Chase	<b>CFPC-VWRPC12TK</b>


## Floor-Mounted Pipe Chase Assemblies

	6" Wide	8" Wide	10" Wide	12" Wide
84" H	<b>VWRPC846</b>	<b>VWRPC848</b>	<b>VWRPC8410</b>	<b>VWRPC8412</b>
96" H	<b>VWRPC966</b>	<b>VWRPC968</b>	<b>VWRPC9610</b>	<b>VWRPC9612</b>



### Electrical Hat Channels with (1) GFI Outlet


Electrical hat channels are used when duplex electric receptacles are required to be flush mounted on both sides of column uprights. Electrical hat channels can be removed easily for access on both sides. These units can be used back to back on the same column upright. Includes (1) GFI-rated receptacle. All outlets 20 AMP. Mounting hardware included. Wiring not included. Color: Black Size: 6½"W x 18"H x 2"D



14" Electrical Insert	VWRHC20146
18" Electrical Insert	VWRHC20188

### Power Bars


Power bar electrical assemblies are used with or in lieu of column mounted electrical outlets. Power Bars include 20 AMP grounded receptacles, an on/off switch with built-in circuit breaker and a 8-ft. power cord. Mounting hardware included. Color: Light Neutral



24" Power Bar with 6-20 Amp Outlets	VWRP20246
30" Power Bar with 6-20 Amp Outlets	VWRP20306
36" Power Bar with 8-20 Amp Outlets	VWRP20368
42" Power Bar with 8-20 Amp Outlets	VWRP20428
48" Power Bar with 8-20 Amp Outlets	VWRP20488
54" Power Bar with 8-20 Amp Outlets	VWRP20548
60" Power Bar with 8-20 Amp Outlets	VWRP20608
66" Power Bar with 8-20 Amp Outlets	VWRP20668
72" Power Bar with 8-20 Amp Outlets	VWRP20728

### Hat Channel Service Chases

Atlas hat channel service chases are used when service lines are too numerous to fit into a standard Atlas upright. Service chases are removable and fixtures can be added as necessary. Mounting hardware included. Color: Black. Size: 6½"W x 3"D x Selected Height




24" Hat Channel	VWRSHC24
36" Hat Channel	VWRSHC36
48" Hat Channel	VWRSHC48
54" Hat Channel	VWRSHC54
60" Hat Channel	VWRSHC60
66" Hat Channel	VWRSHC66
72" Hat Channel	VWRSHC72

### Flush Electrical Inserts with (1) GFI Outlet

Electrical column inserts are used when duplex electric receptacles are required to be flush mounted into column uprights. These inserts are easily removable for access from front. Bottom receptacle is GFI rated. Color: Black Size: 3⅝" W x 24"H x Selected Height.


24" Electrical Insert	VVRE20244
24" Electrical Insert	VVRE20248



### Data Inserts

Data column inserts are used when data connections are required to be flush mounted into column uprights. These inserts are easily removable for access from the front. Data outlets can be RJ-11 or RJ-45 configuration. Color: Black Size: 3⅝" W x 24"H


1-RJ45-CAT5 Outlets	VWRD45241
4-RJ45-CAT5 Outlets	VWRD45244
1-RJ11-CAT3 Outlets	VWRD11241
4-RJ11-CAT3 Outlets	VWRD11244





**Peg Boards**

Atlas epoxy resin, post-mounted peg boards are designed to be adjustable in height, allowing counters to remain clear of clutter. One-inch thick peg boards are sized to mount between standard 36" & 48" center-to-center dimensions of Atlas uprights and are mounted on (4) special brackets, purchased separately. Peg boards include finished backs. Color: Black. Size: As selected.




**Between-Column Tack Boards (Mounting Hardware Included)**

36" Tack Board	18"H x 34¼"W		VWRTP1836
36" Tack Board	24"H x 34¼"W		VWRTP2436
48" Tack Board	18"H x 46¼"W		VWRTP1848
48" Tack Board	24"H x 46¼"W		VWRTP2448
60" Tack Board	18"H x 58¼"W		VWRTP1860
60" Tack Board	24"H x 58¼"W		VWRTP2460
72" Tack Board	18"H x 70¼"W		VWRTP1872
72" Tack Board	24"H x 70¼"W		VWRTP2472

**Bin Bars**

Bin Bars provide another dimension for storage of small parts and pieces. Color: Light Neutral Size: As selected

36" Bin Bar	VWRBB36
42" Bin Bar	VWRBB42
48" Bin Bar	VWRBB48
54" Bin Bar	VWRBB54
60" Bin Bar	VWRBB60
72" Bin Bar	VWRBB72



**LabLink LED LabLight, Task Lighting**

The LabLink LabLight is a streamline LED task light for laboratory workbenches. The specially designed polycarbonate lens directs up to 95% of the illumination from white LED lamps to the work surface. This feature maximizes light dispersion for uniform task lighting in laboratories and similar applications. Environment is Indoor, Dry Locations.


PG LifeLink's LabLight is made from 95% recycled materials and Energy Star rated components. The aluminum housing efficiently manages heat from the high power LEDs and stays cool to the touch. The 50,000 hour LED service life makes the LabLight virtually maintenance free. Expected service life is 20 years @ 50 hours/week.

Power connections on both ends of the fixture allow the "daisy chaining" of up to 10' of light. The LabLight is available in seven models ranging in lengths from 15.125" to 87.125". Add 1.625" to length for each power connection. Length tolerance is ±0.0625".

LabLights are easily attached to metal, plastic or wood furniture using standard magnetic strip mounts or under shelf screw mounts.

C-UL Listed. ARRA Compliant. ANSI C3.

**Ordering Information:** Optional bracket mounts are also available at an additional charge.



Description	Dimensions		Cat. No.
<b>With 24" Daisy Chain</b>			
15" LED LabLight	3.5W x 2.2H x 38.4L cm (1⅜" x 7/8" x 15½")		PG-LL-15-D24
27" LED LabLight	3.5W x 2.2H x 68.9L cm (1⅜" x 7/8" x 27⅞")		PG-LL-27-D24
39" LED LabLight	3.5W x 2.2H x 99.4L cm (1⅜" x 7/8" x 39⅞")		PG-LL-39-D24
51" LED LabLight	3.5W x 2.2H x 129.9L cm (1⅜" x 7/8" x 51⅞")		PG-LL-51-D24
63" LED LabLight	3.5W x 2.2H x 160.3L cm (1⅜" x 7/8" x 63⅞")		PG-LL-63-D24
75" LED LabLight	3.5W x 2.2"H x 190.8L cm (1⅜" x 7/8" x 75⅞")		PG-LL-75-D24
87" LED LabLight	3.5W x 2.2"H x 221.3L cm (1⅜" x 7/8" x 87⅞")		PG-LL-87-D24
<b>With Power Supply</b>			
15" LED LabLight	3.5W x 2.2H x 38.4L cm (1⅜" x 7/8" x 15½")		89207-206
27" LED LabLight	3.5W x 2.2H x 68.9L cm (1⅜" x 7/8" x 27⅞")		89207-208
39" LED LabLight	3.5W x 2.2H x 99.4L cm (1⅜" x 7/8" x 39⅞")		89207-210
51" LED LabLight	3.5W x 2.2H x 129.9L cm (1⅜" x 7/8" x 51⅞")		89207-212
63" LED LabLight	3.5W x 2.2H x 160.3L cm (1⅜" x 7/8" x 63⅞")		89207-214
75" LED LabLight	3.5W x 2.2H x 190.8L cm (1⅜" x 7/8" x 75⅞")		89207-216
87" LED LabLight	3.5W x 2.2H x 221.3L cm (1⅜" x 7/8" x 87⅞")		89207-218
<b>With 12" Daisy Chain</b>			
15" LED LabLight	3.5W x 2.2H x 38.4L cm (1⅜" x 7/8" x 15½")		PG-LL-15-D12
27" LED LabLight	3.5W x 2.2H x 68.9L cm (1⅜" x 7/8" x 27⅞")		PG-LL-27-D12
39" LED LabLight	3.5W x 2.2H x 99.4L cm (1⅜" x 7/8" x 39⅞")		PG-LL-39-D12
51" LED LabLight	3.5W x 2.2H x 129.9L cm (1⅜" x 7/8" x 51⅞")		PG-LL-51-D12
63" LED LabLight	3.5W x 2.2H x 160.3L cm (1⅜" x 7/8" x 63⅞")		PG-LL-63-D12
75" LED LabLight	3.5W x 2.2H x 190.8L cm (1⅜" x 7/8" x 75⅞")		PG-LL-75-D12
87" LED LabLight	3.5W x 2.2H x 221.3L cm (1⅜" x 7/8" x 87⅞")		PG-LL-87-D12
<b>Accessories</b>			
12" Daisy Chain	—		PG-LL-DC12
24" Daisy Chain	—		PG-LL-DC24

### Large Keyboard Tray with Atlas System Mount, Innovative Office Products

Innovative's data-entry arm comprised of our best-selling 7000 monitor arm with a large keyboard tray. Features an articulating arm that suspends the keyboard and offers one finger movement. The large keyboard tray accommodates keyboard and mouse.


This ergonomic keyboard tray RoHS compliant.

<b>Environmentally Preferable Attribute</b>	Product and/or packaging contains recycled/renewable content
<b>Extension Length</b>	Arm extends up to 24"; collapses to 3" when not in use




<b>Description</b>	<b>Cat. No.</b>
Large Keyboard Tray	10147-002

#### Arm Capabilities



#### Keyboard Tray



### Ergonomic Compact Keyboard Tray, Innovative Office Products

Work in comfort and eliminate repetitive stress injuries with a flexible, ergonomic keyboard tray from Innovative Office Products.

This wall mount compact arm includes a 27" keyboard tray to support the keyboard and mouse.

No knob, lever, or ratchet handle required to set tray height, simply lift to desired position and release.

This ergonomic compact keyboard tray is RoHS compliant.



<b>Description</b>	<b>Cat. No.</b>
Ergonomic Compact Keyboard Tray, Black	10799-198

### Ergonomic Extended Reach Keyboard Tray, Innovative Office Products

Work in comfort and eliminate repetitive stress injuries with a flexible, ergonomic keyboard tray.

This wall mount extended reach arm includes a 27" keyboard tray to support the keyboard and mouse.

No knob, lever, or ratchet handle required to set tray height, simply lift to desired position and release.

This ergonomic extended reach keyboard tray is RoHS compliant.



<b>Description</b>	<b>Cat. No.</b>
Ergonomic Extended Reach Keyboard Tray, Black	10805-244

## Ergonomic Dual Monitor Arm with Atlas System Mount, Innovative Office Products

The adjustable wing bracket allows for convenient and ergonomic adjustment of the monitor arms while allowing rotation as well.


This dual monitor arm wall mount is RoHS compliant.




<b>Environmentally Preferable Attribute</b>	Product and/or packaging contains recycled/renewable content
<b>Horizontal Range</b>	28" of horizontal range
<b>Rotation</b>	Rotates 360 degrees at three joints
<b>Vertical Range</b>	18" of vertical range

Description	Capacity	Cat. No.
Dual Monitor Arm Wall Mount	3.5 - 13.5 lbs. per monitor	10147-008

### Arm Capabilities



### Wingback Capabilities



**Max monitor width\* .....23" (58.4 cm)\*\* per monitor**

\* Monitor width is measured from end to end of the monitor as shown here:  
 \*\*This measurement is based on a monitor depth of 2.5" (6.3 cm).  
 The depth of the monitor will have a slight effect on the overall max dimension.

## Ergonomic Laptop Arm with Atlas System Mount, Innovative Office Products

Our flexible laptop arm mount combines our best-selling 7000 monitor arm with our award-winning laptop bracket. Create an instant height-adjustable docking station for your laptop! You'll save valuable space and work in greater ergonomic comfort.


This laptop arm mount is RoHS compliant.





<b>Environmentally Preferable Attribute</b>	Product and/or packaging contains recycled/renewable content
<b>Horizontal Range</b>	24" of horizontal range
<b>Rotation</b>	Rotates 360 degrees at three joints
<b>Vertical Range</b>	16" of vertical range

Description	Capacity	Cat. No.
Flexible Laptop Arm Wall Mount	1 - 14 lbs.	10146-988

### Laptop Holder



### Arm



**For ultra thin laptops/tablet PCs:**  
 The supplied clips can be attached to the clasps to ensure a secure grasp.

## Flexible Data Entry Monitor Arm with Atlas System Mount and Flip-Up Keyboard, Innovative Office Products

This remarkably flexible data entry wall mount features a keyboard tray which can be folded flat when not in use and a slide-out mouse tray which accommodates left- or right-hand users. The wall mount monitor arm includes our spring-assist monitor tilter for effortless monitor motion.

This model is the right choice for space-sensitive environments where the computer is accessed by multiple users in a standing position. Medical and dental, retail and manufacturing environments will find this to be an ideal solution.


This flexible data entry mount is RoHS compliant.




<b>Environmentally Preferable Attribute</b>	Product and/or packaging contains recycled/renewable content
<b>Horizontal Range</b>	27" of horizontal range
<b>Rotation</b>	Rotates 360 degrees at three joints
<b>Vertical Range</b>	18" of vertical range

Description	Capacity	Cat. No.
Flexible Data Entry Arm Wall Mount with Flip-Up Keyboard	6 - 23 lbs.	10147-006

### Arm Capabilities



### Flip-up Keyboard Tray



## Flexible Ergonomic Monitor Arm with Atlas System Mount, Innovative Office Products

The best-selling monitor mount suspends your flat panel monitor, allowing you to quickly reposition with monitor with just one hand as you work in greater ergonomic comfort.


This monitor arm wall mount is RoHS compliant.



<b>Environmentally Preferable Attribute</b>	Product and/or packaging contains recycled/renewable content
<b>Horizontal Range</b>	24" of horizontal range
<b>Rotation</b>	Rotates 360 degrees at three joints
<b>Vertical Range</b>	16" of vertical range

Description	Capacity	Cat. No.
Flexible Monitor Arm Wall Mount	2.5 - 25 lbs.	10147-004

### Arm Capabilities



## Winston Sit/Stand Workstations, Innovative Office Products

Winston Workstation instantly converts any desk into a sit/stand workspace, instantly boosting health and productivity.

Features push-button height adjustment and the industry's first locking cylinder to prevent unwanted movement while standing.

The freestanding base retrofits to any existing desk, while the expansive work surface and storage tray allow user to bring all necessities with them when they stand.

This sit/stand workstation is VESA compatible and RoHS compliant.

**Ordering Information:** For optimal performance, select the model that will support that total weight of your monitor and all items that will be placed on the Winston work surface.



<b>Adjustable Height Range</b>	17" of vertical height adjustment
<b>Color</b>	Silver; Black Work Surface
<b>Environmentally Preferable Attribute</b>	Product and/or packaging contains recycled/renewable content
<b>Work Surface Dimensions</b>	30" wide

Description	Weight Capacity	Cat. No.
Winston Sit/Stand Workstation; One Monitor, Freestanding, Light Duty	9-20 lbs.	10753-160
Winston Sit/Stand Workstation; Dual Monitors, Freestanding, Medium Duty	12-28 lbs.	10753-174

---

## VWR® BENCH SELECTOR GUIDE

---

The VWR Bench Selector allows you to easily design a lab workstation complete with accessories. Follow four (4) easy steps - the configurator walks you through the simple process of selecting a VWR Bench and adding accessories to build a lab workstation that best meets your needs. It has never been easier to design and quote a lab workstation!




[vwr.com/benchselector](https://www.vwr.com/benchselector)

# VWR® Bench Selector Design Tool

The VWR Bench Selector allows users to easily design a lab workspace complete with accessories. Following four easy steps, the tool walks users through the simple process of selecting a VWR bench and adding accessories to build a lab workstation that best meets their needs. Users have the ability to choose from a variety of accessories, including shelving, lighting, etc, all while watching their design take form on the screen. The user can rotate and zoom at any point to better visualize their design.

Once complete, a VWR parts list and drawing will be sent to the user by email.


It has never been easy to design and quote a lab workstation!



## 1

### Laboratory Bench Type

The first step in designing a workstation using the VWR Bench Selector is selecting the bench type. Users can choose between a C-leg bench and a TP 4-leg bench. If a C-leg is chosen, the user must then decide between Standard, Motorized, or Hand-crank height adjustment options.




# 2

## Size & Materials

After selecting a bench type, the user moves on to Step 2: Size & Materials. In this step, the bench size is selected. TP 4-leg frame sizes include 30" x 48", 30" x 60" and 30" x 72". C-leg frames also have a 36" depth option.

The user then must select a worksurface type. Black laminate and phenolic options are available.

\*Please note, choosing between laminate and phenolic will not change the image shown on the right-hand side of the screen, as both are black in color. Once the parts list is provided, the list will reflect the choice in surface.




# 3

## Upright Modules

In Step 3, upright modules are selected. Uprights allow for attachment of accessories to the bench. Modules are available in single bay or double bay, with a double bay allowing for more room for accessories. Upright height is then chosen, with options including 42.3" or 51.2" heights.

Also in Step 3, accessories that are attached to the uprights are selected. These include shelves, bin rails, overhead cabinets and power bars. Lip edges and lights can be added to each shelf after it is added. Accessories can be added to six locations along the length of the upright module.

If a double upright module was chosen, there will be two accessories tabs for adding accessories to each side of the upright module.



# 3

## Upright Modules (continued)


A helpful tool has been added to Step 3 to assist users in adding accessories to their uprights. When the  button is clicked, a window opens that shows users where accessories can be added to the uprights.

\*Note: The  window shows accessories being added to six locations along the length of the upright. These accessories can be relocated anywhere along the upright once they have received the actual bench.



When adding bin rails, users have the option of adding stacking bins to the bin rail. These stacking bins come in four sizes and four colors. The Bench Selector only allows the user to add bins to the rail until the rail reaches its capacity.

Shown right: shelves, shelf lips, LED lights, and a bin rail with bins have been added to a bench.



# 4



## Accessories

Next, additional workstation accessories can added. These accessories include:

- Additional upright accessories such as monitor arms and bottle holders
- Casters
- Transport skates
- CPU holders
- Bag holders
- Footrests
- Keyboard trays
- Drawer units (including Contour™)

\*Drawer mounting brackets required will be included in the PDF part number list.


The VWR® Bench Selector makes adding accessories easy! The tool will not allow users to add accessories to a location where one accessory would conflict with the position of another accessory. Shown right, if two accessories conflict, the first accessory is highlighted in red and a warning window is displayed, giving the user the option to remove one of the accessories.



# 4

## Accessories (continued)

During any step of the configuration process, users can rotate the image by simply clicking and dragging the image. There is also a zoom feature on the bottom right corner of the screen.




## Send for Quotation

Congratulations! Your workstation is complete!

The final step in the process is to submit the design for quotation. After the workstation is complete, click “send for quotation” at the bottom right of the screen. This will open a window with a form. Complete the form with contact information, order quantity, and any additional notes, and click “send for quotation”. Once submitted, the user will receive an email with a PDF attachment that includes the image of their design, along with a details parts list.

The VWR® Bench Selector sends a copy of this email to **furniture@vwr.com**, where the quote request will be fulfilled.


\*Note: The image included in the quotation request PDF will appear exactly as shown when clicking “send for quotation”. Therefore the configured station should be positioned as desired before submitting the form.





## Send for Quotation (continued)

Shown here is an example of the "Send for Quotation" PDF.





**Sales Support:** 1-800-932-5000    furniture@vwr.com

**Inquiry form**

**Configuration reference number:** 198595044  
to discuss your options



**Itemized workbench**

---

**Frame:** C-Leg Bench  
**Model:** Motorized Concept  
**Size:** 30x72  
**Material:** Postformed front edge, laminate (black)  
**Uprights:** 2 x M36

Column 1	Column 2
Steel Shelf 36"x16"	
+Bookend	
+Shelf Lip	None
None	Overhead 'flipper door' cabinet
+Under Lighting	
Steel Shelf 36"x16"	
+Shelf Lip	
+Undershelf Lighting	None
None	None
Bin Rail	None
None	None

**Workbench summary**

Item	Product code	Quantity
C-Leg Bench Motorized	Motorized Concept	1
Postformed front edge, laminate (black)	89524-206	1
30x72	89524-342	1
<b>Uprights</b>		
Upright kit 2xM36	90435002	1
Steel Shelf 36"x16"	89524-498	2
Bookend	89524-506	1
Shelf Lip	89524-512	2
Undershelf Lighting	89524-522	1
Bin Rail	89524-570	1
Overhead 'flipper door' cabinet	89524-596	1
Under Lighting	89524-522	1
<b>Accessories</b>		
Bottle Holder	89524-616	1
LCD Holder Double Arm	89524-560	1
Rollset Casters	89524-482	1
Cpu Holder Barclamp	860551-35	1
Bar for Concept CPU/Footrest	-	1
Foot Rest Bar	89524-582	1
Keyboard tray	89524-556	1
File Drawer	89524-604	1
		x5

**VWR**

---

Thank you for contacting VWR and your interest in VWR Laboratory Benches. We will process your inquiry as soon as possible and will respond back to you with a quotation. (See attached file: quote.pdf)

# designed for your space

## Build Your VWR<sup>®</sup> Contour<sup>™</sup> Chair

Easily get the exact chair you need, specified by you and designed for exceptional performance and ergonomic comfort. Select the options that fit your specifications on the right, and then get a quote for the Contour line chair that's right for you.



Visit [vwr.com/chair-builder](http://vwr.com/chair-builder) to build your custom chair today!

1 Seat
2 Backrest
3 Base
4 Seat Height
5 Control
6 Adjustable Footing
7 Armrests
8 Performance Package
9 Material



**VWR** <sup>®</sup>  
COLLECTION

---

## VWR® LABORATORY BENCHES

---

VWR Laboratory Benches are ergonomically designed to meet the needs of laboratory personnel. There are many possible adjustments, all components fit together seamlessly, and each detail is designed for real world functionality.

The lab bench adapts to the work process and the needs of the user. As a result, all movements of the user can be optimized, which in turn increases performance and productivity and decreases work-related stress.

VWR Lab Benches are height-adjustable in manual, motorized and hand crank versions. All bench frames can be equipped with casters to create a mobile workstation solution. Frame color: Light Neutral. Available with laminate or black phenolic worksurfaces.

\*Benches ship unassembled. Assembled shipping and installation services available upon request. See pages 48-49 for product specs.

- Ergonomic and highly adaptable to the needs of the user
- Flexible, modular design
- Suitable for heavy loads up to 1100 lbs
- Offering a wide range of accessories to meet any need



# Mobile Laboratory Benches



Manual



Max Load: 1100 lbs.

Handcrank



Max Load: 440 lbs.

Motorized



Max Load: 660 lbs (dynamic)/  
1100 lbs. (static)

## VWR® C-Leg Lab Benches

Width	Worksurface	Height	Depth*	Cat. No.
<b>Manual Frame Style</b>				
48"	Phenolic	26.62"-44.34"	30"	89524-198
48"	Phenolic	26.62"-44.34"	36"	89524-832
60"	Phenolic	26.62"-44.34"	30"	89524-220
60"	Phenolic	26.62"-44.34"	36"	89524-834
72"	Phenolic	26.62"-44.34"	30"	89524-200
72"	Phenolic	26.62"-44.34"	36"	89524-836
48"	Laminate	26.62"-44.34"	30"	89524-202
48"	Laminate	26.62"-44.34"	36"	89524-838
60"	Laminate	26.62"-44.34"	30"	89524-204
60"	Laminate	26.62"-44.34"	36"	89524-840
72"	Laminate	26.62"-44.34"	30"	89524-206
72"	Laminate	26.62"-44.34"	36"	89524-842
<b>Handcrank Frame Style</b>				
48"	Phenolic	26.82"-40.60"	30"	89524-270
48"	Phenolic	26.82"-40.60"	36"	89524-904
60"	Phenolic	26.82"-40.60"	30"	89524-268
60"	Phenolic	26.82"-40.60"	36"	89524-906
72"	Phenolic	26.82"-40.60"	30"	89524-272
72"	Phenolic	26.82"-40.60"	36"	89524-906
48"	Laminate	26.82"-40.60"	30"	89524-192
48"	Laminate	26.82"-40.60"	36"	89524-826
60"	Laminate	26.82"-40.60"	30"	89524-194
60"	Laminate	26.82"-40.60"	36"	89524-828
72"	Laminate	26.82"-40.60"	30"	89524-196
72"	Laminate	26.82"-40.60"	36"	89524-830
<b>Motorized Frame Style</b>				
48"	Phenolic	27.80"-43.55"	30"	89524-334
48"	Phenolic	27.80"-43.55"	36"	89524-970
60"	Phenolic	27.80"-43.55"	30"	89524-356
60"	Phenolic	27.80"-43.55"	36"	89524-992
72"	Phenolic	27.80"-43.55"	30"	89524-336
72"	Phenolic	27.80"-43.55"	36"	89524-972
48"	Laminate	27.80"-43.55"	30"	89524-338
48"	Laminate	27.80"-43.55"	36"	89524-974
60"	Laminate	27.80"-43.55"	30"	89524-340
60"	Laminate	27.80"-43.55"	36"	89524-976
72"	Laminate	27.80"-43.55"	30"	89524-342
72"	Laminate	27.80"-43.55"	36"	89524-978

\*Depth with uprights attached are 32.625" for 30"D tables and 38.625" for 36"D tables, due to depth of uprights and notched table top.




**VWR® C-Leg Frame Only**

Width	Height	Depth*	Cat. No
<b>Manual Frame Style</b>			
48"	26.62"-44.34"	30"	10682-750
48"	26.62"-44.34"	36"	10682-780
60"	26.62"-44.34"	30"	10682-776
60"	26.62"-44.34"	36"	10682-782
72"	26.62"-44.34"	30"	10682-778
72"	26.62"-44.34"	36"	10682-784
<b>Handcrank Frame Style</b>			
48"	26.82"-40.60"	30"	10682-786
48"	26.82"-40.60"	36"	10682-752
60"	26.82"-40.60"	30"	10682-788
60"	26.82"-40.60"	36"	10682-754
72"	26.82"-40.60"	30"	10682-790
72"	26.82"-40.60"	36"	10682-756
<b>Motorized Frame Style</b>			
48"	27.80"-43.55"	30"	10682-758
48"	27.80"-43.55"	36"	10682-764
60"	27.80"-43.55"	30"	10682-760
60"	27.80"-43.55"	36"	10682-766
72"	27.80"-43.55"	30"	10682-762
72"	27.80"-43.55"	36"	10682-768

**VWR® C-Leg Manual Frame with Notched Worksurface**

Width	Worksurface	Depth	Cat. No.
<b>Single Bay</b>			
48"	Phenolic	32.625"	10682-568
48"	Phenolic	38.625"	10682-574
60"	Phenolic	32.625"	10682-570
60"	Phenolic	38.625"	10682-576
72"	Phenolic	32.625"	10682-572
72"	Phenolic	38.625"	10682-578
48"	Laminate	32.625"	10682-580
48"	Laminate	38.625"	10682-618
60"	Laminate	32.625"	10682-582
60"	Laminate	38.625"	10682-620
72"	Laminate	32.625"	10682-584
72"	Laminate	38.625"	10682-622
<b>Double Bay</b>			
60"	Phenolic	32.625"	10682-668
60"	Phenolic	38.625"	10682-672
72"	Phenolic	32.625"	10682-670
72"	Phenolic	38.625"	10682-674
60"	Laminate	32.625"	10682-676
60"	Laminate	38.625"	10682-680
72"	Laminate	32.625"	10682-678
72"	Laminate	38.625"	10682-682

Notched work surfaces add 2.625" to depth of table frame. Work surfaces for single bay uprigths include 2 notches (left and right); work surfaces for double bay uprigths include 3 notches (left, center, right)



Shown Above: Single notched top for single bay uprigths

## VWR® C-Leg Handcrank Frame with Notched Worksurface

Width	Worksurface	Depth	Cat. No.
<b>Single Bay</b>			
48"	Phenolic	32.625"	10682-624
48"	Phenolic	38.625"	10682-630
60"	Phenolic	32.625"	10682-626
60"	Phenolic	38.625"	10682-632
72"	Phenolic	32.625"	10682-628
72"	Phenolic	38.625"	10682-634
48"	Laminate	32.625"	10682-636
48"	Laminate	38.625"	10682-642
60"	Laminate	32.625"	10682-638
60"	Laminate	38.625"	10682-644
72"	Laminate	32.625"	10682-640
72"	Laminate	38.625"	10682-646
<b>Double Bay</b>			
60"	Phenolic	32.625"	10682-684
60"	Phenolic	38.625"	10682-688
72"	Phenolic	32.625"	10682-686
72"	Phenolic	38.625"	10682-690
60"	Laminate	32.625"	10682-692
60"	Laminate	38.625"	10682-696
72"	Laminate	32.625"	10682-694
72"	Laminate	38.625"	10682-698

Notched work surfaces add 2.625" to depth of table frame. Work surfaces for single bay uprights include 2 notches (left and right); work surfaces for double bay uprights include 3 notches (left, center, right)

## VWR® C-Leg Motorized Frame with Notched Worksurface

Width	Worksurface	Depth	Cat. No.
<b>Single Bay</b>			
48"	Phenolic	32.625"	10682-650
48"	Phenolic	38.625"	10682-656
60"	Phenolic	32.625"	10682-652
60"	Phenolic	38.625"	10682-658
72"	Phenolic	32.625"	10682-654
72"	Phenolic	38.625"	10682-660
48"	Laminate	32.625"	10682-770
48"	Laminate	38.625"	10682-662
60"	Laminate	32.625"	10682-772
60"	Laminate	38.625"	10682-664
72"	Laminate	32.625"	10682-774
72"	Laminate	38.625"	10682-666
<b>Double Bay</b>			
60"	Phenolic	32.625"	10682-684
60"	Phenolic	38.625"	10682-704
72"	Phenolic	32.625"	10682-702
72"	Phenolic	38.625"	10682-706
60"	Laminate	32.625"	10682-708
60"	Laminate	38.625"	10682-712
72"	Laminate	32.625"	10682-710
72"	Laminate	38.625"	10682-714

Notched work surfaces add 2.625" to depth of table frame. Work surfaces for single bay uprights include 2 notches (left and right); work surfaces for double bay uprights include 3 notches (left, center, right)



Shown Above: Double notched top for double bay uprights



Manual Frame with laminate top



Handcrank Frame with laminate top

**VWR® C-Leg Lab Benches with Single Bay Uprights**

Width	Worksurface	Depth*	Cat. No.
<b>Manual</b>			
48.5"	Phenolic	32.625"	89524-208
48.5"	Phenolic	38.625"	89524-844
60.375"	Phenolic	32.625"	89524-210
60.375"	Phenolic	38.625"	89524-846
72"	Phenolic	32.625"	89524-212
72"	Phenolic	38.625"	89524-848
48.5"	Laminate	32.625"	89524-214
48.5"	Laminate	38.625"	89524-850
60.675"	Laminate	32.625"	89524-216
60.675"	Laminate	38.625"	89524-852
72"	Laminate	32.625"	89524-218
72"	Laminate	38.625"	89524-854
<b>Handcrank</b>			
48.5"	Phenolic	32.625"	89524-274
48.5"	Phenolic	38.625"	89524-910
60.375"	Phenolic	32.625"	89524-276
60.375"	Phenolic	38.625"	89524-912
72"	Phenolic	32.625"	89524-278
72"	Phenolic	38.625"	89524-914
48.5"	Laminate	32.625"	89524-280
48.5"	Laminate	38.625"	89524-916
60.375"	Laminate	32.625"	89524-282
60.375"	Laminate	38.625"	89524-918
72"	Laminate	32.625"	89524-284
72"	Laminate	38.625"	89524-920
<b>Motorized</b>			
48.5"	Phenolic	32.625"	89524-344
48.5"	Phenolic	38.625"	89524-980
60.375"	Phenolic	32.625"	89524-346
60.375"	Phenolic	38.625"	89524-982
72"	Phenolic	32.625"	89524-348
72"	Phenolic	38.625"	89524-984
48.5"	Laminate	32.625"	89524-350
48.5"	Laminate	38.625"	89524-986
60.375"	Laminate	32.625"	89524-352
60.375"	Laminate	38.625"	89524-988
72"	Laminate	32.625"	89524-354
72"	Laminate	38.625"	89524-990

\*Depth with uprights attached are 32.625" for 30"D tables and 38.625" for 36"D tables, due to depth of uprights and notched table top. Width dimensions shown above also include added width from uprights.

**VWR® C-Leg Lab Benches with Double Bay Uprights**

Width	Worksurface	Depth*	Cat. No.
<b>Manual</b>			
60.375"	Phenolic	32.625"	89524-222
60.375"	Phenolic	38.625"	89524-856
72"	Phenolic	32.625"	89524-224
72"	Phenolic	38.625"	89524-858
60.675"	Laminate	32.625"	89524-226
60.675"	Laminate	38.625"	89524-860
72"	Laminate	32.625"	89524-228
72"	Laminate	38.625"	89524-862
<b>Handcrank</b>			
60.375"	Phenolic	32.625"	89524-286
60.375"	Phenolic	38.625"	89524-922
72"	Phenolic	32.625"	89524-288
72"	Phenolic	38.625"	89524-924
60.375"	Laminate	32.625"	89524-290
60.375"	Laminate	38.625"	89524-926
72"	Laminate	32.625"	89524-292
72"	Laminate	38.625"	89524-928
<b>Motorized</b>			
60.375"	Phenolic	32.625"	89524-358
60.375"	Phenolic	38.625"	89524-994
72"	Phenolic	32.625"	89524-360
72"	Phenolic	38.625"	89524-996
60.375"	Laminate	32.625"	89524-362
60.375"	Laminate	38.625"	89524-998
72"	Laminate	32.625"	89524-364
72"	Laminate	38.625"	89525-000

\*Depth with uprights attached are 32.625" for 30"D tables and 38.625" for 36"D tables, due to depth of uprights and notched table top. Width dimensions shown above also include added width from uprights.



11.81"D Steel Shelf with mounting brackets, can be positioned along the height of the upright.

Motorized Frame with phenolic top



\*Depth with uprights attached are 32.625" for 30"D tables and 38.625" for 36"D tables, due to depth of uprights and notched table top. Width dimensions shown above also include added width from uprights.

Motorized Frame with phenolic top

## VWR® C-Leg Lab Benches with Single Bay Uprights and Steel Shelf

Dimensions, W x H	Worksurface	Depth	Steel Shelf	Cat. No.
<b>Manual</b>				
48.5" x 53.5"	Phenolic	32.625"	48"W Shelf	89524-230
48.5" x 53.5"	Phenolic	38.625"	48"W Shelf	89524-864
60.375" x 53.5"	Phenolic	32.625"	60"W Shelf	89524-231
60.375" x 53.5"	Phenolic	38.625"	60"W Shelf	89524-866
72" x 53.5"	Phenolic	32.625"	72"W Shelf	89524-232
72" x 53.5"	Phenolic	38.625"	72"W Shelf	89524-868
48.5" x 53.5"	Laminate	32.625"	48"W Shelf	89524-234
48.5" x 53.5"	Laminate	38.625"	48"W Shelf	89524-870
60.675" x 53.5"	Laminate	32.625"	60"W Shelf	89524-236
60.675" x 53.5"	Laminate	38.625"	60"W Shelf	89524-872
72" x 53.5"	Laminate	32.625"	72"W Shelf	89524-238
72" x 53.5"	Laminate	38.625"	72"W Shelf	89524-874
<b>Handcrank</b>				
48.5" x 53.5"	Phenolic	32.625"	48"W Shelf	89524-294
48.5" x 53.5"	Phenolic	38.625"	48"W Shelf	89524-930
60.375" x 53.5"	Phenolic	32.625"	60"W Shelf	89524-296
60.375" x 53.5"	Phenolic	38.625"	60"W Shelf	89524-932
72" x 53.5"	Phenolic	32.625"	72"W Shelf	89524-298
72" x 53.5"	Phenolic	38.625"	72"W Shelf	89524-934
48.5" x 53.5"	Laminate	32.625"	48"W Shelf	89524-300
48.5" x 53.5"	Laminate	38.625"	48"W Shelf	89524-936
60.375" x 53.5"	Laminate	32.625"	60"W Shelf	89524-302
60.375" x 53.5"	Laminate	38.625"	60"W Shelf	89524-938
72" x 53.5"	Laminate	32.625"	72"W Shelf	89524-304
72" x 53.5"	Laminate	38.625"	72"W Shelf	89524-940
<b>Motorized</b>				
48.5" x 53.5"	Phenolic	32.625"	48"W Shelf	89524-366
48.5" x 53.5"	Phenolic	38.625"	48"W Shelf	89525-002
60.375" x 53.5"	Phenolic	32.625"	60"W Shelf	89524-368
60.375" x 53.5"	Phenolic	38.625"	60"W Shelf	89525-004
72" x 53.5"	Phenolic	32.625"	72"W Shelf	89524-370
72" x 53.5"	Phenolic	38.625"	72"W Shelf	89525-006
48.5" x 53.5"	Laminate	32.625"	48"W Shelf	89524-372
48.5" x 53.5"	Laminate	38.625"	48"W Shelf	89525-008
60.375" x 53.5"	Laminate	32.625"	60"W Shelf	89524-374
60.375" x 53.5"	Laminate	38.625"	60"W Shelf	89525-010
72" x 53.5"	Laminate	32.625"	72"W Shelf	89524-376
72" x 53.5"	Laminate	38.625"	72"W Shelf	89525-012

## VWR® C-Leg Lab Benches with Single Bay Uprights and Two Steel Shelves

Width	Worksurface	Height	Depth	Steel Shelf	Cat. No.
<b>Manual</b>					
48.5"	Phenolic	53.5"	32.625"	48"W Shelves	89524-240
48.5"	Phenolic	53.5"	38.625"	48"W Shelves	89524-876
60.375"	Phenolic	53.5"	32.625"	60"W Shelves	89524-242
60.375"	Phenolic	53.5"	38.625"	60"W Shelves	89524-878
72"	Phenolic	53.5"	32.625"	72"W Shelves	89524-244
72"	Phenolic	53.5"	38.625"	72"W Shelves	89524-880
48.5"	Laminate	53.5"	32.625"	48"W Shelves	89524-246
48.5"	Laminate	53.5"	38.625"	48"W Shelves	89524-882
60.675"	Laminate	53.5"	32.625"	60"W Shelves	89524-248
60.675"	Laminate	53.5"	38.625"	60"W Shelves	89524-884
72"	Laminate	53.5"	32.625"	72"W Shelves	89524-250
72"	Laminate	53.5"	38.625"	72"W Shelves	89524-886
<b>Handcrank</b>					
48.5"	Phenolic	53.5"	32.625"	48"W Shelves	89524-306
48.5"	Phenolic	53.5"	38.625"	48"W Shelves	89524-942
60.375"	Phenolic	53.5"	32.625"	60"W Shelves	89524-308
60.375"	Phenolic	53.5"	38.625"	60"W Shelves	89524-944
72"	Phenolic	53.5"	32.625"	72"W Shelves	89524-310
72"	Phenolic	53.5"	38.625"	72"W Shelves	89524-946
48.5"	Laminate	53.5"	32.625"	48"W Shelves	89524-312
48.5"	Laminate	53.5"	38.625"	48"W Shelves	89524-948
60.375"	Laminate	53.5"	32.625"	60"W Shelves	89524-314
60.375"	Laminate	53.5"	38.625"	60"W Shelves	89524-950
72"	Laminate	53.5"	32.625"	72"W Shelves	89524-316
72"	Laminate	53.5"	38.625"	72"W Shelves	89524-952
<b>Motorized</b>					
48.5"	Phenolic	53.5"	32.625"	48"W Shelves	89524-378
48.5"	Phenolic	53.5"	38.625"	48"W Shelves	89525-014
60.375"	Phenolic	53.5"	32.625"	60"W Shelves	89524-380
60.375"	Phenolic	53.5"	38.625"	60"W Shelves	89525-016
72"	Phenolic	53.5"	32.625"	72"W Shelves	89524-382
72"	Phenolic	53.5"	38.625"	72"W Shelves	89525-018
48.5"	Laminate	53.5"	32.625"	48"W Shelves	89524-384
48.5"	Laminate	53.5"	38.625"	48"W Shelves	89525-020
60.375"	Laminate	53.5"	32.625"	60"W Shelves	89524-386
60.375"	Laminate	53.5"	38.625"	60"W Shelves	89525-022
72"	Laminate	53.5"	32.625"	72"W Shelves	89524-388
72"	Laminate	53.5"	38.625"	72"W Shelves	89525-024

## VWR® C-Leg Lab Benches with Double Bay Uprights and Two Steel Shelves

Width	Worksurface	Depth	Steel Shelf	Cat. No.
<b>Manual</b>				
60.375"	Phenolic	32.625"	30"W Shelves	89524-252
60.375"	Phenolic	38.625"	30"W Shelves	89524-888
72"	Phenolic	32.625"	36"W Shelves	89524-254
72"	Phenolic	38.625"	36"W Shelves	89524-890
60.675"	Laminate	32.625"	30"W Shelves	89524-256
60.675"	Laminate	38.625"	30"W Shelves	89524-892
72"	Laminate	32.625"	36"W Shelves	89524-258
72"	Laminate	38.625"	36"W Shelves	89524-894
<b>Handcrank</b>				
60.375"	Phenolic	32.625"	30"W Shelves	89524-318
60.375"	Phenolic	38.625"	30"W Shelves	89524-954
72"	Phenolic	32.625"	36"W Shelves	89524-320
72"	Phenolic	38.625"	36"W Shelves	89524-956
60.375"	Laminate	32.625"	30"W Shelves	89524-322
60.375"	Laminate	38.625"	30"W Shelves	89524-958
72"	Laminate	32.625"	36"W Shelves	89524-324
72"	Laminate	38.625"	36"W Shelves	89524-960
<b>Motorized</b>				
60.375"	Phenolic	32.625"	30"W Shelves	89524-390
60.375"	Phenolic	38.625"	30"W Shelves	89525-026
72"	Phenolic	32.625"	36"W Shelves	89524-392
72"	Phenolic	38.625"	36"W Shelves	89525-028
60.375"	Laminate	32.625"	30"W Shelves	89524-394
60.375"	Laminate	38.625"	30"W Shelves	89525-030
72"	Laminate	32.625"	36"W Shelves	89524-396
72"	Laminate	38.625"	36"W Shelves	89525-032

\*Depth with uprights attached are 32.625" for 30"D tables and 38.625" for 36"D tables, due to depth of uprights and notched table top. Width dimensions shown above also include added width from uprights.



Motorized Frame with phenolic top

## VWR® C-Leg Lab Benches with Double Bay Uprights and Four Steel Shelves

Width	Worksurface	Depth	Steel Shelf	Cat. No.
<b>Manual</b>				
60.375"	Phenolic	32.625"	30"W Shelves	89524-260
60.375"	Phenolic	38.625"	30"W Shelves	89524-896
72"	Phenolic	32.625"	36"W Shelves	89524-262
72"	Phenolic	38.625"	36"W Shelves	89524-898
60.675"	Laminate	32.625"	30"W Shelves	89524-264
60.675"	Laminate	38.625"	30"W Shelves	89524-900
72"	Laminate	32.625"	36"W Shelves	89524-266
72"	Laminate	38.625"	36"W Shelves	89524-902
<b>Handcrank</b>				
60.375"	Phenolic	32.625"	30"W Shelves	89524-326
60.375"	Phenolic	38.625"	30"W Shelves	89524-962
72"	Phenolic	32.625"	36"W Shelves	89524-328
72"	Phenolic	38.625"	36"W Shelves	89524-964
60.375"	Laminate	32.625"	30"W Shelves	89524-330
60.375"	Laminate	38.625"	30"W Shelves	89524-966
72"	Laminate	32.625"	36"W Shelves	89524-332
72"	Laminate	38.625"	36"W Shelves	89524-968
<b>Motorized</b>				
60.375"	Phenolic	32.625"	30"W Shelves	89524-398
60.375"	Phenolic	38.625"	30"W Shelves	89525-034
72"	Phenolic	32.625"	36"W Shelves	89524-400
72"	Phenolic	38.625"	36"W Shelves	89525-036
60.375"	Laminate	32.625"	30"W Shelves	89524-402
60.375"	Laminate	38.625"	30"W Shelves	89525-038
72"	Laminate	32.625"	36"W Shelves	89524-404
72"	Laminate	38.625"	36"W Shelves	89525-040

\*Depth with uprights attached are 32.625" for 30"D tables and 38.625" for 36"D tables, due to depth of uprights and notched table top. Width dimensions shown above also include added width from uprights.



Motorized Frame with phenolic top

---

## VWR® FOUR-LEGGED BENCHES

---

The classically designed VWR® Basic Four-Legged Workbench is the ideal basic bench, the simplest and the most cost-effective option in our range of workbenches. VWR® Basic Four-Legged Workbenches are standard benches that are suitable for a broad range of tasks and environments when fitted with accessories: testing, packing, repairing, inspecting. VWR® Basic Four-Legged Workbench is guaranteed for regular use with loads up to 660 lbs.

Designed for use in today's production and assembly environments, this definitive modular system, developed over more than 35 years, provides a high quality working environment. Customize your workstation to your own requirements.

### Stepless height adjustment

The height of workbenches in the basic four-legged range can be steplessly adjusted from 25.6" to 35.4", which makes these benches suitable for work while sitting or standing.

- Steel structure
- A cost-effective option
- Easy to equip with a variety of accessories
- Stepless height adjustment with allen-key between 25.6"-35.4"
- Max load: 660 lbs.
- Frames are of epoxy powder coated steel in light neutral
- Upright tubes allow for attachment of a wide range of accessories
- Worksurfaces included (unless otherwise stated)
- See page 50 for product specifications





**VWR® Four-Legged Frame Lab Benches**

Size, D x W*	Worksurface	Cat. No.
30" x 48"	None	10747-874
30" x 60"	None	10747-876
30" x 72"	None	10747-878
30" x 48"	Phenolic	89524-406
30" x 60"	Phenolic	89524-408
30" x 72"	Phenolic	89524-410
30" x 48"	Laminate	89524-412
30" x 60"	Laminate	89524-414
30" x 72"	Laminate	89524-416

**VWR® Four-Legged Frame with Notched Worksurface**

Width	Worksurface	Cat. No.
<b>Single Bay</b>		
48"	Phenolic	10682-716
60"	Phenolic	10682-718
72"	Phenolic	10682-720
48"	Laminate	10682-722
60"	Laminate	10682-724
72"	Laminate	10682-726
<b>Double Bay</b>		
60"	Phenolic	10682-728
72"	Phenolic	10682-730
60"	Laminate	10682-732
72"	Laminate	10682-734



**VWR® Four-Legged Frame Lab Benches with Single Bay Uprights**

Size, D x W*	Worksurface	Cat. No.
32.625" x 48.5"	Phenolic	89524-418
32.625" x 60.375"	Phenolic	89524-420
32.625" x 72"	Phenolic	89524-422
32.625" x 48.5"	Laminate	89524-424
32.625" x 60.375"	Laminate	89524-426
32.625" x 72"	Laminate	89524-428

\*Depth with uprights attached are 32.625" for 30"D tables and 38.625" for 36"D tables, due to depth of uprights and notched table top. Width dimensions shown above also include added width from uprights.



**VWR® Four-Legged Frame Lab Benches with Double Bay Uprights**

Size, D x W*	Worksurface	Cat. No.
32.625" x 60.375"	Phenolic	89524-430
32.625" x 72"	Phenolic	89524-432
32.625" x 60.375"	Laminate	89524-434
32.625" x 72"	Laminate	89524-436

\*Depth with uprights attached are 32.625" for 30"D tables and 38.625" for 36"D tables, due to depth of uprights and notched table top. Width dimensions shown above also include added width from uprights.



**VWR® Four-Legged Frame Lab Benches with Single Bay Uprights and a Steel Shelf**

Size, D x W*	Worksurface	Steel Shelf	Cat. No.
32.625" x 48.5"	Phenolic	48"W Shelf	89524-438
32.625" x 60.375"	Phenolic	60"W Shelf	89524-440
32.625" x 72"	Phenolic	72"W Shelf	89524-442
32.625" x 48.5"	Laminate	48"W Shelf	89524-444
32.625" x 60.375"	Laminate	60"W Shelf	89524-446
32.625" x 72"	Laminate	72"W Shelf	89524-448

\*Depth with uprights attached are 32.625" for 30"D tables and 38.625" for 36"D tables, due to depth of uprights and notched table top. Width dimensions shown above also include added width from uprights.



**VWR® Four-Legged Frame Lab Benches with Single Bay Uprights and Two Steel Shelves**

Size, D x W*	Worksurface	Steel Shelves	Cat. No.
32.625" x 48.5"	Phenolic	48"W Shelves	89524-450
32.625" x 60.375"	Phenolic	60"W Shelves	89524-452
32.625" x 72"	Phenolic	72"W Shelves	89524-454
32.625" x 48.5"	Laminate	48"W Shelves	89524-456
32.625" x 60.675"	Laminate	60"W Shelves	89524-458
32.625" x 72"	Laminate	72"W Shelves	89524-460

\*Depth with uprights attached are 32.625" for 30"D tables and 38.625" for 36"D tables, due to depth of uprights and notched table top. Width dimensions shown above also include added width from uprights.



**VWR® Four-Legged Frame Lab Benches with Double Bay Uprights and Two Steel Shelves**

Size, D x W*	Worksurface	Steel Shelves	Cat. No.
32.625" x 60.375"	Phenolic	30"W Shelves	89524-462
32.625" x 72"	Phenolic	36"W Shelves	89524-464
32.625" x 60.675"	Laminate	30"W Shelves	89524-466
32.625" x 72"	Laminate	36"W Shelves	89524-468


\*Depth with uprights attached are 32.625" for 30"D tables and 38.625" for 36"D tables, due to depth of uprights and notched table top. Width dimensions shown above also include added width from uprights.



**VWR® Four-Legged Frame Lab Benches with Double Bay Uprights and Four Steel Shelves**

Size, D x W*	Worksurface	Steel Shelves	Cat. No.
32.625" x 60.375"	Phenolic	30"W Shelves	89524-470
32.625" x 72"	Phenolic	36"W Shelves	89524-472
32.625" x 60.675"	Laminate	30"W Shelves	89524-474
32.625" x 72"	Laminate	36"W Shelves	89524-476

\*Depth with uprights attached are 32.625" for 30"D tables and 38.625" for 36"D tables, due to depth of uprights and notched table top. Width dimensions shown above also include added width from uprights.



### Upright Component Parts

Description	Cat. No.
Single Slotted Upright, 53"	10486-496
Single Slotted Upright, 63"	10486-498
30"W Horizontal Tube	10486-500
36"W Horizontal Tube	10486-502
48"W Horizontal Tube	10486-504
60"W Horizontal Tube	10486-506
72"W Horizontal Tube	10486-508
Corner Bracket Pair for Upright Modules	10486-510
Straight/Tee Bracket for Upright Modules	10486-512

### Upright Modules

Description	Cat. No.
<b>Single Bay</b>	
48" x 53.5"H	10496-572
48" x 62.3"H	10486-478
60" x 53.5"H	10486-480
60" x 62.3"H	10486-482
72" x 53.5"H	10486-484
72" x 62.3"H	10486-486
<b>Double Bay</b>	
60" x 53.5"H	10486-488
60" x 62.3"H	10486-490
72" x 53.5"H	10486-492
72" x 62.3"H	10486-494

Usable height of uprights = 42.3"/51.18" from top

### Casters

Description	Cat. No.
(1) Caster Set* 3.93" Diameter (adds 5.5")	89524-478
(2) Low Profile Caster Set*	89524-480
(3) Roll & Set Caster Set - 2.75"H (adds 3.5")	89524-482
(4) Transport Skate **	89524-484

\* Set includes two locking and two non-locking casters.

\*\* Adds 1.5" to table height when in use.



#### Transport Skate

With a transport skate, one person can easily move benches. For long-distance moving, you will need two transport skates.



Steel Shelf w/standard bracket (included)



## Shelf Lips

Description	Cat. No.
30" Shelf Lip	89524-510
36" Shelf Lip	89524-512
48" Shelf Lip	89524-514
60" Shelf Lip	89524-516
72" Shelf Lip	89524-518

- Lip height: 1.5"
- Lip height above shelf: .32"



## Bookend Brackets

Description	Cat. No.
11.81" Bookend Bracket	89524-506
15.74" Bookend Bracket	89524-508

## Shelves

Description	Shelf Depth	Max lbs	Cat. No.
30"	11.81"	176	89524-486
36"	11.81"	176	89524-488
48"	11.81"	176	89524-490
60"	11.81"	176	89524-492
72"	11.81"	176	89524-494
30"	15.74"	176	89524-496
36"	15.74"	176	89524-498
48"	15.74"	176	89524-500
60"	15.74"	176	89524-502
72"	15.74"	176	89524-504



Vertical Power Rail



Horizontal Power Rail

## Power Rails

Description	Cat. No.
30" Bin Rail w/6 Outlet Power	10611-600
30" Bin Rail w/6 Outlet Power, Surge	10611-602
36" Bin Rail w/6 Outlet Power	10611-604
36" Bin Rail w/6 Outlet Power, Surge	10611-606
48" Bin Rail w/6 Outlet Power	89524-530
48" Bin Rail w/6 Outlet Power, Surge	89524-532
60" Bin Rail w/6 Outlet Power	89524-534
60" Bin Rail w/6 Outlet Power, Surge	89524-536
72" Bin Rail w/6 Outlet Power	89524-538
72" Bin Rail w/6 Outlet Power, Surge	89524-540
6 Outlet Vertical Power Bar with Clamp*	89524-542
6 Outlet Vertical Power Bar with Clamp, Surge*	89524-544
12 Outlet Vertical Power Bar with Clamp*	89524-546
12 Outlet Vertical Power Bar with Clamp, Surge*	89524-548

\*Can also be mounted underneath the worksurface

Note: All VWR bench accessories with painted steel components match light neutral frame color.



Equipment Holder - C-leg (bar)



A) LCD holder w/ double arm



B/C) Keyboard tray for LCD holder\*



D) Height adjustable LCD monitor arm



Equipment Holder - Four-leg



Equipment Holder - C-leg



E) Bag Holder



F) Bottle Holder



G) Roll Holder

## Holders

Description	Size	Max lbs	Cat. No.
Equipment Holder - C-leg Bar Mount	3.94" x 19.68" x 3.35"	44	89524-550
Equipment Holder - C-leg Leg Mount	3.94"-9.49" x 17.71" x 1.96"	44	89524-552
Equipment Holder - C-leg Motor Leg Mount	3.94"-9.49" x 17.71" x 1.96"	44	89524-554
Equipment Holder - Four-leg	3.94"-9.49" x 17.71" x 1.96"	44	10746-456
LCD Holder with Double Arm, VESA <sup>75/100</sup> (A)	—	33	89524-560
Keyboard Tray for LCD Holder (B)	18.89" x 10.23"	11	89524-562
Keyboard and Mouse Tray for LCD Holder (C)	24.40" x 10.23"	11	89524-564
Height Adjustable Arm with Bracket (D)	—	33	89524-566
Bag Holder for Bench Top Lower Edge (E)	—	—	89524-614
Bottle Holder (3.14" base) (F)	—	—	89524-616
Roll Holder for Upright (Roll ø Height 5.90 inches, Width 6.88 inches) (G)	—	—	89524-618

\*Keyboard trays B and C compatible for LCD holder A only.

## Bin Rails

Description	Size	Cat. No.
30" Bin Rail	29.21" x 3.54"	89524-568
36" Bin Rail	35.23" x 3.54"	89524-570
48" Bin Rail	47.04" x 3.54"	89524-572
60" Bin Rail	58.89" x 3.54"	89524-574
72" Bin Rail	70.66" x 3.54"	89524-576

Bin rail load capacity 110 lbs per rail  
Bin rails may be mounted straight or angled at 15°



Bin rail with brackets. Bins sold separately.

## Accessories, Overhead & Suspended Cabinets



### Footrests

Description	Size	Cat. No.
Footrest, Floor (5-30° tilt)	21.65" x 13.77" x 8.66"	89524-578
Footrest, Floor Height Adjustable (5-30° tilt)	21.65" x 13.77" x 3.54"-11.02"	89524-580
Footrest, Bar Attachment (5-30° tilt)	21.65" x 13.77" x 3.54"-11.02"	89524-582
60" Bin Rail	58.89" x 3.54"	89524-574
72" Bin Rail	70.66" x 3.54"	89524-576



Sliding Door Cabinet



Binder Box



Overhead Flipper Door Cabinet

### Cabinets

Description	Size	Max lbs	Cat. No.
30" Sliding Door Cabinet	28.34" x 14.56" x 15.74"	88	89524-594
Divider for Sliding Door Cabinet	9.84" x 7.87"	—	89524-590
CD/Folder Box, Attached to Upright Tube	5.11" x 11.41" x 12.99"	11	89524-592
30" Overhead "Flipper Door" Cabinet	29.25" x 13" x 15.68"	88	89524-606
36" Overhead "Flipper Door" Cabinet	35.37" x 13" x 15.68"	88	89524-596
Divider for Overhead "Flipper Door" Cabinet	11.80" x 10.40"	—	89524-598



### Bench Drawers

Description	Size	Max lbs	Cat. No.
Bench Drawer Unit 30/15 w/ 1 - 3.94" Drawer	11.81" x 17.20" x 5.91"	22 lbs/drawer	89524-600
Bench Drawer Unit 30/22 w/ 2 - 3.94" Drawers	11.81" x 17.20" x 8.74"	22 lbs/drawer	89524-602
Bench Drawer Unit 38 w/ 2 - 100% Extendable Drawers	14.96" x 20.47" x 10.23"	132 lbs/drawer	89524-604

### Bench Drawer Attachment Brackets

Description	Cat. No.
Drawer Attachment Bracket for 30"D x 48"W Manual C-leg Bench	10820-620
Drawer Attachment Bracket for 30"D x 60"W Manual C-leg Bench	10820-622
Drawer Attachment Bracket for 30"D x 72"W Manual C-leg Bench	10820-624
Drawer Attachment Bracket for 30"D Handcrank C-leg Bench	89524-608
Drawer Attachment Bracket for 36"D Manual/Handcrank C-leg Bench	10746-454
Drawer Attachment Bracket for 30"/36"D Motorized C-leg Bench	89524-610
Drawer Mounting Bracket for 30" x 48" Four-leg Bench	10611-638
Drawer Mounting Bracket for 30" x 60" Four-leg Bench	10611-640
Drawer Mounting Bracket for 30" x 72" Four-leg Bench	10611-642
Drawer Mounting Bracket for Mounting a Second Drawer Unit to Four-leg Bench	10769-372

### Contour Cabinet Attachment Brackets


Description	Cat. No.
Bracket for C-leg Manual 30" x 48", 36" x 48"	10611-644
Bracket for C-leg Manual 30" x 60", 36" x 60"	10611-646
Bracket for C-leg Manual 30" x 72", 36" x 72"	10611-648
Bracket for C-leg Motor 30" x 48", 36" x 48"	10611-650
Bracket for C-leg Motor 30" x 60", 36" x 60"	10611-652
Bracket for C-leg Motor 30" x 72", 36" x 72"	10611-654
Drawer Mounting Bracket for 30" x 48" Four-leg Bench	10611-638
Drawer Mounting Bracket for 30" x 60" Four-leg Bench	10611-640
Drawer Mounting Bracket for 30" x 72" Four-leg Bench	10611-642
Drawer Mounting Bracket for Mounting a Second Contour Cabinet to Four-leg Bench	10769-372

Contour™ cabinets are not able to attach to C-leg Handcrank frames.


# VWR® Contour™ Suspended Cabinets




All suspended height cabinets are 19" or 25" high, 22" deep and 18" wide. Cabinets are designed to hang beneath table frames. Cabinets are phosphate-coated with a baked, chemical-resistant, powder-coat finish. Color: light neutral.



## 19" Drawer Cabinets



Three Drawer w/ Pull Board  
18" W CFD-5002-B




One Drawer w/ Pull Board  
One File Drawer w/ Hangers  
18" W CFD-5852-B




Swings Right  
18" W CFC-5032-A  
Swing Right w/ Pull Board  
18" W CFC-5032-B

Swings Left  
18" W CFC-5032-L  
Swing Left w/ Pull Board  
18" W CFC-5032-BL


## 25" Drawer Cabinets



Three Drawers  
18" W CFD-4002-A  
Three Drawers w/ Pull Board  
18" W CFD-4002-B




Two File Drawers with Hangers  
18" W CFD-4812-A




Three Drawers w/ Pull Board  
One File Drawer w/ Hangers  
18" W CFD-4852-B

## 25" Cupboard Cabinets




Swings Right  
18" W CFC-4032-A  
Swing Left  
18" W CFC-4032-L




Swings Right w/ Pull Board  
18" W CFC-4032-B  
Swings Left w/ Pull Board  
18" W CFC-4032-BL

## 25" Cupboard and Drawer Cabinets



Swings Right  
18" W CFE-4402-A  
Swings Left  
18" W CFE-4402-L



Swings Right w/ Pull Board  
18" W CFE-4402-B  
Swings Left w/ Pull Board  
18" W CFE-4402-BL

## Worktable Frames

- VWR® Worktable Frame consist of a pair of legs, rear case, and front rail.
- VWR® Frames are available in depths: 30" and 36"
- VWR® Frames are available in widths: 48", 60" and 72"
- VWR® Table Frames can be easily equipped with any 12mm stem caster.

## Steel Thickness:

Steels used in the construction of the frames are 12-, 14-, and 16-gauge CRS steel, except as follows:


- Leg attachment parts: 7-gauge CRS steel
- Cable rail: 20-gauge CRS steel

## Tools needed for C-leg frame assembly:

- 19mm wrench
- 13mm wrench
- 5mm Allen wrench
- #2 Phillips head screwdriver

## C-Leg Manual

- Leg pair can support 1100 lb load, evenly loaded.
- Height adjustment of VWR® Worktable Frame is 26.62" - 44.34" with table top (¾") and adjustment glide. Height adjustment is made by adjusting telescoping legs and tightening friction screws on the back and sides of the legs.



## Tools needed for Four-leg frame assembly:

- 10mm wrench
- 8mm Allen wrench

## Recommended table top assembly:

- VWR® Laminate Tops assembled with wood screws provided
- VWR® Phenolic Tops can be assembled in 2 ways:
  - Mechanical fasteners (screws) - Pre-drill a pilot hole in the bottom of the work surface, taking care to keep from penetrating through the top. The pre-drilled hole should be the same length as the screw shaft. The screw should be slightly larger in diameter than the pilot hole. Use a coarse threaded screw, wood or machine.
  - Adhesive - Use 100% silicon. Apply several "spots" or a continuous bead to the top of the table frame. Position the work surface as centered on the frame and let cure for 12 hours at room temperature.


	30" x 48"	36" x 48"	30" x 60"	36" x 60"	30" x 72"	36" x 72"
A	47.24"	47.24"	59.06"	59.06"	70.86"	70.86"
B	29.07"	35.09"	29.07"	35.09"	29.07"	35.09"
C	40.78"	40.78"	52.60"	52.60"	64.41"	64.41"
D	40.11"	40.11"	51.93"	51.93"	63.74"	63.74"
E	22.93"	28.88"	22.93"	28.88"	22.93"	28.88"
F	11.43"	11.43"	11.43"	11.43"	11.43"	11.43"
G	29.52"	29.52"	29.52"	29.52"	29.52"	29.52"
H	17.57"	17.57"	17.57"	17.57"	17.57"	17.57"
I	11.12"	11.12"	11.12"	11.12"	11.12"	11.12"
J	26.82-44.34"	26.82-44.34"	26.82-44.34"	26.82-44.34"	26.82-44.34"	26.82-44.34"
K	5.5"	5.5"	5.5"	5.5"	5.5"	5.5"
*K1	2.87"	2.87"	2.87"	2.87"	2.87"	2.87"

\*K1 - low profile caster option

### C-Leg Handcrank

- Pair of legs of crank frame can support 440 lbs dynamic load, and 1,100 lbs static load, evenly loaded.
- VWR® Frames are available in depths: 30" and 36"
- VWR® Frames are available in widths: 48", 60" and 72"
- Height adjustment of crank frame is 27.4"-41.2" with table top (3/4").

- The crank handle can be installed on the left or right side of the frame. Standard table configuration is right-hand crank.
- Height adjustment change one inch per each full rotation of the crank handle.
- The crank folds away under the frame when not in use.




	30" x 48"	36" x 48"	30" x 60"	36" x 60"	30" x 72"	36" x 72"
A	47.24"	47.24"	59.06"	59.06"	70.86"	70.86"
B	29.07"	35.09"	29.07"	35.09"	29.07"	35.09"
C	40.78"	40.78"	52.60"	52.60"	64.41"	64.41"
D	40.11"	40.11"	51.93"	51.93"	63.74"	63.74"
E	16.44"	16.44"	16.44"	16.44"	16.44"	16.44"
F	11.43"	11.43"	11.43"	11.43"	11.43"	11.43"
G	29.52"	29.52"	29.52"	29.52"	29.52"	29.52"
H	18.28"	18.28"	18.28"	18.28"	18.28"	18.28"
I	9.26"	9.26"	9.26"	9.26"	9.26"	9.26"
J	26.82-40.60"	26.82-40.60"	26.82-40.60"	26.82-40.60"	26.82-40.60"	26.82-40.60"
K	5.5"	5.5"	5.5"	5.5"	5.5"	5.5"
*K <sub>1</sub>	2.87"	2.87"	2.87"	2.87"	2.87"	2.87"
L	8.06"	8.06"	8.06"	8.06"	8.06"	8.06"

\*K<sub>1</sub> - low profile caster option

### C-Leg Motorized

- Leg pair of VWR® Motorized Frame can support a static load of 1,100 lbs, evenly loaded on a worktop and a dynamic load of 660 lbs, evenly loaded.
- VWR® Frames are available in depths: 30" and 36"
- VWR® Frames are available in widths: 48", 60" and 72"
- Height adjustment range of VWR® Motorized Frame is 27.80" - 43.55" with table top (3/4") and adjustment glides. Adjustment is made with 2 linear actuators, using an up/down memory desk switch.
- The speed of the frame is .70" / sec without load.

- The speed of the frame is .59" / sec with load.
- The desk switch can be installed on the front, left or right side of frame.
- The desk switch has multi-position memory system to be able to set the frame to three customer based heights and up-down buttons for random settings.
- All units have soft-start and soft-stop circuitry.
- The table can be manually driven over the memory positions.
- Silent operation (45 db (A) DS/EN ISO 3746
- Capable of operations in ambient temperature 41 to 104°F.



	30" x 48"	36" x 48"	30" x 60"	36" x 60"	30" x 72"	36" x 72"
A	47.24"	47.24"	59.06"	59.06"	70.86"	70.86"
B	29.07"	34.57"	29.07"	34.57"	29.07"	34.57"
C	38.12"	38.12"	49.94"	49.94"	61.75"	61.75"
D	40.11"	40.11"	51.93"	51.93"	63.74"	63.74"
E	14.88"	14.88"	14.88"	14.88"	14.88"	14.88"
F	11.43"	11.43"	11.43"	11.43"	11.43"	11.43"
G	29.52"	29.52"	29.52"	29.52"	29.52"	29.52"
H	18.28"	18.28"	18.28"	18.28"	18.28"	18.28"
I	9.26"	9.26"	9.26"	9.26"	9.26"	9.26"
J	31"	31"	43.53"	43.53"	55.33"	55.33"
K	5.5"	5.5"	5.5"	5.5"	5.5"	5.5"
*K <sub>1</sub>	2.87"	2.87"	2.87"	2.87"	2.87"	2.87"

\*K<sub>1</sub> - low profile caster option


## TP Four-Leg Bench

- VWR® TP Worktable Frame consists of a pair of legs, rear case, and front rail.
- VWR® TP Frames are available in 30" depth.
- VWR® TP Frames are available in widths: 48", 60" and 72"
- Table Frames can be easily equipped with any 12mm stem caster.
- Weight Capacity: 660 lbs.

### Steel Thickness:

Steel used in the construction of the frames are 12,- 14- and 16-gauge CRS steel, except as follows:

- Leg attachment parts: 7-gauge CRS steel
- Cable rail: 20-gauge CRS steel




	30" x 48"	30" x 60"	30" x 72"
A	47.2"	59.1"	70.9"
B	27.6"	27.6"	27.6"
C	25.6-35.4"	25.6-35.4"	25.6-35.4"
D	24.4"	24.4"	24.4"

## Slotted Uprights


- Uprights are available in heights: 53.5" and 62.3" (usable height above the work surface is 42.3" and 51.18")
- Depth of the uprights is 2.36".
- Width of the uprights is 1.18".
- Uprights can be attached to each table frame to enable the use of accessories
- Upright module consist of a minimum of 2 uprights, which are installed with the correct spacing for the desired module size (see table).
- Upright modules can easily be extended to create double modules (3 uprights needed).
- Each upright can support loads equal to the maximum load limited by table frame or accessories in use.
- For accessories the height adjustment steps are 1.1".
- Steel Thickness: Steel used in the construction of the upright is 14-gauge CRS steel in light neutral.

	48"	60"	72"	2" x 30"	2" x 36"
A	2.36"	2.36"	2.36"	2.36"	2.36"
B	1.18"	1.18"	1.18"	1.18"	1.18"
C	46.06"	57.87"	69.92"	28.35"	34.37"
D	47.2"	59.1"	71.1"	19.7"	19.7"
E	48.42"	60.23"	72.28"	20.86"	20.86"
F	—	—	—	59.06"	71.11"
G	—	—	—	60.23"	72.28"

Note: Use of uprights adds an additional 2.625" to depth of frame.



## Reagent Rack Hole Mounting Info




# Alpha Overhead Service Carriers

## Alpha Overhead Service Carriers

Designed with VWR Furniture Supplies technological expertise, the Alpha Overhead Service Carrier meets the needs of the modern laboratory. This flexible and adaptable system provides integrated plumbing, power, data, and accessories, and can be mounted overhead, along a wall, or under cabinets.

**Contact your VWR Furniture Specialist for assistance in configuring a service carrier to meet your requirements.**



## Alpha Overhead Service Enclosures

These enclosures can integrate multiple plumbing, power, lighting, data, and localized exhaust elements, and feature a suspended carrier design. Available in custom lengths and configurations to accommodate all overhead service requirements, Alpha Overhead Service Enclosures can be pre-plumbed and pre-wired for easy single point connection at the job site.

- Pre-plumbed lab services, including high purity
- Pre-wired power and data ports



## Alpha Ceiling Panels

Panels are designed to fit into a ceiling grid to economically deliver plumbing, power, and data to today's lab bench or equipment.

- Integrated in-ceiling grid
- Plumbing fittings
- Power outlets
- Data ports

**Keyed Quick Connect Body, WaterSaver**

The quick connect keying system is based upon the inner and outer diameters of the sleeves on the body and plug. Quick connect bodies and plugs color coded to designate the service and to facilitate matching the correct body and plug.

CSA certified for natural gas.



**Keyed Quick Connect Plug, WaterSaver**

Plugs have an internal shut-off valve that closes when the valve is disengaged. Internal valve components are stainless steel to assure compatibility with service to assure durable performance.

CSA certified for natural gas.



Description	Color	Cat. No.
Flanged Keyed Quick Connect	Black	89176-950
Flanged Keyed Quick Connect	Blue	89176-958
Flanged Keyed Quick Connect	Brown	89176-964
Flanged Keyed Quick Connect	Green	89176-954
Flanged Keyed Quick Connect	Orange	89176-952
Flanged Keyed Quick Connect	Purple	89176-962
Flanged Keyed Quick Connect	White	89176-960
Flanged Keyed Quick Connect	Yellow	89176-956

Description	Color	Cat. No.
Quick Connect Plug	Black	89176-966
Quick Connect Plug	Blue	89176-974
Quick Connect Plug	Brown	89176-980
Quick Connect Plug	Green	89176-970
Quick Connect Plug	Orange	89176-968
Quick Connect Plug	Purple	89176-978
Quick Connect Plug	White	89176-976
Quick Connect Plug	Yellow	89176-972

**LabLink Surface Mount Raceway, Versa-Duct, PG LifeLink**

- Extruded from #6063-T5 aluminum alloy
- Single channel assembly (1 3/4" x 2 3/4" x 48"-120")
- Double channel assembly (2 1/4" x 5 1/4" x 48"-120")
- Single circuit configuration
- 12" pigtails are provided for electrical connections
- Clear satin anodized finish is durable and corrosion resistant
- 12" removable cover plates for easy access in the field
- 20 amp (NEMA-5-20R) grey duplex receptacles
- Grey communication holder with 2 Cat 5e data jacks (not wired)
- C-UL Listed



LabLink Surface Mount Raceway (Versa-Duct), made from anodized extruded aluminum, offers a superior design to combine electrical and/or data into one unified system. Factory pre-cut and pre-wired to standard lengths. Single channel raceway is configured to include power receptacles and dual channel raceway will include power and data. The raceway is provided in a clear satin anodized finish. Versa-Duct multi-outlet assemblies save you time and money by reducing on-site labor costs by up to 70%. Whether designing a raceway system for your laboratory work bench or a surface mounted system for an entire building, Versa-Duct is the raceway for you. Cover with two end plates.

cUL Listed.

Description	Dimensions	No. of Outlets	Cat. No.
Single Channel Raceway, 4'	1 3/4" x 2 3/4" x 48"	2 Cover Mounted Spec Grade Duplex Receptacles	89422-924
Single Channel Raceway, 5'	1 3/4" x 2 3/4" x 60"	2 Cover Mounted Spec Grade Duplex Receptacles	89422-926
Single Channel Raceway, 6'	1 3/4" x 2 3/4" x 72"	3 Cover Mounted Spec Grade Duplex Receptacles	89422-928
Single Channel Raceway, 7'	1 3/4" x 2 3/4" x 84"	3 Cover Mounted Spec Grade Duplex Receptacles	89422-930
Single Channel Raceway, 8'	1 3/4" x 2 3/4" x 96"	4 Cover Mounted Spec Grade Duplex Receptacles	89422-932
Single Channel Raceway, 9'	1 3/4" x 2 3/4" x 108"	4 Cover Mounted Spec Grade Duplex Receptacles	89422-934
Single Channel Raceway, 10'	1 3/4" x 2 3/4" x 120"	5 Cover Mounted Spec Grade Duplex Receptacles	89422-936
Double Channel Raceway, 4'	2 1/4" x 5 1/4" x 48"	2 Cover Mounted Spec Grade Duplex Receptacles and 1 Data Holder with 2 Cat5 Jacks, Not Wired	89422-938
Double Channel Raceway, 5'	2 1/4" x 5 1/4" x 60"	2 Cover Mounted Spec Grade Duplex Receptacles and 1 Data Holder with 2 Cat5 Jacks, Not Wired	89422-940
Double Channel Raceway, 6'	2 1/4" x 5 1/4" x 72"	3 Cover Mounted Spec Grade Duplex Receptacles and 1 Data Holder with 2 Cat5 Jacks, Not Wired	89422-942
Double Channel Raceway, 7'	2 1/4" x 5 1/4" x 84"	3 Cover Mounted Spec Grade Duplex Receptacles and 1 Data Holder with 2 Cat5 Jacks, Not Wired	89422-944
Double Channel Raceway, 8'	2 1/4" x 5 1/4" x 96"	4 Cover Mounted Spec Grade Duplex Receptacles and 2 Data Holders with 2 Cat5 Jacks, Not Wired	89422-946
Double Channel Raceway, 9'	2 1/4" x 5 1/4" x 108"	4 Cover Mounted Spec Grade Duplex Receptacles and 2 Data Holders with 2 Cat5 Jacks, Not Wired	89422-948
Double Channel Raceway, 10'	2 1/4" x 5 1/4" x 120"	5 Cover Mounted Spec Grade Duplex Receptacles and 2 Data Holders with 2 Cat5 Jacks, Prewired	89422-950

### LabLink Power Strips, PG LifeLink

- Extruded from #6063-T5 aluminum alloy
- Dimensions 1¾" x 2¾" x 18" - 60"
- Ideal for mobile benches
- Single circuit configuration
- Clear satin anodized finish is durable and corrosion resistant
- 20 amp (NEMA 5-20R) grey duplex receptacles
- 8', 3 conductor, black power cord with straight blade plug (NEMA5-20P)

LabLink Power Strips are made of high quality extruded aluminum and can be permanently mounted to your work station. Power strips are available in lengths of 18", 24", 36", 48" and 60", providing more functionality for today's modular lab environment. Cover with 2 end plates.

C-UL Listed.

**Ordering Information:** Not to be used in conjunction with an extension cord or another power strip.



Description	Dimensions	No. of Outlets	Cat. No.
Power Strip, 1 Channel, 1.5'	1¾" x 2¾" x 18"	2 Cover Mounted Spec Grade Duplex Receptacle	89422-952
Power Strip, 1 Channel, 2'	1¾" x 2¾" x 24"	2 Cover Mounted Spec Grade Duplex Receptacle	89422-954
Power Strip, 1 Channel, 3'	1¾" x 2¾" x 36"	3 Cover Mounted Spec Grade Duplex Receptacles	89422-956
Power Strip, 1 Channel, 4'	1¾" x 2¾" x 48"	4 Cover Mounted Spec Grade Duplex Receptacles	89422-958
Power Strip, 1 Channel, 5'	1¾" x 2¾" x 60"	4 Cover Mounted Spec Grade Duplex Receptacles	89422-960

### Economy Shelf Bins, Quantum Storage Systems

- Edges are reinforced for added strength
- Designed for 12", 18" and 24" deep shelf units, racks, and standard shelving
- Wide hopper front allows complete access to the contents of the bin
- Available in six colors in the 4" height; blue, green, ivory, red, yellow, and black
- Available in red, blue or yellow in the 6" high bin
- Bin height available of 4" or 6"

Durable bins feature a mold-in label holder and built-in rear hanglock which allows bins to tilt out for complete access when on shelving. Seven divider slots.



Description	Exterior Dimensions	Interior Dimensions	Color	Cat. No.	Unit
4" High Shelf Bin	11½" x 4½" x 4"	10¼" x 3" x 4"	Ivory	10752-526	Case of 36
4" High Shelf Bin	11½" x 6½" x 4"	10¼" x 5½" x 4"	Ivory	10752-538	Case of 30
4" High Shelf Bin	11½" x 6½" x 4"	10¼" x 5½" x 4"	Red	10752-542	Case of 30
4" High Shelf Bin	11½" x 6½" x 4"	10¼" x 5½" x 4"	Yellow	10752-544	Case of 30
4" High Shelf Bin	11½" x 6½" x 4"	10¼" x 5½" x 4"	Black	10752-546	Case of 30
4" High Shelf Bin	11½" x 8¾" x 4"	10¼" x 7½" x 4"	Blue	10752-548	Case of 20
4" High Shelf Bin	11½" x 11½" x 4"	10¼" x 10" x 4"	Green	10752-562	Case of 8
4" High Shelf Bin	11½" x 11½" x 4"	10¼" x 10" x 4"	Yellow	10752-568	Case of 8
4" High Shelf Bin	17½" x 4½" x 4"	16½" x 3" x 4"	Green	10752-574	Case of 20
4" High Shelf Bin	17½" x 4½" x 4"	16½" x 3" x 4"	Ivory	10752-576	Case of 20
4" High Shelf Bin	17½" x 4½" x 4"	16½" x 3" x 4"	Red	10752-578	Case of 20
4" High Shelf Bin	17½" x 11½" x 4"	16½" x 10" x 4"	Black	10752-634	Case of 8
4" High Shelf Bin	11½" x 4½" x 4"	10¼" x 3" x 4"	Yellow	10752-530	Case of 36
4" High Shelf Bin	11½" x 8¾" x 4"	10¼" x 7½" x 4"	Green	10752-550	Case of 20
4" High Shelf Bin	11½" x 8¾" x 4"	10¼" x 7½" x 4"	Yellow	10752-556	Case of 20
4" High Shelf Bin	11½" x 4½" x 4"	10¼" x 3" x 4"	Blue	10752-522	Case of 36
4" High Shelf Bin	11½" x 4½" x 4"	10¼" x 3" x 4"	Red	10752-528	Case of 36
4" High Shelf Bin	11½" x 8¾" x 4"	10¼" x 7½" x 4"	Red	10752-554	Case of 20
4" High Shelf Bin	11½" x 11½" x 4"	10¼" x 10" x 4"	Blue	10752-560	Case of 8
4" High Shelf Bin	11½" x 11½" x 4"	10¼" x 10" x 4"	Red	10752-566	Case of 8
4" High Shelf Bin	17½" x 4½" x 4"	16½" x 3" x 4"	Black	10752-598	Case of 20
4" High Shelf Bin	17½" x 6½" x 4"	16½" x 5½" x 4"	Green	10752-602	Case of 20

## Economy Shelf Bins, Quantum Storage Systems (cont.)

- Edges are reinforced for added strength
- Designed for 12", 18" and 24" deep shelf units, racks, and standard shelving
- Wide hopper front allows complete access to the contents of the bin
- Available in six colors in the 4" height; blue, green, ivory, red, yellow, and black
- Available in red, blue or yellow in the 6" high bin
- Bin height available of 4" or 6"

Durable bins feature a mold-in label holder and built-in rear hanglock which allows bins to tilt out for complete access when on shelving. Seven divider slots.



Description	Exterior Dimensions	Interior Dimensions	Color	Cat. No.	Unit
4" High Shelf Bin	17 <sup>7</sup> / <sub>8</sub> " x 8 <sup>3</sup> / <sub>8</sub> " x 4"	16 <sup>1</sup> / <sub>2</sub> " x 7 <sup>1</sup> / <sub>8</sub> " x 4"	Green	10752-614	Case of 10
4" High Shelf Bin	17 <sup>7</sup> / <sub>8</sub> " x 8 <sup>3</sup> / <sub>8</sub> " x 4"	16 <sup>1</sup> / <sub>2</sub> " x 7 <sup>1</sup> / <sub>8</sub> " x 4"	Ivory	10752-616	Case of 10
4" High Shelf Bin	17 <sup>7</sup> / <sub>8</sub> " x 8 <sup>3</sup> / <sub>8</sub> " x 4"	16 <sup>1</sup> / <sub>2</sub> " x 7 <sup>1</sup> / <sub>8</sub> " x 4"	Black	10752-622	Case of 10
4" High Shelf Bin	17 <sup>7</sup> / <sub>8</sub> " x 11 <sup>1</sup> / <sub>8</sub> " x 4"	16 <sup>1</sup> / <sub>2</sub> " x 10" x 4"	Ivory	10752-628	Case of 8
4" High Shelf Bin	17 <sup>7</sup> / <sub>8</sub> " x 6 <sup>5</sup> / <sub>8</sub> " x 4"	16 <sup>1</sup> / <sub>2</sub> " x 5 <sup>1</sup> / <sub>2</sub> " x 4"	Yellow	10752-664	Case of 20
4" High Shelf Bin	11 <sup>5</sup> / <sub>8</sub> " x 4 <sup>1</sup> / <sub>8</sub> " x 4"	10 <sup>3</sup> / <sub>4</sub> " x 3" x 4"	Black	10752-532	Case of 36
4" High Shelf Bin	11 <sup>5</sup> / <sub>8</sub> " x 8 <sup>3</sup> / <sub>8</sub> " x 4"	10 <sup>1</sup> / <sub>4</sub> " x 7 <sup>1</sup> / <sub>8</sub> " x 4"	Ivory	10752-552	Case of 20
4" High Shelf Bin	11 <sup>5</sup> / <sub>8</sub> " x 8 <sup>3</sup> / <sub>8</sub> " x 4"	10 <sup>1</sup> / <sub>4</sub> " x 7 <sup>1</sup> / <sub>8</sub> " x 4"	Black	10752-558	Case of 20
4" High Shelf Bin	11 <sup>5</sup> / <sub>8</sub> " x 11 <sup>1</sup> / <sub>8</sub> " x 4"	10 <sup>1</sup> / <sub>4</sub> " x 10" x 4"	Ivory	10752-564	Case of 8
4" High Shelf Bin	11 <sup>5</sup> / <sub>8</sub> " x 11 <sup>1</sup> / <sub>8</sub> " x 4"	10 <sup>1</sup> / <sub>4</sub> " x 10" x 4"	Black	10752-570	Case of 8
4" High Shelf Bin	17 <sup>7</sup> / <sub>8</sub> " x 6 <sup>5</sup> / <sub>8</sub> " x 4"	16 <sup>1</sup> / <sub>2</sub> " x 5 <sup>1</sup> / <sub>2</sub> " x 4"	Blue	10752-600	Case of 20
4" High Shelf Bin	17 <sup>7</sup> / <sub>8</sub> " x 8 <sup>3</sup> / <sub>8</sub> " x 4"	16 <sup>1</sup> / <sub>2</sub> " x 7 <sup>1</sup> / <sub>8</sub> " x 4"	Red	10752-618	Case of 10
4" High Shelf Bin	17 <sup>7</sup> / <sub>8</sub> " x 8 <sup>3</sup> / <sub>8</sub> " x 4"	16 <sup>1</sup> / <sub>2</sub> " x 7 <sup>1</sup> / <sub>8</sub> " x 4"	Yellow	10752-620	Case of 10
4" High Shelf Bin	17 <sup>7</sup> / <sub>8</sub> " x 6 <sup>5</sup> / <sub>8</sub> " x 4"	16 <sup>1</sup> / <sub>2</sub> " x 5 <sup>1</sup> / <sub>2</sub> " x 4"	Ivory	10752-644	Case of 20
4" High Shelf Bin	17 <sup>7</sup> / <sub>8</sub> " x 6 <sup>5</sup> / <sub>8</sub> " x 4"	16 <sup>1</sup> / <sub>2</sub> " x 5 <sup>1</sup> / <sub>2</sub> " x 4"	Red	10752-662	Case of 20
4" High Shelf Bin	17 <sup>7</sup> / <sub>8</sub> " x 6 <sup>5</sup> / <sub>8</sub> " x 4"	16 <sup>1</sup> / <sub>2</sub> " x 5 <sup>1</sup> / <sub>2</sub> " x 4"	Black	10752-666	Case of 20
4" High Shelf Bin	17 <sup>7</sup> / <sub>8</sub> " x 11 <sup>1</sup> / <sub>8</sub> " x 4"	16 <sup>1</sup> / <sub>2</sub> " x 10" x 4"	Blue	10752-624	Case of 8
4" High Shelf Bin	11 <sup>5</sup> / <sub>8</sub> " x 4 <sup>1</sup> / <sub>8</sub> " x 4"	10 <sup>3</sup> / <sub>4</sub> " x 3" x 4"	Green	10752-524	Case of 36
4" High Shelf Bin	17 <sup>7</sup> / <sub>8</sub> " x 11 <sup>1</sup> / <sub>8</sub> " x 4"	16 <sup>1</sup> / <sub>2</sub> " x 10" x 4"	Green	10752-626	Case of 8
4" High Shelf Bin	17 <sup>7</sup> / <sub>8</sub> " x 4 <sup>1</sup> / <sub>8</sub> " x 4"	16 <sup>1</sup> / <sub>2</sub> " x 3" x 4"	Blue	10752-572	Case of 20
4" High Shelf Bin	17 <sup>7</sup> / <sub>8</sub> " x 8 <sup>3</sup> / <sub>8</sub> " x 4"	16 <sup>1</sup> / <sub>2</sub> " x 7 <sup>1</sup> / <sub>8</sub> " x 4"	Blue	10752-612	Case of 10
4" High Shelf Bin	17 <sup>7</sup> / <sub>8</sub> " x 11 <sup>1</sup> / <sub>8</sub> " x 4"	16 <sup>1</sup> / <sub>2</sub> " x 10" x 4"	Red	10752-630	Case of 8
4" High Shelf Bin	11 <sup>5</sup> / <sub>8</sub> " x 6 <sup>3</sup> / <sub>8</sub> " x 4"	10 <sup>3</sup> / <sub>4</sub> " x 5 <sup>1</sup> / <sub>2</sub> " x 4"	Green	10752-536	Case of 30
4" High Shelf Bin	17 <sup>7</sup> / <sub>8</sub> " x 11 <sup>1</sup> / <sub>8</sub> " x 4"	16 <sup>1</sup> / <sub>2</sub> " x 10" x 4"	Yellow	10752-632	Case of 8
4" High Shelf Bin	11 <sup>5</sup> / <sub>8</sub> " x 6 <sup>5</sup> / <sub>8</sub> " x 4"	10 <sup>3</sup> / <sub>4</sub> " x 5 <sup>1</sup> / <sub>2</sub> " x 4"	Blue	10752-534	Case of 30
4" High Shelf Bin	17 <sup>7</sup> / <sub>8</sub> " x 4 <sup>1</sup> / <sub>8</sub> " x 4"	16 <sup>1</sup> / <sub>2</sub> " x 3" x 4"	Yellow	10752-580	Case of 20
6" High Shelf Bin	11 <sup>5</sup> / <sub>8</sub> " x 4 <sup>1</sup> / <sub>8</sub> " x 6"	10 <sup>3</sup> / <sub>8</sub> " x 3" x 6"	Red	10752-638	Case of 36
6" High Shelf Bin	17 <sup>7</sup> / <sub>8</sub> " x 6 <sup>3</sup> / <sub>8</sub> " x 6"	16 <sup>1</sup> / <sub>2</sub> " x 5 <sup>1</sup> / <sub>2</sub> " x 6"	Blue	10752-714	Case of 20
6" High Shelf Bin	17 <sup>7</sup> / <sub>8</sub> " x 6 <sup>5</sup> / <sub>8</sub> " x 6"	16 <sup>1</sup> / <sub>2</sub> " x 5 <sup>1</sup> / <sub>2</sub> " x 6"	Yellow	10752-718	Case of 20
6" High Shelf Bin	11 <sup>5</sup> / <sub>8</sub> " x 4 <sup>1</sup> / <sub>8</sub> " x 6"	10 <sup>3</sup> / <sub>8</sub> " x 3" x 6"	Yellow	10752-640	Case of 36
6" High Shelf Bin	11 <sup>5</sup> / <sub>8</sub> " x 6 <sup>3</sup> / <sub>8</sub> " x 6"	10 <sup>3</sup> / <sub>8</sub> " x 5 <sup>1</sup> / <sub>2</sub> " x 6"	Yellow	10752-670	Case of 30
6" High Shelf Bin	11 <sup>5</sup> / <sub>8</sub> " x 8 <sup>3</sup> / <sub>8</sub> " x 6"	10 <sup>3</sup> / <sub>8</sub> " x 7 <sup>1</sup> / <sub>8</sub> " x 6"	Yellow	10752-676	Case of 20
6" High Shelf Bin	11 <sup>5</sup> / <sub>8</sub> " x 11 <sup>1</sup> / <sub>8</sub> " x 6"	10 <sup>3</sup> / <sub>8</sub> " x 10" x 6"	Yellow	10752-682	Case of 8
6" High Shelf Bin	17 <sup>7</sup> / <sub>8</sub> " x 8 <sup>3</sup> / <sub>8</sub> " x 6"	16 <sup>1</sup> / <sub>2</sub> " x 7 <sup>1</sup> / <sub>8</sub> " x 6"	Blue	10752-720	Case of 10
6" High Shelf Bin	11 <sup>5</sup> / <sub>8</sub> " x 6 <sup>5</sup> / <sub>8</sub> " x 6"	10 <sup>3</sup> / <sub>8</sub> " x 5 <sup>1</sup> / <sub>2</sub> " x 6"	Blue	10752-642	Case of 30
6" High Shelf Bin	11 <sup>5</sup> / <sub>8</sub> " x 6 <sup>3</sup> / <sub>8</sub> " x 6"	10 <sup>3</sup> / <sub>8</sub> " x 5 <sup>1</sup> / <sub>2</sub> " x 6"	Red	10752-668	Case of 30
6" High Shelf Bin	11 <sup>5</sup> / <sub>8</sub> " x 8 <sup>3</sup> / <sub>8</sub> " x 6"	10 <sup>3</sup> / <sub>8</sub> " x 7 <sup>1</sup> / <sub>8</sub> " x 6"	Red	10752-674	Case of 20
6" High Shelf Bin	11 <sup>5</sup> / <sub>8</sub> " x 11 <sup>1</sup> / <sub>8</sub> " x 6"	10 <sup>3</sup> / <sub>8</sub> " x 10" x 6"	Red	10752-680	Case of 8
6" High Shelf Bin	17 <sup>7</sup> / <sub>8</sub> " x 8 <sup>3</sup> / <sub>8</sub> " x 6"	16 <sup>1</sup> / <sub>2</sub> " x 7 <sup>1</sup> / <sub>8</sub> " x 6"	Yellow	10752-724	Case of 10
6" High Shelf Bin	17 <sup>7</sup> / <sub>8</sub> " x 11 <sup>1</sup> / <sub>8</sub> " x 6"	16 <sup>1</sup> / <sub>2</sub> " x 10" x 6"	Yellow	10752-736	Case of 8
6" High Shelf Bin	17 <sup>7</sup> / <sub>8</sub> " x 8 <sup>3</sup> / <sub>8</sub> " x 6"	16 <sup>1</sup> / <sub>2</sub> " x 7 <sup>1</sup> / <sub>8</sub> " x 6"	Red	10752-722	Case of 10
6" High Shelf Bin	17 <sup>7</sup> / <sub>8</sub> " x 11 <sup>1</sup> / <sub>8</sub> " x 6"	16 <sup>1</sup> / <sub>2</sub> " x 10" x 6"	Blue	10752-728	Case of 8
6" High Shelf Bin	17 <sup>7</sup> / <sub>8</sub> " x 4 <sup>1</sup> / <sub>8</sub> " x 6"	16 <sup>1</sup> / <sub>2</sub> " x 3" x 6"	Blue	10752-708	Case of 20
6" High Shelf Bin	17 <sup>7</sup> / <sub>8</sub> " x 4 <sup>1</sup> / <sub>8</sub> " x 6"	16 <sup>1</sup> / <sub>2</sub> " x 3" x 6"	Yellow	10752-712	Case of 20
6" High Shelf Bin	11 <sup>5</sup> / <sub>8</sub> " x 8 <sup>3</sup> / <sub>8</sub> " x 6"	10 <sup>3</sup> / <sub>8</sub> " x 7 <sup>1</sup> / <sub>8</sub> " x 6"	Blue	10752-672	Case of 20
6" High Shelf Bin	17 <sup>7</sup> / <sub>8</sub> " x 6 <sup>5</sup> / <sub>8</sub> " x 6"	16 <sup>1</sup> / <sub>2</sub> " x 5 <sup>1</sup> / <sub>2</sub> " x 6"	Red	10752-716	Case of 20
6" High Shelf Bin	17 <sup>7</sup> / <sub>8</sub> " x 4 <sup>1</sup> / <sub>8</sub> " x 6"	16 <sup>1</sup> / <sub>2</sub> " x 3" x 6"	Red	10752-710	Case of 20
6" High Shelf Bin	17 <sup>7</sup> / <sub>8</sub> " x 11 <sup>1</sup> / <sub>8</sub> " x 6"	16 <sup>1</sup> / <sub>2</sub> " x 10" x 6"	Red	10752-734	Case of 8
6" High Shelf Bin	11 <sup>5</sup> / <sub>8</sub> " x 4 <sup>1</sup> / <sub>8</sub> " x 6"	10 <sup>3</sup> / <sub>8</sub> " x 3" x 6"	Blue	10752-636	Case of 36
6" High Shelf Bin	11 <sup>5</sup> / <sub>8</sub> " x 11 <sup>1</sup> / <sub>8</sub> " x 6"	10 <sup>3</sup> / <sub>8</sub> " x 10" x 6"	Blue	10752-678	Case of 8



# The VWR Furniture Catalog

**VWR Furniture** | designed • delivered • installed

Whether you're looking to build a new lab or renovate an existing one—check out VWR's Furniture Catalog!

Find innovative storage ideas, modular benches, fixed casework, and functional seating and accessories – all from one quick, convenient source. VWR offers a full line of quality lab furniture and components, along with a complete turn-key service offering that can help you with lab planning, design, and installation.



Go to [vwr.com/furniture](http://vwr.com/furniture) today to request your free printed copy or access the online, interactive version.

**vwr.com**  
**1.800.932.5000**

## VWR Furniture

designed • delivered • installed

Mobile Laboratory Benches

Worksurfaces

Suspended Cabinets

Shelves

Overhead Cabinets

Fixtures/Power Bars

Technical Sheets

Accessories



**1.800.932.5000 | vwr.com**



Prices and product details are current when published; subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks, of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). | Visit vwr.com to view our privacy policy, trademark owners and additional disclaimers. ©2018 VWR International, LLC. All rights reserved.

11505\_LM\_JS\_2017

1217 3M Lit. No. 930231