



CERTIFICATE OF ANALYSIS

CiDehol[®]70

ISOPROPYL ALCOHOL SOLUTION

Lot # **A09042004E**

Manufactured September 2020

Expiration September 2022

Decon 8432 32 oz. Size

Meets Current USP Specifications

TEST	SPECIFICATION	UNIT	TEST RESULT
Appearance	Clear, colorless liquid	N/A	Pass
Isopropyl alcohol, % v/v	68.0 – 72.0	%	70.4
Specific gravity at 20°C	0.872 - 0.883	N/A	0.876
Limit of non-volatile residues	< 5	mg	< 5
Acidity	<0.02 meq/50ml	meq/50ml	<0.01

The above material has been filtered through a Hydrophobic Cartridge filter, pore size 0.2 µm.

Reviewed and Approved by:

Date: September 2020

A handwritten signature in blue ink, appearing to read 'Karen Wend'.

Karen Wend
Manager, Quality and Technical Services

Updated November 16, 2011- USP note

BSE/TSE Statement: There are no materials of animal origin in the raw materials or in the process aides used to manufacture any Decon products.

460 Glennie Circle • King of Prussia, PA 19406 • (610) 755-0800 • Fax (610) 270-8905